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mail: techmed@techmed.pl

Poland, Europe

**USER MANUAL
INFRARED THERMOMETER NON-CONTACT
MEASUREMENT**

MODEL.: TMB-70 EXP

NOTE: THE USER MUST CHECK THAT THE EQUIPMENT FUNCTIONS SAFELY AND SEE THAT IT IS IN PROPER WORKING CONDITION BEFORE BEING USED.

CONTENTS

- I. Features.....1
- II. Intended use..... 1
- III. Safety Precautions..... 1
- IV. How The Non-contact Infrared Thermometer Works.....2
- V. Normal Temperatures According To Age.....2
- VI. How To Take A Temperature.....2
- VII. Basic Instrument.....2
- VIII. Instructions..... 3
- IX. Setting And Function Of Menu..... 3
- X. Battery Low Voltage Detection..... 4
- XI. Changing the Battery.....4
- XII. Technical Specifications.....4
- XIII. Maintenance of The Product..... 5
- XIV. Accessories..... 5
- XV. Guidelines..... 5
- XVI. Troubleshooting..... 5
- XVII. Explanation of symbols 6
- XVIII. EMC Declaration.....6

THE MANUFACTURER RESERVES THE RIGHT TO ALTER THE SPECIFICATIONS OF THE PRODUCT WITHOUT PRIOR NOTIFICATION.

- I. Features**
1. Special design to take the Human Body Temperature with a 3~5cm (1.2~2 in) distance from forehead.
 2. Reliable and stable measurement, thanks to the advantage Infrared Detection System.
 3. Audible alarm if temperature is more than 38°C (100.4°F).
 4. Memorize the last 32 temperature measurements.
 5. Three color backlits LCD digital display screen.
 6. Automatic power-off (< 30 secs) to conserve energy.
 7. Longevity use (100,000 readings) .
 8. Practical, easy to use.
 9. Use direct mode to measure

II. Intended use

The device is an infrared thermometer intended to measure forehead temperature of infants and adults without contacting human body. It can be used by consumers in household environment and doctor in clinic as reference.

- III. Safety Precautions**
- Follow the maintenance advice stipulated in this instruction manual.
 - This device may be used for professional purposes or for personal home use.
 - This device must only be used for the purposes described in this instruction manual.
 - This device must only be used in an ambient temperature range of between 10°C and 40°C.
 - This device must be always kept in a clean, dry area.
 - Do not expose this thermometer to electric shocks.
 - Do not expose this thermometer to extreme temperature conditions of >55°C or <-20°C.
 - Do not use this device in relative humidity higher than 85%.
 - Do not touch the glass of the infrared sensor with your fingers.
 - Clean the glass with a cotton bud lightly moistened with 95° alcohol.
 - Do not expose the thermometer to sunlight or to water.
 - Never drop the device.
 - Should a problem occur with your device, please contact your retailer. Do not attempt to repair this device yourself.
 - The device is not suitable for use in the presence of flammable anesthetic mixtures with oxygen or nitrous oxide.
 - The device requires no calibration.
 - The device is not repairable and contains no user serviceable parts.
 - Manufacturer will provide circuit diagrams, component part lists, descriptions, calibration instruction to assist to SERVICE PERSONNEL in parts repair.
 - Modification of this equipment is not allowed.
 - Remove batteries if equipment is not likely to be used for some time.
 - The user must check that the equipment functions safely and see that it is in proper working condition before being used.
 - The manufacturer does not require such preventive inspections by other persons.
 - Contact the manufacturer for assistance or reporting.
 - Keep the equipment away from children and pets.
 - Do not dispose of electrical appliances as unsorted municipal waste, use separate collection facilities.
- Contact you local government for information regarding the collection systems available. If electrical appliances are disposed of in landfills or dumps, hazardous substances can leak into the groundwater and get into the food chain, damaging your health and well-being.

It is essential to use the Non-contact Infrared Thermometer Model TMB-70 EXP. You are therefore advised to read this instruction manual and the safety precautions carefully before use.

The TMB-70 EXP Non-contact Infrared Thermometer is pre-set at the factory.

It is not necessary to calibrate the device when starting it up. However, in order to obtain reliable results, you are advised to allow the thermometer to acclimate to room/ambient temperature for 15 to 20

minutes before using each time there is a significant change in environmental temperature. It is also important to allow a 3 to 5 seconds intervals between measurements.

IV. How the Non-contact Infrared Thermometer Model TMB-70 EXP Works.

All objects, solid, liquid or gas, emit energy by radiation. The intensity of this energy depends on the temperature of the object. The TMB-70 EXP infrared thermometer is therefore able to measure the temperature of a person by the energy the person emits. This measurement can be taken thanks to an external temperature probe on the device which permanently analyses and registers the ambient temperature. Therefore, as soon as the operator holds the thermometer near the body and activates the radiation sensor, the measurement is taken instantly by detection of the infrared heat generated by the arterial blood flow. Therefore body heat can be measured without any interference from the heat of the surrounding environment.

Use the Non-contact Infrared Thermometer Model TMB-70 EXP:

- When you have reason to believe you are sick.
- When you are recovering from surgery, illness or exhaustion.
- For monitoring body temperatures.
- When exercising, hiking, or doing strenuous physical activity.
- When traveling or when you are under abnormal stress.


V. Normal Temperatures According To Age

Age	°C	°F
0-2 years	36.4-38.0	97.5-100.4
3-10 years	36.1-37.8	97.0-100.0
11-65 years	35.9-37.6	96.6-99.7
> 65 years	35.8-37.5	96.4-99.5

Manufacturer does not take responsibility for the results or actions of non- medical personnel using the Non-contact infrared thermometer Model TMB-70 EXP.

VI. How to Take a Temperature



Aim at the FOREHEAD, over the right temporal region, from a distance of 3cm~5cm (1.2in~2in), press the scan button to display temperature instantly. (as fig.1)

Warning: The reliability of the measurement cannot be guaranteed if the temperature is measured over another part of the body (e.g. arm, torso...)

VII. Basic Instrument

The type BF applied part: Sensor.

- VII. Instructions**
1. Install battery.
 2. For the first use or when inserting new battery, wait from 10~15 minutes for the warm-up of the unit. This will allow the unit to become acclimated to the temperature of the room.
 3. Press "On/Scan button to turn on the unit, aim towards the desired location, keeping the unit at a distance of 3cm~5cm (1.2in~2in) from forehead, and press the "On/Scan" button. Results will be displayed in 1 second.
 4. Before taking the temperature, make sure to remove hair and perspiration from the forehead. You can also take the temperature

IX. Setting And Function of Menu

1. Switch on the device

Press the "On/Scan" button, the screen panel in full display, it will enter the standby mode with the sign "—"°C".When then press the "On/Scan" button will get the measuring result in 1 second. But if there is no key to press, it will turn off in 30 seconds automatically.

- 2. Mode choice**
- a) The TMB-70 EXP is specially designed to take the body temperature of a human being. For this use the BODY mode. Measurement range for BODY mode: 32.0°C - 43.0°C (89.6°F - 109.4°F).
 - b) You can use the TMB-70 EXP to measure the temperature of an area or an object, a food, a liquid. For this, use the SURFACE TEMP mode. Measurement range for SURFACE TEMP mode: 0°C-60°C (32°F - 140°F).
 - c) You also can use the TMB-70 EXP to measure the temperature of a room. For this use the ROOM mode. Measurement range for ROOM mode: 0°C - 40°C (32°F - 104°F).
- In the switch on state, press the "MODE" button. You can freely select "BODY", "SURFACE TEMP" or "ROOM" mode.
- ※ **Note:** The TMB-70 EXP is default setting as BODY mode.
- ※ **Note:** The surface temperature is different from the internal human temperature. To obtain the internal human temperature always use the BODY mode.

3. F1: Choosing the temperature unit

In the switch on state, Press "MODE" button for 3 seconds, the screen will display "F-1", then press "MEM" button for Celsius and Select "On/Scan" button for Fahrenheit, Confirm by press "MODE" button. And press "MODE" button three times to exit the Settings.

4. F2: Alarm setup

In the switch on state, Press "MODE" button for 3 seconds, the screen will display "F-1", then press "MODE" button once, the screen will display "F-2", press "On/Scan" button to add 0.1°C. Select the "MEM" button to minus 0.1°C. Confirm by pressing "MODE" button. And press "MODE" button two times to exit the Settings.

Note: For temperature unit Celsius, the alarm range is 37.0°C to 39.0°C.

Note: High temperature alarm default setting is 38°C (100.4°F).

5. Recalibration of device via the F4 MENU

When there is the difference between TMB-70 EXP and mercury thermometer, and you believe mercury thermometer from its temperature but it is not convenient to use. You can use recalibration function to adjust the TMB-70 EXP to make it the same test result with mercury thermometer after recalibration.

Besides, when you use TMB-70 EXP for the people with different skin color (For example: the yellow race, the white race, black people and so on) you can use recalibration too.

Instructions for recalibration:

In the switch on state, Press "MODE" button for 3 seconds, the screen will display "F-1", then press "MODE" button twice, the screen will display "F-4", Press "On/Scan" button to increase 0.1°C, press "MEM" button to reduce 0.1°C, Confirm by pressing "MODE" button. And press "MODE" button one time to exit the Settings.

Note: Default setting as 0.0°C.

Note: For temperature unit Celsius, the alarm range is 3.0°C to -3.0°C.

6. F5: Buzzer ON/OFF

In the switch on state, Press "MODE" button for 3 seconds, the screen will display "F-1", then press "MODE" button thrice, the screen will display "F-5", Press "On/Scan" button to open the beeper (a symbol "🔊" flashes on the screen and shows ON); press the "MEM" button to close the beeper (a symbol "🔊" disappear and shows OFF), Confirm by pressing "MODE" button and exit the settings.

7. Date memory

In the switch on state, press "MEM" button, which will display the last temperature, and allows for a view of the last 32 measurements taken.

8. Delete the memory

In the switch on state, press the "MODE" and the "MEM" buttons at the same time for 5 seconds, will delete all memory data. Then press "MEM" button again will display "CLr".

6. High Temperature Alarm

In BODY mode, if the temperature is less than the setting value (the default value is 38°C), The buzzer sends out 1 sound. When the temperature is more than or equal to the alarm temperature, the buzzer sends out 5 sounds.


When the ambient temperature is less than 10°C or more than 40°C, measuring in BODY mode or SURFACE TEMP mode, the LCD displays ERR1, at the same time, sends out 2 sounds.

X. Battery Low Voltage Detection

The TMB-70 EXP Non-contact Infrared Thermometer will display the flashed symbol "🔋" to remind you to replace the battery.

XI. Changing the Battery

This TMB-70 EXP is with 2 AAA alkaline batteries for 20,000 times use. When the LCD screen displays "🔋", the battery should be replaced. Open the lid and change the battery, taking great care with the correct positioning. A mistake with this could cause damage to the unit and compromise the guarantee of your Non-contact Infrared Thermometer. (as fig.2)

XII. Technical Specifications

- Normal using condition
Ambient temperature: 10°C ~ 40°C (50°F ~ 104°F)
Relative humidity: ≤ 85%
Pressure altitude: 700hPa ~ 1060hPa
- Storage and shipping condition
Ambient temperature: -20°C ~ 55°C (-4°F ~ 131°F)
Relative humidity: ≤ 95%
Pressure altitude: 700hPa ~ 1060hPa
- Batteries: DC 3V
- Unit size: 170x42x39mm (6.70x1.65x1.54in) L x W x H
- Unit weight: 80g (without Batteries)
- Temperature display resolution: 0.1°C (0.1°F)
- Measuring range:

In body mode: 32.0°C ~ 43.0°C (89.6°F ~ 109.4°F)
Under body mode, there is three color backlights:
Green color backlight: ≤37.3°C (99.1°F), means normal temperature.
Orange color backlight: 37.4°C~37.9°C (99.3°F ~ 100.2°F), means low fever.
Red color backlight: ≥38°C (100.4°F), means high fever.
In surface temp mode: 0°C ~ 60°C (32°F ~ 140°F)
In room mode: 0°C ~ 40°C (32°F ~ 104°F)

8. Precision:

32.0°C ~ 34.9°C (89.6°F ~ 94.8°F) ±0.3°C(±0.6°F)
35.0°C ~ 42.0°C (95°F ~ 107.6°F) ±0.2°C(±0.4°F)
42.1°C ~ 43.0°C (107.8°F ~ 109.4°F) ±0.3°C(±0.6°F)

9. Consumption: ≤ 300mW

10. Measuring distance: 3cm ~ 5 cm (1.2 ~ 2 in)

11. Automatic power-off: < 30 secs.

12. Memory: 32 sets

※Note: The Non-contact Infrared Thermometer Model TMB-70 EXP can take temperature readings below 32°C or above 43.0°C (89.6°F to 109.4°F) but precision is not guaranteed outside of this range.

Longevity of the Product:

The Non-contact Infrared Thermometer Model TMB-70 EXP was conceived for an intense and professional use. Its longevity is guaranteed for 100,000 takings.

XIII. Maintenance of the Product

- The infrared sensor is the most important and fragile part of the thermometer. Take care when using, storing or cleaning.
- Clean infrared sensor with cotton fabric dampened with a 95% alcohol solution.
- Do not throw battery in fire. Please dispose of battery properly.
- Remove the battery when thermometer is not used for an extended period of time.
- Do not expose the thermometer to water or excessive sunlight.
- An impact may damage the product.

XIV. Accessories

- User manual: 1 pc
- AAA Alkaline batteries: 2 pcs

XV. Guidelines

This device complies with the EU Directive 93/42/EEC concerning medical products, the ISO 80601-2-56 and the European Standard EN60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility.


XVI. Troubleshooting

If you have problems while using your thermometer, please refer to this guide to help resolve the problem. If the problem persists, please contact our customer service.

To take the surface temperature, press the "MODE" button and set to the reading called "Body". If the device is in Surface Temp Mode, the 32°C (89.6°F) temperature displayed is showing the external temperature of your body, rather than the internal.


THE SCREEN DISPLAYS THE MESSAGE HI

When using the Non-Contact Infrared Thermometer Model TMB-70 EXP the message "HI" can show on the screen. In this case, the temperature is either above the measurement range selected, or above 43.0°C (109.4°F) in Body Mode.

THE SCREEN DISPLAYS THE MESSAGE LO

When using the Non-contact Infrared Thermometer Model TMB-70 EXP, the message "Lo" can show on the screen. In this case, the temperature analyzed is either under the measuring range selected, or less than 32°C (89.6°F) in Body Mode.

This message displays for various reasons. Please find below a list of the main issues:

Reasons for "Lo" message display	Advice
Temperature reading hampered by hair or perspiration.	Make sure there is no obstruction or dampness prior to taking temperature.
Temperature hampered by an air draft or dramatic change in ambient temperature.	Make sure there is no air blowing in the area of use; this could affect the infrared reading.
Temperature readings are too close together, and the thermometer did not have time to reboot.	Pause for 3~5 seconds minimum between readings; a 15 seconds pause is recommended.
The measuring distance is too far.	Take measurements at the recommended distance (app. 3cm~5cm; 1.2in~2.0 in).

XVII. Explanation of symbols

Symbol	Reference
	IEC 60417-5333, Type BF applied part
	IEC 60417-5031 Direct current
IP22	Protected against access to hazardous parts with a finger and against vertically falling water drops when enclosure tilted up to 15°
	Refer to instruction manual / booklet
	DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.
	This symbol shall be accompanied by the name and the address of the manufacturer
SN	Specifies serial number

XVIII. EMC Declaration

Guidance and manufacturer's declaration – electromagnetic immunity

The TMB-70 EXP is intended for use in the electromagnetic environment specified below. The customer or the user of the TMB-70 EXP should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines	Not Applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not Applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	Not Applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the TMB-70 EXP requires continued operation during power interruptions, it is recommended that the TMB-70 EXP be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity			
The TMB-70 EXP is intended for use in the electromagnetic environment specified below. The customer or the user of the TMB-70 EXP should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not Applicable	Portable and mobile RF should communications be used no closer equipment to part cables, of the than the TMB-70 EXP, recommended including any separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1.167\sqrt{P}$ $d=1.167\sqrt{P}$ 80MHz to 800MHz $d=2.333\sqrt{P}$ 800MHz to 2.5 GHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b Interference may occur in the vicinity of equipment marked with the following symbol:
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a	Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio To broadcast assess and the TV electromagnetic broadcast cannot environment be predicted due to fixed theoretically RF with transmitters, accuracy. an electromagnetic site survey should be considered. If the measured field strength in normal the compliance location operation. level in If which above, abnormal the the TMB-70 EXP performance TMB-70 EXP is should exceeds observed, be the additional observed applicable to measure verify fly s may be necessary, such as reorienting or relocating the TMB-70 EXP. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.		
b			

Recommended separation distances between portable and mobile RF communications equipment and the Medical TMB-70 EXP

The TMB-70 EXP is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Medical TMB-70 EXP can help prevent electromagnetic interference by maintaining a minimum and distance the TMB-70 EXP between as portable and recommended mobile below, RF according communications to the equipment maximum output (transmitters) power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=1.167\sqrt{P}$	80 MHz to 800 MHz $d=1.167\sqrt{P}$	800 MHz to 2,5 GHz $d=2.333\sqrt{P}$
0,01	0.117	0.117	0.233
0,1	0.369	0.369	0.738
1	1.167	1.167	2.333
10	3.690	3.690	7.377
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Disposal

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

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