
# BLOOD PRESSURE MONITOR ARDOR

## MANUAL ECB004

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Guangdong Richjinn Medical Technology Co., Ltd.  
201, Building 38, Xi Tian Second Industrial Park, GongMing Street,  
Guangming District, Shenzhen, 518107 Guangdong, P.R.China

0197

## BLOOD PRESSURE MONITOR

### WHAT IS BLOOD PRESSURE?


The blood pressure is the pressure the blood is exerting on the walls of the arteries. The blood pressure in the arteries is subject to continuous variation during a heart cycle. The highest pressure within this cycle is called 'systolic pressure' (SYS), the lowest pressure is called 'diastolic pressure' (DIA). With these two blood pressure values, a doctor can determine a patient's blood pressure. Your blood pressure may be influenced by many different factors as physical effort, anxiety or the time of the day etc. Generally, blood pressure is low in the morning and raises in the afternoon towards the evening. In summer, blood pressure is lower than in winter.

### WHY IS IT A GOOD THING TO MEASURE BLOOD PRESSURE AT HOME

Having one's blood pressure measured by a doctor in a hospital or clinic environment tends to stimulate nervousness in a patient and may lead to a higher blood pressure. Blood pressure varies in accordance with a variety of conditions and so judgment is not always possible on the basis of a single measurement. Blood pressure measured first thing in the morning after getting up, before eating, and keeping still, is known as the fundamental blood pressure. In practice it is rather difficult to record the fundamental blood pressure but measurement at home in the right environment aids this.


### VARIATIONS OF BLOOD PRESSURE

Individual blood pressure is subject to variations on a day-to-day basis or in relation to the time of the year. These variations are more perceptible in patients with high blood pressure. Usually, blood pressure is higher during work and is lowest in the night when sleeping. The figure shows the variations of the blood pressure in the course of a day with measurements taken every 5 minutes. The bold line represents the sleep phase. The increase in blood pressure at 4 pm (shown as letter 'A' in the figure) and at midnight (shown as letter 'B' in the figure) mark pain events of the patient.

### BLOOD PRESSURE CLASSIFICATION OF THE WHO

The standards for the determination of high or low blood pressure independently from the patient's age – as shown on the table – were determined by the World Health Organization (WHO).


## IMPORTANT INSTRUCTIONS OF THE DEVICE

- Blood pressure should be judged by a medical doctor or other medical specialist who is familiar with your case history. Through the regular use of the blood pressure gauge you can continuously record your pressure values. These records can be used for therapeutic purposes.
- Tightly wrap the cuff around your upper arm. The cuff must be positioned on the level of the heart.
- The monitor must not be moved or shaken during measuring as this would lead to inaccurate results.
- Carry out measuring in a still and relaxed position.
- Do not wrap the cuff around the sleeve of a jacket, shirt which could make measuring impossible.
- Keep in mind that blood pressure naturally varies in the course of the day and is influenced by many different factors as, e.g., smoking, alcohol, drugs and physical effort. We recommend not to smoke, eat or do physical exercises for about 30 minutes before measuring.
- To the people with following condition which leads to circulatory disorders (diabetes, kidney diseases, hardening of the arteries), the values measured with this device may be lower than those measured with a gauge taking the blood pressure at the upper arm.
- Consult your doctor to make sure that measuring blood pressure at your upper arm leads to accurate measuring values.
- The measuring device is designed for the measuring of blood pressure in adults only.
- Do not use the device if your upper arm is injured.
- If the automatic inflation function of the cuff does not stop automatically it must be opened and taken off immediately!

**Note:** This device is not suitable for invasive blood pressure measuring in adults or children. It may only be used for non-invasive blood pressure monitoring and only in adults!


## PARTS DESCRIPTION

- A. Cuff jack
- B. Type-C jack
- C. LCD display
- D. SET button
- E. MEM button
- F. START/STOP button

## THE SYMBOLS ON THE LCD DISPLAY

1. Systolic blood pressure
2. Diastolic blood pressure
3. WHO classification symbol ( blood pressure level)
4. Pulse display
5. Low battery symbol
6. Date/Time display
7. Average value symbol (AVG)
8. Memory symbol (M)
9. User 1
10. User 2
11. Irregular heartbeat symbol
12. Heartbeat symbol (flashes during measurement)
13. Number of measurement memory
14. Light indicator of WHO classification

## Light indicator color / Blood pressure level related to blood pressure value

| Light indicator of WHO classification | WHO classification symbol ( blood pressure level ) | Blood pressure value   |
|---------------------------------------|----------------------------------------------------|------------------------|
| Green                                 | 1                                                  | SYS<120 & DIA<80       |
| Green                                 | 2                                                  | 120≤SYS<130 80≤DIA<85  |
| Yellow                                | 3                                                  | 130≤SYS<140 85≤DIA<90  |
| Yellow                                | 4                                                  | 140≤SYS<160 90≤DIA<100 |

|     |   |                                                         |
|-----|---|---------------------------------------------------------|
| Red | 5 | $160 \leq \text{SYS} < 180$ $100 \leq \text{DIA} < 110$ |
| Red | 6 | $\text{SYS} > 180$ $\text{DIA} > 110$                   |

## USING THIS DEVICE FOR THE FIRST TIME

### Selecting the power source

You can optionally operate the device by a AC adapter with type C connector (output : DC, 5V, 1A) or 1.5V "AA" alkaline batteries. Note : AC adapter is not included in the set,

### Using the AC adapter

- Use only the CE approved AC adapter to avoid damaging the device.
- Ensure that the AC adapter and cable are not damaged.
- Plug the adapter cable into the AC adapter jack on the right side of the device.
- Plug the adapter into your mains outlet. When the AC adapter is connected, no battery current is consumed.

### Battery Installation

- Lift the hook on the bottom of the battery cover and lift the cover off in the direction of the arrow.
- Install 3 "AA" size batteries so the + (positive) and - (negative) polarities match the polarities of the battery compartment, replace the battery cover. Make sure that the battery cover is securely in position.

### Battery replacement

#### Low Battery Indicator

- When the Low Battery Indicator appears on the display, turn the device off and remove all the batteries. Replace with 3 new batteries at the same time. Long-life alkaline batteries are recommended.
- To prevent the damage of device from leaked battery fluid, please take out of battery if the device unused in a long time (generally more than 3 months). If battery fluid should get in your eyes, immediately rinse with plenty of clean water. Contact a physician immediately.
- Dispose of the device, components and optional accessories according to applicable local regulations. Unlawful disposal may cause environmental pollution.

## SYSTEM SETTINGS

Before setting, ensure that all batteries are loaded correctly or AC adapter connected correctly.

### Time setting

Long press the SET button for 3 seconds, the year starts flashing, You can press MEM button to increase continuously 1 by 1 until 2049 and then return to the original year. After year set is done, press SET button to go to Month/Date setting.

**Note : When the set value blinks continuously for 7 seconds without any other operation, the system automatically switches off and confirms the previously set value.**

Initial Month/Date is 1/01. When the month value is flashing, press MEM button and it will increase 1 by 1, press SET button to go to date setting, and do in the same way to set the date.

After date set is done, press SET button to go to Hour/Minute setting.

Initial Hour/Minute is 01/01, When the hour display is flashing, press the MEM button, the hour will increase 1 by 1, press SET button to go to minute setting, and do in the same way to set minute.

### Pressure Unit Setting

The blood pressure value is displayed in mmHg/KPa and the default unit is mmHg. In the power off state, press and hold the START/STOP button for 10 seconds to enter the unit setting interface. Press MEM button to select mmHg or KPa. Press START/STOP button to save the selected unit automatically.

### User ID(1 or 2) setting:

Short press the SET button, the user ID shows up, You can press SET button again to choose user ID 1 or user ID 2.

Note : When the set user ID shows continuously for 12 seconds without any other operation, the system automatically switches off and confirms the user setting. Or press START/STOP button after the user ID set immediately, system automatically save the settings.

## MEASUREMENT PROCEDURE


### Before measurement:

- Avoid eating and smoking as well as all forms of exertion directly before measurement. These factors influence the measurement result. Find time to relax by sitting in an armchair in a quiet atmosphere for about ten minutes before taking a measurement.
- Remove any garment that fits closely to your wrist.
- Always measure on the same arm (normally left).

### Fitting the Cuff

- Insert the connector of rubber tube into the cuff jack. Make sure that they are properly connected.
- Wrap the cuff around your left upper arm. The rubber tube should be on the inner side of your arm extending downward to your hand. Make sure the lower edge of cuff lies approximately 2 to 3 cm above the elbow. Important! The ARTERY Mark on the lower edge of the cuff must lie over your artery which runs down the inner side of your arm.

- To fix the cuff, wrap it around your arm and secure the cuff with the hook and loop adhesive so that it lies comfortably and not too tight. There should be a little free space between your arm and the cuff. It should be able to fit 2 fingers between your arm and the cuff. Cuffs that don't fit properly result in false measurement values.
- Lay your arm on a table (palm upwards) so the cuff rests at about the same height as the heart. Make sure that the tube is not kinked.

### MEASURE YOUR BLOOD PRESSURE

The device is designed to take measurements and store the measurement values in memory for two persons using User ID 1 and User ID 2.

- Press the START/STOP button. The pump begins to inflate the cuff. The increasing pressure is continually displayed.
- After automatically reaching an individual pressure, the pump stops and the pressure slowly falls. The pressure value is displayed during the measurement.
- When the device has detected your pulse, the heart symbol in the display begins to blink.
- When the measurement has been concluded, the measured systolic and diastolic blood pressure values, as well as the pulse value will be displayed.
- According to the measurement result, the WHO symbol will also display on left corner of LCD to show at which level of WHO is your measurement result. And the WHO color indicator will also light up correspondingly.
- The measurement results are displayed until you switch the device off. If no button is pressed for 60 seconds, the device switches off automatically.

„NOTE:

The Movement Error Symbol ( Er 2 ) is displayed if you move your body during the measurement. Please remove the cuff, and wait 2-3 minutes. Reapply the cuff and take another measurement.”

### IRREGULAR HEARTBEAT DETECTOR

An irregular heartbeat rhythm is defined as a rhythm that varies by 25% higher or lower than average rhythm detected while the device is measuring the systolic and diastolic blood pressure.

The irregular heartbeat symbol displaying on LCD, which indicates that certain pulse irregularities were detected during the measurement. When irregular heartbeat is detected, the symbol will display after the measurement. If the symbol appears more frequently (e.g. several times per week on measurements performed daily) or if it suddenly appears more often than usual, we recommend the patient to seek medical advice. The instrument does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

### MEMORY

Each device stores 90 sets measurements for 2 users, totally 180 sets (User A and B). Measurements for each user are stored separately. Be certain that you are viewing the measurements for the correct user.

#### View the memory of current user ID


- Press the MEM button. The monitor will display the current measured user ID and an average value of the last 3 times measurements stored in the device. If measurements are less than 3 sets, it directly displays the first set. Each time you press the MEM button, it will display the memory value from the latest to the oldest in turn.

#### view the memory of another user ID

- Short press the SET button, the user ID shows up. You can press SET button again to choose another user ID. The selected user icon will show on LCD display. Press the START/STOP button to save the changes of user ID automatically. Do the same operations as above showed for viewing the memory of current user ID.

#### Delete memory

In memory viewing mode, long press the MEM button for 3 seconds, then it will delete all measurements for the current user.

Note: Take out the battery won't lead to a memory value missing.

### CARE AND MAINTENANCE

- Wash hands after each time measurement. If one device is used by different patients, wash hands before and after each use.
- Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- The cuff contains a sensitive air-tight bubble. Handle this cuff carefully and avoid all types of stress through twisting or buckling.
- Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. The cuff with bladder must not be washed in a dishwasher, clothes washer, or submerged in water.
- Handle the tube carefully. Do not pull on it. Do not allow the tubing to kink and keep it away from sharp edges.
- Do not drop the monitor or treat it roughly in any way. Avoid strong vibrations.
- Never open the monitor! This invalidates the manufacturer's warranty.
- Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

### ERROR INDICATES

| SYMBOL | CAUSE                                                                                             | CORRECTION                                                                                   |
|--------|---------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|
| Er 1   | No pulse detected                                                                                 | Wear the cuff correctly to keep the air tube smooth and do not remove the air tube.          |
| Er 2   | Too much interference                                                                             | Be quiet during the measurement. Do not talk nor move; Do not squeeze the cuff and air tube. |
| Er 6   | Air tube blockage                                                                                 | Keep the air tube smooth and do not remove the air tube.                                     |
| Er 7   | Sensor abnormal                                                                                   | Sensor defective, return to factory to repair.                                               |
| ErH    | When the inflation pressure is higher than 295mmHg, it will display an error and quickly deflate. | Be quiet during the measurement. Do not talk nor move; Do not squeeze the cuff and air tube. |
| ErU    | Air leak causing inflatable                                                                       | Make sure the cuff and air tube are properly connected to the main unit.                     |

### TECHNICAL SPECIFICATIONS

LCD Display size: 52 \* 70mm

Product size: 126 (W) \* 130 (L) \* 50 (H) mm

Operating temperature: 5° to 40°

Operating humidity: 15% to 80% RH

Storage and shipping temperature: -25° to 55°

Storage and shipping humidity: 10% - 93%RH

Atmospheric pressure range: 86kPa ~ 106kPa

Measuring method: Oscillometric

Pressure sensor: Resistive

Measuring range: DIA: 40-190mmHg; SYS: 60-250mmHg

Pulse: 40 to 180 per minute

Cuff pressure display range: <290mmHg

Memory: Automatically stores 90 measurements for 2 users (total 180 measurements)

Measuring resolution: 1 mmHg

Accuracy: Pressure within  $\pm 3$  mmHg / pulse  $\pm 5$  % of the reading

Power source A: 4\*AA batteries, 1.5 V (not included)

Power source B: AC adapter (not included): Input: 220-240V, 50/60 Hz. Output: 5V DC 1A

USB cable: Input port-type C output port-USB 50cm

Cuff: 22-32cm

Automatically power off: 60 seconds without operation

## WARRANTY

1. This warranty is issued by EDC Poterek Sp. Jawna with address:  
ul. Poznańska 129/133, 05-850 Ożarów Mazowiecki, Phone: 022-7213571 w.10 – called hereinafter EDC Poterek Sp. Jawna.
2. This warranty confers upon you specific rights and remedies against EDC Poterek Sp. Jawna. Please read this document carefully, by making a claim under this warranty, you will be deemed to have understood and accept its terms.
3. EDC Poterek Sp. Jawna warrants that any EDC Poterek Sp. Jawna product together with attached accessories is free from defects in materials and workmanship under conditions of its normal use.
4. If you make a claim under this warranty, EDC Poterek Sp. Jawna will, at its option, repair or replace any defective product into new one by using new or use factory-refurbished parts. In case of using factory-refurbished parts they will be equal to new products.
5. Defective product should be returned together with receipt, commercial invoice or any other proof of purchase to the place of purchase. It should include all accessories attached in original packing. EDC Poterek Sp. Jawna has right to reject warranty if any accessories will be missed.
6. EDC Poterek Sp. Jawna covers the costs of delivery of the goods complained about to the place indicated in the guarantee.
7. EDC Poterek Sp. Jawna is committed to respond to the complaint within 14 calendar days of its notification in the point of sale. EDC Poterek Sp. Jawna will repair the defective product in the shortest possible time. If will be need to import from abroad of any parts which are necessary to repair this time may be extended.
8. If in the performance of its duties EDC Poterek Sp. Jawna provided the entitled person with a warranty instead of a defective item free from defects or made significant repairs to the item covered by the warranty, the warranty period runs anew from the moment of delivery of the item free from defects or returning the repaired item. If EDC Poterek Sp. Jawna mentioned a part of things, the above provision shall apply accordingly to the part mentioned. In other cases, the warranty period is extended by the time during which, as a result of a defect in the item covered by the warranty, the holder of the warranty could not use it.
9. In case of replacing defective product onto new one, EDC Poterek Sp. Jawna has right to replace it with newer model with the same technical parameters if the old one model was replaced by new one during warranty period.
10. If you make claim under this warranty, you are deemed to have agreed that repair and (upon availability) replacement, as applicable, under the warranty services described herein, is your sole and exclusive remedy vis-a-vis EDC Poterek Sp. Jawna.
11. The warranty is valid in Poland.

## EXCLUSIONS:

### This warranty does not extend to:

1. Any product not manufactured by EDC Poterek Sp. Jawna
2. Any product that has been damaged or rendered defective as a result of:
  - use it in other way than normal intended use
  - service or repair by anyone other than EDC Poterek Sp. Jawna service
  - improper store, transportation or packing
  - improper installation of product or software
  - mechanical, physical, chemical etc. damages
3. In no event shall EDC Poterek Sp. Jawna be liable for any damages, loss of business, profits, data or loss based on use EDC Poterek Sp. Jawna products with other devices.
4. Products with removed warranty and other labels and products spent or have hard use marks.

All products are covered by EDC Poterek Sp. Jawna 2-year manufacturer's warranty while product is purchased by the consumer (the consumer is a person who buys product for the purpose of not professional or commercial activity). In the case of purchase for the company (for purposes relating to the activities of professional or commercial) warranty is 1 year from date of purchase. EDC Poterek Sp. Jawna products have all the necessary safety certificates required in the EU.

According to Directive 2012/19/EU, disposal of waste electrical and electronic equipment together with household waste is prohibited. Their components need to be recycled or disposed separately. Otherwise, they can cause the release of hazardous and noxious substances which pose a threat to the environment and human health. The law requires each consumer to return unnecessary or waste electrical and electronic equipment to the producers, sales agents or waste treatment facilities. Detailed conditions are governed by national legislation. The above symbol placed on the product, in the user manual or on the packaging refers to this requirement. Thanks to proper segregation and compliance with the provisions relating to waste disposal each user contributes significantly to environmental protection.

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**EC Certificate**  
**Directive 93/42/EEC Annex V**  
**Production Quality Assurance**  
**Medical Devices**

**Registration No.:** DD 60148122 0001

**Report No.:** 17054607 002

**Manufacturer:** DONGGUAN RICHJINN MEDICAL  
TECHNOLOGY CO., LTD.  
Room 233, Building #1, No.12 Industrial  
West Road, Songshan Lake Park  
523808 Dongguan, Guangdong,  
P.R. China

**Products:** Electronic Thermometers and Blood Pressure Monitors

**Expiry Date:** 2024-05-26

The Notified Body hereby declares that the requirements of Annex V of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex V, section 4 of the aforementioned directive. For placing on the market of class IIb and class III devices covered by this certificate an EC type-examination certificate according to Annex III is required.

**Effective Date:** 2020-07-27

**Notified Body**

**Date:** 2020-07-27

Jason Pan

**TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg**  
TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC  
concerning medical devices with the identification number 0197.

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*GUANGDONG RICHJINN MEDICAL TECHNOLOGY CO., LTD.  
201, Building 38, Xi Tian Second Industrial  
Park, GongMing Street, Guangming District,  
Shenzhen, 518107, Guangdong,  
P.R. China*

Contact

Tel. +49 911 655-5225  
Mail: [medical\\_products@de.tuv.com](mailto:medical_products@de.tuv.com)  
Date May 11, 2024

### **Notified Body Confirmation Letter**

Reference. : RICHJ\_PLA\_2024-03-14; order#10924253

To whom it may concern,

**Confirmation of the status of a formal application, written agreement, and appropriate surveillance in the framework of Regulation EU 2023/607 amending Regulations (EU) 2017/745 and (EU) 2017/746 as regards the transitional provisions for certain medical devices and in vitro diagnostic medical devices.**

This letter confirms that **TÜV Rheinland LGA Products GmbH**, a Notified Body (NB) designated against Regulation (EU) 2017/745 (MDR) and identified by the number **0197** on NANDO, has received a formal application in accordance with Section 4.3, first subparagraph of Annex VII of MDR and has signed a written agreement in accordance with Section 4.3, second subparagraph of Annex VII of MDR with the following manufacturer:

GUANGDONG RICHJINN MEDICAL TECHNOLOGY CO., LTD.  
201, Building 38, Xi Tian Second Industrial  
Park, GongMing Street, Guangming District,  
Shenzhen, 518107, Guangdong,  
P.R. China  
SRN Number (if available): CN-MF-000040778

The devices covered by the formal application and the written agreement mentioned above are identified in the tables below. Table 1 identifies the devices for which an MDR application has been received, written agreement concluded and for which the NB is also responsible for appropriate surveillance under the applicable Directive. Table 2 identifies the devices for which an MDR application has been received and a written agreement concluded, but the NB has not yet taken the responsibility for appropriate surveillance of the corresponding devices under the applicable Directive.

In the case of devices covered by certificates issued under Directive 90/385/EEC (AIMDD) or Directive 93/42/EEC (MDD) that expired after May 26, 2021 but before March 20, 2023 without having been withdrawn, this letter also confirms that the manufacturer either signed the written agreement under MDR by the date of MDD/AIMDD certificate expiry; or provided evidence that a competent authority of a Member State had granted a derogation or exemption from the applicable conformity assessment procedure in accordance with Article 59(1) of MDR or Article 97(1) of the MDR respectively, by March 20, 2023 for the relevant devices.

TÜV Rheinland  
LGA Products GmbH

Am Grauen Stein  
51105 Köln  
Germany

Headquarter

Tillystraße 2  
90431 Nuremberg

Phone. +49 911 655 5225  
Fax +49 911 655 5226  
[service@de.tuv.com](mailto:service@de.tuv.com)  
[www.tuv.com/safety](http://www.tuv.com/safety)

Board of Management

Dipl.-Ing.  
Thomas Weigand, Spokesman

Dipl.-Kfm.  
Dr. Jörg Schlösser

Nuremberg HRB 26013  
VAT No.: DE 811835490

Chairman of the  
Supervisory Board

Dr.-Ing. Michael Fühl

The transition timelines that apply to the devices covered by this letter, subject to the manufacturer's continued compliance to the other conditions specified in Article 120.3c of MDR (as amended by (EU) 2023/607), are shown below:

- May 26, 2026 for Class III custom-made implantable devices
- December 31, 2027 for Class III devices and Class IIb implantable devices excluding Well-established technologies (WET - sutures, staples, dental fillings, dental braces, tooth crowns, screws, wedges, plates, wires, pins, clips and connectors)
- December 31, 2028 for other Class IIb devices, Class IIa, Class I devices placed on the market in sterile condition or have a measuring function
- December 31, 2028 for devices not requiring the involvement of a notified body under MDD but requiring it under MDR (e.g., class I devices that qualify as re-usable surgical instruments)

On behalf of the Notified Body

Samuel Qin

Certification body

**Table 1: Devices covered by this letter and for which the NB is also responsible for appropriate surveillance of the corresponding devices under the applicable Directive:**

| Device name or Basic UDI-DI (under MDR application)                                                                                      | MDR Device classification (as proposed by the manufacturer and verified at the pre-application stage) | If the MDR device is a substitute device, identification of the corresponding MDD/AIMDD device | MDD/AIMDD Certificate Reference(s) of the devices under MDR application, and the NB Identification |
|------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|
| <b>Digital Blood Pressure Monitor (Arm Type)</b><br><b>Model:</b><br>RJA-001, RJA-002<br><br><b>Basic UDI-DI:</b><br>697545553RJA001MM   | Class IIa                                                                                             | N/A                                                                                            | Certificate # DD<br>60148122 0001<br><br>NB# 0197                                                  |
| <b>Digital Blood Pressure Monitor (Wrist Type)</b><br><b>Model:</b><br>RJW-001, RJW-002<br><br><b>Basic UDI-DI:</b><br>697545553RJW001SF | Class IIa                                                                                             | N/A                                                                                            | Certificate # DD<br>60148122 0001<br><br>NB# 0197                                                  |
| <b>Digital Thermometer</b><br><b>Model:</b><br>RJT-001, RJT-002<br><br><b>Basic UDI-DI:</b><br>697545553RJT001RS                         | Class IIa                                                                                             | N/A                                                                                            | Certificate # DD<br>60148122 0001<br><br>NB# 0197                                                  |

**Table 2: Devices covered by this letter and for which the NB is NOT responsible for appropriate surveillance of the corresponding devices under the applicable Directive:**

| Device name or Basic UDI-DI (under MDR application) | MDR Device classification (as proposed by the manufacturer and verified at the pre-application stage) | If the MDR device is a substitute device, identification of the corresponding MDD/AIMDD device | MDD/AIMDD Certificate Reference(s) of the devices under MDR application, and the NB Identification |
|-----------------------------------------------------|-------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|
| N/A                                                 | N/A                                                                                                   | N/A                                                                                            | N/A                                                                                                |

#### **Confirmation Letter Revision History**

| Date       | NB internal reference traceable to each version of the letter | Action        |
|------------|---------------------------------------------------------------|---------------|
| 2024-05-11 | RICHJ_CL607_2024-05-11                                        | Initial issue |

## DECLARATION OF CONFORMITY

We, the manufacturer:

**Guangdong Richjinn Medical Technology Co.,Ltd.**

201, Building 38, Xi Tian Second Industrial Park, GongMing Street, Guangming District, Shenzhen, 518107 Guangdong, P.R.China

**We confirm that this product has CE.**

Herewith, Declares that the product:

**Product description: Blood Pressure Monitor**

**Model Number: RJA-002 (ECB004)**

The products is compatible with the following norms/standards:

**Directive 93/42/EEC on Medical Devices (MDD), Annex V**

**Issued by TUV Rheinland LGA Products GmbH, Notified body: 0197, EC Certificate No.: DD 60148122 0001**

**This certificate is valid from 27 July 2020 until 26 May 2024.**

Manufacturers/Authorized representative

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2023-09 SHENZHEN

## MODEL IDENTITY DECLARATION

We, the manufacturer:

**Guangdong Richjinn Medical Technology Co.,Ltd.**

201, Building 38, Xi Tian Second Industrial Park, GongMing Street, Guangming District, Shenzhen, 518107 Guangdong, P.R.China

Herewith, Declares that the product:

Product description: **Blood Pressure Monitor**

Model Number: **RJA-002**

Is technically identical with the following product imported from our company by EDC Poterek Sp. Jawna company:

Product description: **Blood Pressure Monitor**

Model Number: **ECB004**

Trademark: **ESPERANZA**

Manufacturers/Authorized representative

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2023-09 SHENZHEN