FIZJOLOVE

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Importer:
Mydecor sp. z o.o.
ul. Ostródzka 112a
03-289 Warszawa
NIP: 5242759978

The acupressure mat with pillow has been designed for daily relaxation, body recovery, and support of natural relaxation processes. Regular use may help reduce muscle tension, improve well-being, and enhance the quality of rest.

1. Package Contents

- Acupressure mat
- Acupressure pillow
- User manual

2. Principle of Operation

The mat and pillow are equipped with special stimulating elements (rosettes with spikes) that act on the body's acupressure points. Gentle pressure:

- stimulates blood circulation,
- supports muscle relaxation,
- promotes a feeling of deep relaxation.

3. Preparation for Use

- Lay the mat on a flat, stable surface (floor, exercise mat, bed).
- Place the pillow in a position suited to the selected body area (most commonly under the neck or lower back).
- For initial sessions, it is recommended to place a thin layer of clothing or a towel between the body and the mat.

4. Methods of Use

Lying on your back

- Slowly lie down on the mat.
- Place the pillow under your neck.
- Breathe calmly and deeply.

Lying on your stomach

- Recommended for more advanced users.
- Helps stimulate the abdominal area.

Sitting

- Place the mat on the backrest of a chair or armchair.
- An ideal solution for office work.

Feet

- Stand or sit with your feet resting on the mat.
- Helps relax the entire body.

5. Duration of Use

- Beginners: 5-10 minutes
- Intermediate users: 10-20 minutes
- Advanced users: up to 30-40 minutes
- Recommended frequency: 1-2 times per day

6. Safety and Contraindications

Do not use the product if you:

- have blood clotting disorders,
- suffer from skin conditions at the contact area,
- are pregnant (without consulting a doctor),
- have open wounds or inflammation,
- experience severe discomfort or pain.

If in doubt, consult a doctor.

7. Cleaning and Care

- Clean by hand with a damp cloth and mild detergent.
- Do not machine wash.
- Do not tumble dry.
- Store in a dry place, away from sunlight.

8. Storage

- Roll the mat loosely, without excessive bending.
- Do not place heavy objects on the product.

We wish you a pleasant and relaxing experience!