
FIZJOLOVE

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The basic-type exercise bike is a compact rehabilitation device designed for exercising both the lower and upper limbs. It can be used at home, in rehabilitation centers, and in physiotherapy clinics.

The device helps improve circulation, joint mobility, and overall physical fitness.

2. SAFETY RULES


WARNING: Before using the device, carefully read this manual.

- ▲ The device is intended solely for rehabilitation and recreational exercises performed in a seated position.
- ▲ Do not stand on the bike or use it while standing.
- ▲ Always use the safety straps on the pedals to prevent your feet or hands from slipping.
- ▲ Perform exercises on a stable, flat surface.
- ▲ Do not use the device if any parts are damaged or loose.
- ▲ Children may use the device only under adult supervision.
- ▲ Stop exercising immediately if you experience pain, dizziness, or discomfort.
- ▲ Maximum allowable user weight: 100 kg.

3. PACKAGE CONTENTS

The set includes:

- Basic exercise bike frame
- 2 stabilizers
- 2 pedals with straps
- LCD display
- Screws, nuts, and wrench
- User manual

4. DEVICE ASSEMBLY

Stabilizer Assembly

- Position the stabilizers at the lower part of the main frame.
- Secure them using the supplied screws.
- Tighten all connections with the included wrench, ensuring they are firmly fastened.

Pedal Assembly

- Identify the pedals marked L (left) and R (right).
- Install the L pedal into the left crank arm by turning it counterclockwise.
- Install the R pedal into the right crank arm by turning it clockwise.
- Ensure both pedals are securely tightened to prevent loosening during exercise.

LCD Display Installation

- Attach the display bracket to the frame of the unit.
- Connect the display cable to the corresponding socket in the frame.
- Ensure the display is securely mounted and powers on correctly.

Device Adjustment

- Place the exercise bike on a flat and stable surface.
- Check that all components have been installed correctly.
- Set the resistance adjustment knob to the initial position, then adjust the resistance level according to your needs.

Pre-Use Inspection

- Check that all screws and connections are properly tightened.
- Ensure the pedals rotate smoothly.
- Turn on the LCD display and verify that the displayed information is correct.

Completed

The exercise bike is correctly assembled and ready for use.

Remember to regularly inspect the technical condition of the device to ensure safe and effective training.

5. DISPLAY FUNCTIONS

To start the display, press the red button or begin pedaling.

The unit will automatically shut off after approximately 4 minutes of inactivity.

During training, all functions activate automatically and are deactivated when pedaling stops.

The display mode can be changed by pressing the red button; the indicator (arrow) shows the currently selected function.

Available Display Modes:

- **CNT (Count)** - Displays the total number of revolutions performed.
- **TIME** - Shows the duration of the current workout.
- **CAL (Calories)** - Displays an approximate number of calories burned.
- **RPM (Revolutions Per Minute)** - Indicates the number of revolutions per minute.
- **SCAN** - Automatic mode in which the display cycles through the functions in the following order:
TIME → CNT → RPM → CAL

6. EXERCISE AND TRAINING RECOMMENDATIONS

Leg Exercises

Starting Position:

Sit upright with a straight posture. Place your feet on the pedals and rest your hands on your thighs or on the armrests.

Exercise Execution:

Begin pedaling forward with a smooth, continuous motion.

Direction Change:

To engage different muscle groups, periodically pedal backward during the exercise.

www.fizjoloove.pl

8. ENVIRONMENTAL PROTECTION

Compliance with regulations

The Rotor rehabilitation bike has been designed and manufactured with respect for environmental protection principles and in accordance with applicable European Union regulations, including the WEEE Directive (2012/19/EU) on waste electrical and electronic equipment.

The crossed-out trash bin symbol on the product indicates that it must not be disposed of with household waste.

Disposal and recycling

- After the end of its use, the bike should be taken to an electrical and electronic equipment collection point or an authorized service center.
- Do not dispose of the device in regular trash bins.
- Sorting and recycling materials help reduce environmental impact, save energy, and decrease the amount of waste in landfills.

Batteries and electrical components

- If your bike model contains batteries or accumulators:
 - Remove them according to local regulations for hazardous waste.
 - Take batteries to designated battery and accumulator collection points.
- Improper handling of batteries can pose risks to the environment and human health.

Packaging

The bike's packaging is made from recyclable materials (cardboard, plastic, paper).

After unpacking, take the packaging to recycling points according to local waste sorting regulations.

9. IMPORTER

The importer is Mydecor sp. z o.o., headquartered at: Ostródzka 112a, 03-289 Warszawa NIP:5242759978

The product complies with the requirements of applicable EU directives and safety standards.

Arm Exercises

Starting Position:

Place the exercise bike on a stable, flat surface, such as a table, at an appropriate height.

Grip:

Hold the pedals with your hands, keeping your wrists in a neutral, straight position.

Movement:

Perform circular hand motions, maintaining a steady tempo and smooth exercise flow.

Training Recommendations

Consistency:

Exercise regularly, starting with 10–15 minutes per day, then gradually increase the duration and intensity of your workouts.

Breathing:

Maintain calm and even breathing during exercises to ensure proper oxygenation of the body.

Hydration:

Drink water before starting your workout and after finishing to stay properly hydrated.

7 . MAINTENANCE AND CLEANING

- After each workout, clean the device with a dry, soft cloth.
- Do not use harsh chemicals.
- Regularly check the condition of the pedals and straps.
- Store in a dry place.