FIZJOLOVE

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Szanowny Kliencie,

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CUPPING THERAPY SET

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Dear Customer,

Thank you for choosing a product from the **Fizjolove** brand. We are delighted that you have placed your trust in us and we believe that your purchase will meet your expectations. Our mission is to support health and well-being, which is why we carefully select our partners and ensure the quality of the products we offer. As a responsible and trustworthy company, we care about our customers and the community in which we operate.

Before using the product for the first time, please read the user manual carefully – it contains important information and practical tips that will help you use the device correctly and safely.

With kind regards,
The **Fizjolove** Team

1. WARNINGS

This product has been designed with safety in mind; however, before using it, be sure to read this manual carefully.

Failure to follow the recommendations may cause skin irritation or damage. The manufacturer is not responsible for the consequences of improper use.

- Keep the cups out of reach of children.
- Do not apply the cups on wounds, swelling, birthmarks, varicose veins, or areas where the skin is damaged or highly sensitive.
- Do not place the cups near the heart, neck, chest, or the abdomen of pregnant women.
- Avoid using the cups immediately after a meal.
- Always monitor the suction strength – excessive suction may cause bruising or broken capillaries.

2. CONTRAINDICATIONS

Cupping is not recommended in the following cases:

- Skin conditions (e.g., fungal infections, eczema, wounds)
- Blood clotting disorders, varicose veins, high blood pressure, heart disease, or cancer
- Fever weakened health or acute infections

If you have any doubts about performing the treatment, consult a doctor who can assess whether cupping therapy is safe for you.

3. GENERAL INFORMATION

Cupping massage improves blood circulation, enhances skin elasticity, and supports tissue regeneration. Regular use may help reduce muscle tension and improve the appearance of the skin.

For best results, it is recommended to perform the massage 2–3 times per week, and after completing a series, allow at least 24 hours before the next session.

4. HOW TO USE

1. Preparation:

Clean and dry the skin, then apply oil or cream to make it easier to move the cup.

2. Applying the cup:

Choose the appropriate cup size. Lift the valve to allow airflow.

Then place the cup on the skin and attach the pump. Pull the pump handle until the cup gently suctions the skin.

3. Adjustment:

Check that the cup is securely attached. If the suction is too strong, slightly lift the valve to release some air.


4. Removal:

To remove the cup, hold it with your hand and gently lift the valve to let air in.

Note: During the treatment, you may feel warmth and mild skin tension – this is a normal reaction. If you experience pain or severe redness, stop the massage immediately.


5. MAINTENANCE AND STORAGE

- After each use, wash the cups in warm water with a mild detergent, then dry them thoroughly.
- Store the set in a dry, clean place, away from heat sources and direct sunlight.
- Do not sterilize the cups in boiling water or in a microwave.

1. Choosing the cup


Select a cup of the appropriate size. Lift the valve to allow air to flow freely.

2. Suction


Attach the cup to the pump nozzle and place it on the chosen area of the body.

Gently pull the pump handle to create suction and attach the cup. Do not pump too strongly to avoid excessive pressure.

3. Using the tube

If you need more flexibility, connect the pump to the cup using the flexible tube.

4. Removal

To remove the cup, hold it with your hand and gently lift the valve to release the air. Once the pressure is equalized, the cup will detach from the skin on its own.