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Device and controls

- 1 On/off key
- 2 Temperature level key
- 3 Display and indicator light
- 4 Operating time key

Explanation of symbols

 Do not use the heated cushion when it is folded!

 Do not puncture the heated cushion

 Not suitable for children under 3 years!

 Only use the heated cushion indoors!

 The heat pad can be cleaned in a gentle wash cycle at max. 40°C!

 Do not bleach!

 Do not dry the heated cushion in a tumble drier!

 Do not iron the heated cushion!

 Do not dry clean!

 **IMPORTANT!**
Non-observance of these instructions can result in serious injury or damage to the device.

 **WARNING**
These warning notes must be observed to prevent any injury to the user.

 **CAUTION**
These notes must be observed to prevent any damage to the device..

 **NOTE**
These notes provide useful additional information on the installation or operation.

 Protection category II

LOT LOT number

 Manufacturer

 Recycling symbols/codes: These are used to provide information about the material and its proper use and recycling.

GB Safety instructions

 **Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use. Should you give this device to another person, it is vital that you also pass on these instructions for use.**

- Carefully check the heat pad before each use, if you notice wear or damage.
- Do not use if you notice wear, damage or signs of improper use on the heat pad, the switch or the cables or if the device is not working.
- Before connecting the device to your power supply, please ensure that the supply voltage stated on the rating plate is compatible with your mains supply.
- Do not fold or crease the heat pad during use.
- Neither pierce the pad with, nor attach safety pins or other spiky or sharp objects to the pad.
- Do not use the heat pad on children or persons who are disabled, asleep or insensitive to heat (persons who cannot react on overheating).
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision and instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children shall not play with the appliance.
- Cleaning and user maintenance shall not be made by children without supervision.
- Children below the age of 3 years are not allowed to use this device, because they are not able to react on overheating.
- The heat pad must not be used by young children over the age of 3 years, unless the control unit has been adjusted accordingly by its parents or legal guardians and unless the child has been well instructed how to use the control unit properly.
- The appliance can be used by children older than 3 years and younger than 8 years under supervision, but the control unit must always be set to the lowest level.
- The unit's cable and control unit can cause tangling, strangulation, tripping and other hazards if improperly arranged. The user must ensure that excess cable fixings and electrical wiring are safely routed.
- Never use the pad unsupervised.
- The device is only for private use and is not to be used in hospitals.
- Do not fall asleep while the heat pad is switched on.
- Move the switch to position 1 for any continuous use.
- Prolonged use of the heated cushion at a high setting could lead to skin burns.
- Do not cover the pad with any other cushions.
- The control unit should not be covered or placed on or below the pad while the unit is in operation.
- Do not sit on the heat pad. Instead, rest the pad on or against the part of the body to be treated.
- Never touch a device that has fallen into water. Unplug the unit from the mains outlet immediately.
- Keep the mains lead away from hot surfaces.
- Never carry, pull or turn the heat pad by the mains lead and do not let the lead get tangled.
- Do not use the heat pad when wet and only use it in dry environments (not in the bathroom or similar environments).
- The controls and leads must not be exposed to any kind of moisture.
- If a fault occurs, do not attempt to repair the unit yourself. Repairs must only be carried out by an authorized specialist dealer or other suitably qualified personnel.
- The heat pad may only be operated with the control unit (HP 630) supplied.

- If the power supply cable is damaged, it may only be replaced by **medisana**, an authorised dealer or by qualified personnel in order to avoid hazards.
- A damaged power supply cable must be replaced with a power cable of the same type.
- When storing the heat pad, allow it to cool down before laying out flat in a clean and dry place.
- Do not store the device with any objects covering or resting on top of it in order to avoid sharp crinkles.
- Should you have any concerns about health issues, consult your doctor before using the heat pad.
- Do not use the heat pad on parts of the body that are swollen, injured or inflamed.
- Should you experience persistent pains in the muscles and joints, please inform your doctor. Persistent pains can be a symptom of serious illness.
- If you wear a pacemaker, please be aware that the electromagnetic fields emanating from this electric product may, under certain circumstances, interfere with the functions of your pacemaker. Therefore, ask your doctor and the manufacturer of the pacemaker for advice before using this product.
- If you experience any pain or discomfort while using the unit, stop using it immediately.

Scope of supply

Please check first of all that the device is complete and is not damaged in any way. In case of doubt, do not use the it and contact your dealer or your service centre.

The following parts are included:

- 1 **medisana** heat pad for shoulders and back **HP 630** with removable control unit
- 1 Instruction manual

If you notice any transport damage during unpacking, please contact your dealer without delay.

 **WARNING**
Please ensure that the polythene packing is kept away from the reach of children!
Risk of suffocation!

Operating

Always make sure that the heat pad is not bent or folded when in use.

- Establish a power connection. The indicator light on display ③ lights up.
- Press the on/off key ① to switch on the appliance. The control light goes out and a „1“ (for temperature level 1) appears in the display ③.
- You can set one of 4 temperature levels by repeatedly pressing the temperature level key ②, with level 4 being the one with the highest heating.
- To preset an operating time, use the operating time key ④. By repeatedly pressing it, you can preset a time (30 / 60 / 90 minutes). When you press the operating time key ④, the currently set time only appears for a few seconds on the display ③ together with the indicator light that lights up. Then the temperature level (without indicator light) is displayed again.
- After the preset operating time has elapsed, the appliance switches off automatically. However, you can also switch it off manually at any time by pressing the on/off key ①.

Unplug the unit from the mains outlet once you have finished using it.

The particular shape is especially designed for use on shoulders, neck and back. The heat pad features an electric temperature control that regulates the temperature according to the chosen setting. The two snap buttons on the front part and the practical belt keep the heat pad for shoulders and back **HP 630** in its position.

Cleaning and maintenance

- Before cleaning the pad, unplug it from the mains outlet and let it cool down for at least ten minutes.
- The heat pad is equipped with a removable control unit (HP 630). Unplug the control unit from the heat pad by pulling out the connection cable from the plug on the heat pad.
- The heat pad can be cleaned dry (with a soft brush) or as described in the instructions on the product label.
- Never use strong detergents or hard brushes.
- To dry, lay the pad flat on a moisture-absorbent towel or mat.
- Only use the heat pad again once it is completely dry.
- Reconnect the control unit to the heat pad in the correct manner (the arrow markings are pointing towards each other).
- You should untangle the cable from time to time.
- Store the unit in a clean and dry place, laid out flat with nothing covering or resting on top of it.

Disposal

 This product must not be disposed together with the domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Consult your municipal authority or your dealer for information about disposal.

Technical data

Name and model:	medisana Heat pad for shoulders and back HP 630
Power supply:	230V~, 50 Hz
Heat output:	100 watts
Automatic switch-off:	after approx. 30/60/90 minutes
Dimensions approx.:	55 x 65 cm
Weight approx.:	0.71 kg
Operating conditions:	Only use in dry rooms according to this user manual
Storage conditions:	Laid out flat, clean and dry
Article number:	61158
EAN number:	40 15588 61158 2

In accordance with our policy of continual product improvement, we reserve the right for technical and design changes without notice.

The current version of this instruction manual can be found under www.medisana.com

Warranty and repair terms

Your statutory warranty rights are not limited by our warranty set out below. Please contact your dealer or the service centre in case of a claim under the warranty. If you have to return the unit, please enclose a copy of your receipt and state what the defect is.

The following warranty terms apply:

- The warranty period for **medisana** products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
- Defects in material or workmanship will be removed free of charge within the warranty period.
- Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
- The following is excluded under the warranty:
 - All damage which has arisen due to improper treatment, e.g. non-observance of the user instructions.
 - All damage which is due to repairs or tampering by the customer or unauthorised third parties.
 - Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
 - Accessories which are subject to normal wear and tear.
- Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

 **medisana GmbH**
Carl-Schurz-Str. 2
41460 NEUSS
GERMANY

The service centre address is shown on the attached leaflet.

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