medivon.

User manual

Cosy Air Pro

PL | EN | DE | FR | ES | PT | IT | CZ | SK
SL | EE | LT | LV | HU | RO | RU | NB | FI
EL | HR | SV | DA | NL | TR | BG | AR

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1. Configuration and usage

- Attach the mat to the chair with the elastic straps
- Connect Mat connector (1) with the power supply connector (2)
- Plug the home power adapter into an electrical outlet
- Turn on the device with the remote control (instructions later)
- Place the RC in the bag on the side of the massage mat

Note: The removable flap (3) or cover for the neck massage section allows you to adjust the intensity between a more gentle or more intense massage. If the massage is too strong, place a blanket or towel between you and the mat (4)

If you want to massage your neck, shoulder, lower back, but the shiatsu massage heads cannot reach these areas, you can adjust the height of the massage by placing a pillow on the sitting area (5).

To avoid overloading the muscles, we recommend that you do not exceed 15 minutes for one massage session. Do not massage yourself more than 3 times during the day. (in the full range of device operation - 15 minutes)

2. Care and cleaning

- Do not machine wash or submerge the product or adapter in water.
- Always turn off the power before cleaning and remove the adapter from the power supply with clean, dry hands.
- Wipe the device with a slightly damp cloth and dry it thoroughly.
- Do not wipe the device with preparations containing alcohol or harsh chemicals such as thinner or bleach.
- The neck cover is removable and washable.
- Do not stretch the elastic straps excessively on the back.

- Store away from heat, fire or direct sunlight.
- Do not put excessive pressure on the device and the massage heads. Applying a direct force in excess of 49 kg may damage the massager and void the warranty.
- Store the massager in a dry place inaccessible to children.
- Avoid contact with sharp edges or sharp objects that could cut or puncture the material of the device.

3. Remote controller description

- (6) - press to turn on the massager vibration, you can switch between 3 levels of intensity.
- (7) - press to change the 2D / 3D massage mode. In the 3D option, it additionally takes into account the exact area of the spine.
- (8) - press to keep the shiatsu function in one area. Press the arrow buttons to move the shiatsu heads up or down.
- (9) - press to set the intensity of air massage on low, medium or high level.
- (10) - press to enable or disable rear air massage.
- (11) - press to designate the area with shiatsu massage.
- (12) - press to turn the device on or off.
- (13) - press to turn the heating on or off. The heating cannot be turned on if the back shiatsu function is not turned on.
- (14) - press to enable or disable seat air massage.
- (15) - press to turn on the roll massage function. Repeated pressing changes modes.
- (16) - press to start shiatsu massage on the neck. The arrows change the position of the massage heads in the mat.
- (17) - a quick demonstration of all the functions of the massage mat

4. Safety instructions

- The device should never be left unattended when connected to the power supply. Before assembling or disassembling, disconnect the device from the power supply.
- Do not use the device under a cover where the ventilation opening may be blocked. The device may overheat.
- This device requires close supervision by adults when it is used by children or people with disabilities.
- Do not use with other accessories not recommended by the manufacturer.
- Never use the appliance if it has a damaged cord or plug, if it does not work properly, if it has been dropped, damaged or dropped into water. In such cases, contact the service immediately.
- Do not carry the device with the power cord or use it as a handle.
- Never use the device with blocked air openings.
- Never insert any objects into any openings.
- Do not use outdoors.
- When using, do not wear jewelry or loose clothing.
- Longer hair must be tied back before use so that it does not come into

contact with the device.

- Do not use the warm-up function in places where you have skin disorders, open wounds, swelling or inflammation.
- If the plugs or cables are damaged, they must be replaced by the customer service
- In case of failure, disconnect the device from the mains immediately.
- Do not use in places if you have skin abnormalities, open wounds, swelling or inflammation.
- Improper use or incorrect use excludes all liability for damages.
- Do not use this device while driving or sleeping.
- We have equipped the device with overheating protection for safe use.
- Do not use the device during pregnancy.
- Do not use if you have pacemakers or are under treatment for diseases or other ailments.

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Used electrical and electronic equipment and batteries and accumulators must not be mixed with other waste. This is indicated by the crossed-out bin symbol on the product packaging. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution. Proper handling of used electrical and electronic equipment and batteries and accumulators in the household contributes positively to the reuse and recovery, including recycling, of used equipment. Used electrical and electronic equipment and batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the used equipment is of the same type and has performed the same functions as the equipment sold. A distributor, when delivering household equipment to a purchaser, is obliged to take back used household equipment free of charge at the place of delivery of such equipment, provided that the used equipment is of the same type and has performed the same functions as the equipment delivered. A distributor running a retail outlet with a sales area within the meaning of art. 2 item 19 of the Act of 27 March 2003 on spatial planning and development of at least 400 m² devoted to the sale of household equipment, is obliged to accept, free of charge, in that outlet or in its immediate vicinity, used household equipment, none of the external dimensions of which exceeds 25 cm, without the need to purchase new household equipment. The product was introduced to the market after 13 August 2005.

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Short Declaration of Conformity

The manufacturer declares under its sole responsibility that the Cosy Air Pro product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility, Directive 2014/35/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits, and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU as regards the list of restricted substances.

The full text of the Declaration of Conformity is available at the following website: www.medivon.pl/deklaracje

Importer: Armare S.A.
ul Grzybowska 87,
00-844 Warszawa, Poland
Made in China
Marka: Medivon
E-mail: serwis@medivon.pl
www.medivon.pl
SN:01265908275917137

Producent: XIAMEN COMFIER TECHNOLOGY CO.,LTD
2-5/F, No.65 Sunban South Road, Jimei District, Xiamen, China.
E-mail: sales@comfier.com

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Shortened warranty terms

This Guarantee is provided by Armare S.A., a company governed by Polish law, with its registered office in Warsaw, ul. Grzybowska 87, 00-844, Poland. You can find us in the official Polish register of companies under the KRS number: 0000730590, NIP: 5213833614 and REGON: 380721860. We guarantee that the equipment you have purchased from us (hereinafter referred to as the "Product") will work properly. If something goes wrong for a reason inherent in the Product, this warranty protects you. The warranty starts from the day you receive the Product and lasts: 24 months, if you are a consumer or a natural person running a business and the purchase of the Product is not of a professional nature for you. Otherwise, the warranty period is 12 months. If the law in your area requires a longer warranty, the warranty period will be extended to meet local law requirements.

In addition, if the Product does not comply with the contract, you have the right to legal protection from and at the expense of the seller. This warranty does not affect these rights - it does not limit, suspend or exclude them.

Detailed warranty conditions, how to submit complaints, what happens when you find a defect, what are exceptions to the warranty, when the warranty expires or is extended, and what are the rules regarding the processing of personal data, you can find on our website: <http://medivon.pl/warranty>