medivon.

User manual

Medivon Nova

PL | EN | DE | FR | ES | PT | IT | CZ | SK
SL | EE | LT | LV | HU | RO | RU | NB | FI
EL | HR | SV | DA | NL | TR | BG | AR

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Device manual

The massage mat is ready to work immediately after purchase. Just connect the cable to the power supply and turn on the ON / OFF button after lying comfortably on the mat. The device will automatically start working on medium intensity and program 1. Use the remote control to select your favorite massage area according to your preferences or start the heating mode. You can change the intensity of the massage at any time.

The shiatsu massage pillow can be used while massaging on the massage mat, or completely independently, all you need to do is unzip the pillow. Connect the power to the massage cushion and press the ON/OFF button located on the cushion to start it in shiatsu mode with heating. Press the button again to change the mode to massaging without heating. To turn off the pillow, hold the button for at least 3 seconds. Shiatsu massage is a very deep and effective massage, so it is not recommended to massage one area for more than 15 minutes in one session to avoid strains and injuries.

Fastening straps will allow you to lay the mat on an armchair or sofa

Warnings and Safety Instructions

Particular care should be taken by those who have recently undergone surgery or are undergoing treatment. People with hypertension, heart disease, cerebrovascular atherosclerosis and cerebral thrombosis. People with high fever, cancer and osteoporosis. Women who are pregnant or menstruating. Using a pacemaker or foot implant. This product has a surface heating function and should only be used by people who are sensitive to heat. The product is not suitable for children to play with. This product can be used by children 12 years or older. When using this product, people with disabilities or people with limited behaviors should be under the supervision of their guardians. The product must not be used on children or unconscious people. Improper connection of the plug may result in electric shock. If you have any doubts about the safety of this product or the outlet is not suitable for the plug of this product, please contact customer service. If the product gets wet, it cannot be reused.

Do not use the product in a humid environment such as a bathroom or near a water source to avoid malfunction or danger to life and health. Do not use this product at a room temperature higher than 40 degrees Celsius or lower than 5 degrees Celsius, because the thickening of the lubricating oil may increase the wear of the machine transmission parts and affect the service life. Please do not use in the open or with prolonged exposure to hot sun, which may damage the varnish layer. Do not use the device under a cover where the ventilation opening may be blocked. Overheating of the device may occur. Do not carry the device by the power cord or use it as a handle. Never use the device with the air holes blocked. Never insert any objects into any openings. Do not use outdoors. Before use, long hair should be clipped so that it does not come into contact with the device. Do not use the warming function on areas where you have skin disorders, open wounds, swelling or inflammation. Improper use excludes any liability for damage. Do not use this device while driving or sleeping.

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Used electrical and electronic equipment and batteries and accumulators must not be mixed with other waste. This is indicated by the crossed-out bin symbol on the product packaging. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution. Proper handling of used electrical and electronic equipment and batteries and accumulators in the household contributes positively to the reuse and recovery, including recycling, of used equipment. Used electrical and electronic equipment and batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the used equipment is of the same type and has performed the same functions as the equipment sold. A distributor, when delivering household equipment to a purchaser, is obliged to take back used household equipment free of charge at the place of delivery of such equipment, provided that the used equipment is of the same type and has performed the same functions as the equipment delivered. A distributor running a retail outlet with a sales area within the meaning of art. 2 item 19 of the Act of 27 March 2003 on spatial planning and development of at least 400 m² devoted to the sale of household equipment, is obliged to accept, free of charge, in that outlet or in its immediate vicinity, used household equipment, none of the external dimensions of which exceeds 25 cm, without the need to purchase new household equipment. The product was introduced to the market after 13 August 2005.

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Short Declaration of Conformity

The manufacturer declares under its sole responsibility that the Medivon Nova product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility, Directive 2014/35/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits, and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU as regards the list of restricted substances.

The full text of the Declaration of Conformity is available at the following website: www.medivon.pl/deklaracje

Importer: Armare S.A.
ul Grzybowska 87,
00-844 Warszawa, Poland
Made in China
Marka: Medivon
E-mail: serwis@medivon.pl
www.medivon.pl
SN:01265904119281334

Producent: XIAMEN COMFIER TECHNOLOGY CO.,LTD
2-5/F, No.65 Sunban South Road, Jimei District, Xiamen, China.
E-mail: sales@comfier.com

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Shortened warranty terms

This Guarantee is provided by Armare S.A., a company governed by Polish law, with its registered office in Warsaw, ul. Grzybowska 87, 00-844, Poland. You can find us in the official Polish register of companies under the KRS number: 0000730590, NIP: 5213833614 and REGON: 380721860. We guarantee that the equipment you have purchased from us (hereinafter referred to as the "Product") will work properly. If something goes wrong for a reason inherent in the Product, this warranty protects you. The warranty starts from the day you receive the Product and lasts: 24 months, if you are a consumer or a natural person running a business and the purchase of the Product is not of a professional nature for you. Otherwise, the warranty period is 12 months. If the law in your area requires a longer warranty, the warranty period will be extended to meet local law requirements.

In addition, if the Product does not comply with the contract, you have the right to legal protection from and at the expense of the seller. This warranty does not affect these rights - it does not limit, suspend or exclude them.

Detailed warranty conditions, how to submit complaints, what happens when you find a defect, what are exceptions to the warranty, when the warranty expires or is extended, and what are the rules regarding the processing of personal data, you can find on our website: <http://medivon.pl/warranty>