medivon.

User manual

PL | EN | DE | FR | ES | PT | IT | CZ | SK
SL | EE | LT | LV | HU | RO | RU | NB | FI
EL | HR | SV | DA | NL | TR | BG | AR

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User Manual

The device is ready for immediate use after purchase. Make sure the device is charged, then adjust the position of the massager using the top and side knobs to your liking. When the massager is turned off, press and hold the power button briefly to turn it on. The first - default massage mode will be launched. Press the button to cyclically change the operating modes according to your preferences.

Below is a diagram of the cyclic sequence of modes and functions:

- 1) Vibration massage + Air massage + Heating massage
- 2) Vibration massage + Air massage
- 3) Vibration massage + Heating massage
- 4) Air massage
- 5) Heating massage

If the device is discharged, connect the device to a power source to recharge the battery. A full battery will be indicated by a blue icon. The supply voltage is 5V / 1A (USB-C).

To play music, you need to connect to your mobile device via bluetooth. To do this, open the bluetooth settings and search for our head massager in the list of devices. After connecting, we can control the playback using the buttons on the massager. A short press of the buttons plays the next or previous track, and a longer press increases or decreases the volume of the music being played.

The massager is ideal for people who require deep relaxation, want to improve the condition of their hair or improve the circulation of the scalp. With systematic use of the massage, you will reduce the symptoms of fatigue, insomnia or chronic headaches. The device must not be used by persons with injuries or who have undergone surgery in the last 3 months. Do not use the device in high temperatures and damp rooms. Absolutely keep away from water. Also, people with hypertension, hypotension or severe strokes cannot use the device. People with skin problems should consult a doctor before use.

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Used electrical and electronic equipment and batteries and accumulators must not be mixed with other waste. This is indicated by the crossed-out bin symbol on the product packaging. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution. Proper handling of used electrical and electronic equipment and batteries and accumulators in the household contributes positively to the reuse and recovery, including recycling, of used equipment. Used electrical and electronic equipment and batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the used equipment is of the same type and has performed the same functions as the equipment sold. A distributor, when delivering household equipment to a purchaser, is obliged to take back used household equipment free of charge at the place of delivery of such equipment, provided that the used equipment is of the same type and has performed the same functions as the equipment delivered. A distributor running a retail outlet with a sales area within the meaning of art. 2 item 19 of the Act of 27 March 2003 on spatial planning and development of at least 400 m² devoted to the sale of household equipment, is obliged to accept, free of charge, in that outlet or in its immediate vicinity, used household equipment, none of the external dimensions of which exceeds 25 cm, without the need to purchase new household equipment. The product was introduced to the market after 13 August 2005

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The manufacturer hereby declares that the Medivon Hebe device is in compliance with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65/UE of the European Parliament and of the Council EU of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, as well as Commission Delegated Directive (EU) 2015/863 of March 31, 2015 amending Annex II to Directive 2011/65/EU.

The full text of the declaration of conformity is available at the following website:
www.medivon.pl/deklaracje

Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska

Model: Medivon Hebe | SN: 06255904119281433 | Marka: Medivon | email: serwis@medivon.pl