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The temperature in the forehead and temple area differs from the internal temperature, which is taken orally or rectally. Vasoconstriction, an effect which constricts the blood vessels and cools the skin, can occur during the early stages of a fever. In this case, the temperature measured by the A202 Infrared thermometer may be unusually low. If the measurement therefore does not match the patient's own perception or is unusually low, repeat the measurement every 15 minutes. As a reference, you can also measure the internal body temperature using a conventional oral or rectal thermometer. Body temperature can vary from one individual/person to next. It also varies by location on the body and time of day. Below shows the statistical normal ranges from different sites. Please keep in mind that temperatures measured from different sites, even at the same time, should not be directly compared. Fever indicates that the body temperature is higher than normal. This symptom may be caused by infection, overexposure or immunization. Some people may not experience fever even when they are ill. These include, but are not limited to, infants younger than 3 months old, persons with compromised immune systems, persons taking antibiotics, steroids or antipyretics (aspirin, ibuprofen, acetaminophen), or persons with certain chronic illnesses. Please consult your physician when you feel ill even if you do not have fever.

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RO

Infrared thermometer Model: A202 USER MANUAL

PLEASE NOTE:

THIS MEDICAL DEVICE MUST BE USED
ACCORDING TO INSTRUCTIONS
TO ENSURE ACCURATE READINGS.

Thank you for purchasing the Infrared Thermometer model A202. Please read this instruction manual first, so you can use this thermometer safely and correctly. Please keep this instruction manual for future reference. Batch code on the packaging. This innovative medical device uses advanced infrared (IR) technology to measure temperature instantly and accurately on the forehead or object.

IMPORTANT SAFETY INSTRUCTIONS READ BEFORE USE

The following basic safety precautions should always be taken.

1. Close supervision is necessary when the thermometer is used by, on, or near children, handicapped persons or invalids.
2. Use the thermometer only for the intended use described in this manual.
3. Do not use the thermometer if it is not working properly, or if it has suffered any damage.

KEEP THESE INSTRUCTIONS AT A SAFE PLACE

BEFORE YOU BEGIN

Cautions and Warnings

- As with any thermometer, proper technique is crucial to obtaining accurate temperature readings. Please read this manual thoroughly and carefully before using.
- Always operate the thermometer in an operating temperature range 5°C to 40°C (41°F to 104°F), and relative humidity 15 to 93%.
- Always store the thermometer in a cool and dry place -20°C to 70°C (-4°F to 158°F) and relative humidity 15% to 93%.
- The device requires no calibration.
- The device contains no user serviceable parts.
- The user must check that the equipment functions safely and see that it is in proper working condition before being used.
- The manufacturer does not require such preventive inspections by other persons.
- No modification of this equipment is allowed.
- The device is not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.
- Manufacturer will provide circuit diagrams, component part lists, descriptions, calibration instructions to assist to SERVICE PERSONNEL in parts repair.
- Do not clean or carry out maintenance of the device when it is in use.
- Avoid direct sunlight.
- Avoid dropping the thermometer, if it happens and you think the thermometer may be damaged, please contact customer services immediately.
- Do not touch the lens.
- Do not touch the patient's skin surface.
- Do not disassemble the thermometer.
- Basic safety precautions should always be observed, especially when the thermometer is used on or near children and disabled persons.

- This thermometer is not intended to substitute for a consultation with your physician.
- This thermometer and the subject must remain in a stable environment for at least 30 minutes before measuring the temperature.
- When the measured temperature falls within the fever temperature range of $\geq 38.0^{\circ}\text{C}$ (100.4°F) and $< 43^{\circ}\text{C}$ (109.4°F), as indicated by the red LCD on display, please consult with your physician immediately.

RESTRICTIONS OF USE

This thermometer is clinically proven to produce accurate temperature measurements. However, please be advised that the accuracy can not be ensured when the thermometer is not clean. Check that the probe is clean before taking a measurement.

INTENDED USE

A202 Infrared Thermometer is intended for the periodic measurement and monitoring of human body temperature from forehead measurement at home, clinics and hospital. The thermometer is intended for non-contact measurement.

A control measurement using a conventional thermometer is recommended in the following cases:

1. If the reading is surprisingly low.
2. For new-born infants, up to 100 days old.
3. For children under three years of age who have a weakened immune system or who react unusually in the presence or absence of fever.

HOW DOES IT WORK

The thermometer measures the infrared heat generated by the surface of the skin over the vessel and its surrounding tissue. The thermometer then converts it into a temperature value.


NOTE: The thermometer does not emit any infrared energy.

HIGHLIGHTED FEATURES

- Measurement that does not require probe cover, thereby saving cost of replacement.
- Automatically power off if left idle for 10 seconds.
- Memory function allows you to recall previous results up to 25 previous results.
- Easy to read LCD with green backlight in a dark environment.
- Color for visualization of fever (red or yellow) and measurement in progress (green).


Meter Overview

1. Infrared Sensor
2. LED Light
3. Display Screen
4. ON Button / Measurement / Memory / Setup Button
5. Mode Switch
6. Battery Cover

Display Screen

1. Body mode indication
2. Object mode indication
3. Memory indication
4. Temperature reading
5. Temperature unit
6. Battery indication

DISPLAY MODE

Two modes can be selected.

1. Body Mode

This mode is used to measure the forehead temperature.

2. Object Mode

This mode is used to measure the object temperature.

SELECT THE TEMPERATURE UNIT

This meter provides two measurement units used for indicating the body/object temperature, $^{\circ}\text{C}$ or $^{\circ}\text{F}$, for your preferred selection.

Be sure the thermometer is OFF before selecting the temperature unit.

Long-press the Measurement Button for 8 seconds until $^{\circ}\text{C}$ or $^{\circ}\text{F}$ displayed on the LCD panel as figure shown. Re-press the Measurement Button to select $^{\circ}\text{C}$ or $^{\circ}\text{F}$ unit. Then long-press the Measurement Button for 3 seconds, it will automatically save your choice and turn off.

Fahrenheit / $^{\circ}\text{F}$

Celsius / $^{\circ}\text{C}$

TURN ON THE LED LIGHT

This meter provides Led Light to help users placing the thermometer at the correct position.

Be sure the thermometer is OFF before turn ON/OFF the LED Light.

1. Long-press the Measurement Button for 10 seconds until Signal – Led displayed on the LCD panel.

2. Then short press the measurement to turn on or off LED light – when LED is turned on, the blue light will indicate, while it is turned off there is no light. Then long-press the Measurement Button for 3 seconds, it will automatically save your choice and turn off.

BEEP SETTING

Make sure the thermometer is turned off before turning the beep on / off.

1. Press and hold the measurement button for approximately 12 seconds until OFF or ON appears on the LCD panel.

2. Press the measurement button again to turn the beep ON or OFF. Confirm the selected option by holding down the measurement button until the LCD screen turns off completely.

REPLACING THE BATTERY

The thermometer comes with two 1.5 V AAA alkaline batteries. The meter will display to alert you when the meter power is getting low, please follow the steps below to replace new batteries immediately.


1. Remove the battery cover according to the arrow direction.

2. Remove the old batteries and replace with two 1.5V AAA size alkaline batteries. Taking care to match the Positive (+) and Negative (-) indications.

3. Close the battery cover according to the arrow direction.

NOTE:

- Although the thermometer works when appearing, we still recommend that you change the batteries to obtain an accurate result.
- Remove the batteries if stored for a long period of time.

WARNING!

The batteries should be kept away out of children's reach. Risk of swallowing batteries! If they are swallowed, promptly see a doctor for help. Before using the device, make sure that it is not damaged. Do not use the device if it is damaged as this may cause injury. The battery replacement should be carried out only by adults. Be sure to observe the polarity marked in the battery holder. Avoid the occurrence of short circuits of the batteries' terminals. Do not install used batteries; there is a risk of electrolyte leakage and damage to the device. Do not dispose of batteries in normal household waste and places not intended for this. Do not recharge used batteries, as they may explode.

DETAILED INFORMATION

About Normal Body Temperature & Fever

The temperature in the forehead and temple area differs from the internal temperature, which is taken orally or rectally. Vasoconstriction, an effect which constricts the blood vessels and cools the skin, can occur during the early stages of a fever. In this case, the temperature measured by the A202 Infrared thermometer may be unusually low. If the measurement therefore does not match the patient's own perception or is unusually low, repeat the measurement every 15 minutes. As a reference, you can also measure the internal body temperature using a conventional oral or rectal thermometer. Body temperature can vary from one individual/person to next. It also varies by location on the body and time of day. Below shows the statistical normal ranges from different sites. Please keep in mind that temperatures measured from different sites, even at the same time, should not be directly compared. Fever indicates that the body temperature is higher than normal. This symptom may be caused by infection, overdressing or immunization. Some people may not experience fever even when they are ill. These include, but are not limited to, infants younger than 3 months old, persons with compromised immune systems, persons taking antibiotics, steroids or antipyretics (aspirin, ibuprofen, acetaminophen), or persons with certain chronic illnesses. Please consult your physician when you feel ill even if you do not have fever.

Table*1 Normal Temperature Range of various body sites

Oral	0.6°C (1°F) or more above or below 37°C (98.6°F)
Rectal/ear	0.3°C to 0.6°C (0.5°F to 1°F) higher than oral temperature
Axillary (armpit)	0.3°C to 0.6°C (0.5°F to 1°F) lower than oral temperature

PERFORM THE TEST

As a Body Measurement Thermometer

Press "ON" Button to turn on the thermometer first.

Push the Mode switch to select Body mode. The temperature unit flashes.

Move the probe close to the forehead and take measurements. Make sure the probe is flat and close to the forehead, not at an angle. Perform a forehead measurement with a distance within 3 cm.

Read the result. The measurement result will be done in 1 second. The reading is shown together with LED lighting and one long beep informs about the temperature measurement, and after a while another short beep, confirming the saving of the result to the memory and readiness for the next measurement.

Press the "ON" button to turn off the unit, Or leave it idle for 10 seconds, the unit will switch off automatically.

NOTE:

- As the forehead measurement temperature is likely to be affected by sweat, oil and the surroundings, the reading shall be taken as a reference only.
- If the probe is placed at an angle close to the forehead for measurement, the reading will be affected by surrounding temperature. Babies' skin reacts very quickly to the ambient temperature. Therefore, do not take their temperature with the A202 Infrared thermometer during or after breastfeeding, because the skin temperature may then be lower than the internal body temperature.
- If the measured temperature is $< 32^{\circ}\text{C}$ (89.6°F), the display will show with LO symbol.
- If the reading is $\geq 38.0^{\circ}\text{C}$ (100.4°F) and $< 43.0^{\circ}\text{C}$ (109.4°F), the display will show together with red LED and six short beeps.
- The thermometer will automatically turn off if left idle for 10 seconds.

As an Object Measurement Thermometer


Push the Mode switch to select Object mode.
The machine after switching on emits two short beeps indicating that the object mode is turned on and the temperature unit flashes.

Short-press "ON" Switch to turn on the thermometer first.

Move the probe close to the object and take measurements.
Make sure the probe is flat and close to the object, not at an angle. Perform a measurement with a distance within 3 cm. When ready, push measurement button to take measurements.

Read the result.
The measurement result will be done in 1 second. The reading is shown together with LED lighting and one long beep informs about the temperature measurement, and after a while another short beep, confirming the saving of the result to the memory and readiness for the next measurement.

Press the "ON" button to turn off the unit, Or leave it idle for 10 seconds the unit will switch off automatically.

MEMORY

Recalling the Memory

This thermometer stores 25 most recent readings

Be sure the thermometer is OFF before recalling this memory.

Press  for 2 seconds to enter the memory mode.

Each time you press the Memory Button, a result will be displayed in the order of dates (latest result shown first), together with "MEM" and number (from 1 to 25). LED in green, yellow or red will appear according to the memory reading.
When the memory is full, the oldest result is deleted as the new one added. When the last record displayed in the display, press Memory Button again to return the first record.

Exit the memory.

Long press ON Button to exit the memory or keep the meter in idle for 10 seconds to switch off automatically.

Clear the memory

Removing one of the batteries, keep long pressing the measurement button, then reload the battery until "—" on the LCD screen.

MAINTENANCE

Care & Cleaning

- The probe is not waterproof. Please wipe with a clean and dry cotton swab to clean the probe on the inside.
- The body of the thermometer is not water-resistant. Never put the thermometer under a running tap or submerge it into water. Use a soft and dry cloth to clean it. Do not use abrasive cleaners.
- Store the thermometer in a cool and dry location. Free from dust and away from direct sunlight.

FAULT INDICATION		
FAULT OR FAULT SYMBOL	FAULT DESCRIPTION	CORRECTIVE MEASURE
No display On the LCD panel	The battery has run out. Incorrect battery polarity.	Replace the battery. Please note: The (+) side of the battery must face upwards.
Measurement not possible (or an abnormal value is displayed)	The thermometer is not ready.	Wait until the °C symbol is displayed.
An abnormal temperature value is displayed.	The probe tip is dirty or damaged. Did you hear the beep after pressing the ON button?	Clean the probe tip or get it repaired. Wait until you hear the beep before removing the thermometer from the ear or forehead.
LO or HI symbol is displayed	The temperature measured is outside the measuring range. LO-temperature < 32°C (89.6°F) HI-temperature ≥ 43.0°C (109.4°F)	Check that the probe tip is clean and that the thermometer is properly placed on the forehead.
 Symbol is displayed	The battery has run out.	Replace the battery.
 Symbol is displayed	The ambient temperature is outside the operating temperature range or is changing too fast.	To ensure accurate measurement, let the thermometer rest at operating temperature for 30 minutes prior to use.

SYMBOL INFORMATION	
SYMBOL	REFERENT
	Information of manufacture, including name and address
	Serial number
	Catalogue number
	For use follow instructions
	Keep dry

SYMBOL	REFERENT
	Information of authorized representative, including name and address
	Batch code
	Type BF Equipment
	Caution, consult accompanying documents

SPECIFICATIONS	
Model	A202
Measurement range	Body/Forehead: 32.0–43°C (89.6°F–109.4°F) Object: 0.0°C to 99.9°C (32°F to 211.8°F)
Display resolution	0.1°C / 0.1°F
Accuracy	For Forehead mode: ±0.2°C (±0.4°F) from 36.0°C (96.8°F) to 39.0°C (102.2°F) ±0.3°C (±0.5°F) : out of the range For object mode: ±1°C (±2°F) from 0°C (32°F) to 60°C (140.0°F) ±4°C (±7.2°F) : out of the range
Indicator light	Green light for temperature lower than 37.5°C (99.5°F) Yellow light for temperature equal or between 37.5°C (99.5°F) to 38.0°C (100.4°F) Red light for temperature higher than 38.0°C (100.4°F)
Voice	Power on and ready for working: a short beep. Measurement finished: 1 long beep equal or lower than 38.0°C (100.4°F) 6 short beeps higher than 38.0°C (100.4°F)
Silent mode	1. Be sure the infrared thermometer is off 2. Hold the measurement button for about 12 seconds after "ON" shows on the LCD panel. (Note: keep holding the button when "°C/°F" is displayed on the LCD panel). 3. Press the Measurement button again to turn ON/OFF the beep. 4. Then long-press the Measurement Button for 3 seconds, it will automatically save your choice and turn off.
Memory	25 groups
Operating conditions	5 °C to 40 °C (41 °F to 104 °F), Humidity: 15 to 93 % RH (non-condensing)
Air pressure	700 hPa - 1060 hPa (70Kpa-106Kpa)
Storage and transport environment	Temperature: -20°C to 70°C (-4°F to 158°F) Storage humidity: 15% to 93% RH (non-condensing)
Auto shut-off	About 10 seconds after no using
Battery	2pcs 1.5V AAA Alkaline Battery
Dimension (LxWxH)	142.5mmx38mmx40mm
Weight	79g (Net: 57g)
Life time	2 years

DISPOSAL OF DEVICE

Adhere to the applicable regulations when disposing of the device. This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Please remove the batteries before disposing of the device/unit. Do not dispose of old batteries with your household waste, but at a battery collection station at a recycling site or in a shop.

REFERENCE OF STANDARDS

Device Standards:

Device corresponds to the requirements of the standard for infrared thermometers
IEC 60601-1-2: 2014
IEC 60601-1: 2012
EN ISO 80601-2-56: 2017

Classification:

Anti-electric Shock Type: Internally powered equipment

Applied part: Type BF

Mode of operation: Continuous Operation

EMC: type B class I

Electromagnetic Compatibility:

Device fulfills the stipulations of the standard IEC 60601-1-2

The stipulations of EU-Directive 93/42/EEC for Medical Devices Class Ila have been fulfilled.

* Technical alterations reserved!

Software identity no.: A202 V1.0.0

MANUFACTURER'S DECLARATION OF THE EMC

1. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected service life.
2. Guidance and manufacturer's declaration - electromagnetic emissions and Immunity.

Table 1 - Guidance and manufacturer's declaration - electromagnetic emissions

Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not application
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not application

Table 2 - Guidance and manufacturer's declaration - electromagnetic Immunity

Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	Power supply lines ±2 kV input/output lines ±1 kV	Not application
Surge IEC 61000-4-5	line(s) to line(s) ±1 kV. line(s) to earth ±2 kV. 100 kHz repetition frequency	Not application

Immunity Test	IEC 60601-1-2 Test level	Compliance level
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0 0% 300 cycle	Not application
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC 61000-4-6	150kHz to 80MHz 3Vrms 6Vrms (in ISM and amateur radio bands) 80% AM at 1 kHz	Not application
Radiated RF IEC 61000-4-3	10V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz	10V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz

NOTE Ur is the a.c. mains voltage prior to application of the test level.

Table 3 - Guidance and manufacturer's declaration - electromagnetic Immunity

Radiated RF IEC 60904-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)							
Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)		Distance (m)	Immunity test Level (V/m)
385	380-390	TETRA 400	Pulse Modulation 18 Hz	1,8	0,3	27	
450	430-470	GMRS 460, FRS 460	FM ±5kHz deviation 1 kHz sine	2	0,3	28	
710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0,3	9	
745							
780							
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0,3	28	
870							
930							
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3,4, 25; UMTS	Pulse modulation 217 Hz	2	0,3	28	
1845							
1970							
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0,3	28	
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0,3	9	

Radiated RF IEC 61000-4-3
(Test specifications for ENCLOSURE PORT IMMUNITY to RF-wireless communications equipment)

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