medivon.

User manual

PL | EN | DE | FR | ES | PT | IT | CZ | SK
SL | EE | LT | LV | HU | RO | RU | NB | FI
EL | HR | SV | DA | NL | TR | BG | AR

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User Manual

The device is ready for immediate use after purchase. Make sure the device is charged, then start the massager and select the appropriate operating mode according to your preferences. When the massager is turned off, short press the POWER button to turn it on. The first slow speed mode will start. To turn off the device, press the POWER button again.

If the device is discharged, the battery indicator will turn on. Connect the device to a power source to charge the battery. A full battery will be indicated by a green LED. The supply voltage is 5V / 2A (USB-C).

The massage icon pressed cyclically changes the intensity of the air massage. Three massage strengths are available.

Pressing the heating button activates the heating mode, pressing it again disables the function. Using the gear icon, we cycle through one of the three operating modes of the device. They have been automatically prepared for the best effects of hand massage.

The clock icon is used to set the timer of the device. The available time is: 5/10/15 min.

With systematic use of the massage, you will reduce the symptoms of hand fatigue, insomnia or chronic nervous pain. The massager is also ideal for people who spend a lot of time working at the computer.

The strength of the massage should be adjusted to your own preferences. The device must not be used by persons with injuries or who have undergone surgery in the last 3 months. Do not use the device in high temperatures and damp rooms. Absolutely keep away from water. Also, people with hypertension, hypotension or severe strokes cannot use the device. People with skin problems should consult a doctor before use.

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Waste electrical and electronic equipment as well as batteries and accumulators cannot be mixed with other waste. This is indicated by the crossed-out bin symbol on the product packaging. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage may result in health damage or environmental pollution.

Proper household handling of waste electrical and electronic equipment and batteries and accumulators positively contributes to the reuse and recovery, including recycling, of waste equipment.

Waste electrical and electronic equipment as well as batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the equipment is of the same type and performed the same functions as the equipment sold. When delivering equipment intended for households to the buyer, the distributor is obliged to collect used equipment from households free of charge at the place of delivery of the equipment, provided that the used equipment is of the same type and performed the same functions as the delivered equipment. A distributor running a retail unit with a sales area within the meaning of Art. 2 point 19 of the Act of March 27, 2003 on spatial planning and development of at least 400 m² devoted to the sale of equipment intended for households, is obliged to accept waste equipment from households, free of charge, in this unit or in its immediate vicinity, none of the external dimensions exceeds 25 cm, without the need to purchase new household equipment. The product was introduced to the market on August 13, 2005

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The manufacturer declares, under its sole responsibility, that the Medivon Hand Lux device complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, Directive 2014/35/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to the making available on the market of electrical equipment intended for use within specified voltage limits, and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances. The full content of the declaration of conformity is available at www.medivon.pl/deklaracje

Importer: Armare S.A.
ul Grzybowska 87, 00-844 Warszawa, Poland.
Made in China,
Marka: Medivon
Model: Medivon Hand Lux
email: serwis@medivon.pl
SN:06255904119281570

Producent: Shenzhen Forrest Health
Technology Co., Ltd
Room 201, Building 1, No.30, QinFu
Road, Liu Yue South Community, Heng-
Gang Street, LongGang District, Shen-
Zhen, China.
email: jason@forrestmassager.com