**medivon.**

## **User manual**

PL | EN | DE | FR | ES | PT | IT | CZ | SK  
SL | EE | LT | LV | HU | RO | RU | NB | FI  
EL | HR | SV | DA | NL | TR | BG | AR

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## User Manual

The device is ready for immediate use after purchase. Make sure the device is charged, then start the massager and select the appropriate operating mode according to your preferences. When the massager is turned off, short press the POWER button to turn it on. The first default operating mode will be launched. To turn off the device, hold down the POWER button again. When the device is turned on, short cyclical presses of the power button change the predefined operating modes with various functions of the device.

If the device is discharged, the battery indicator will turn red. Connect the device to a power source to charge the battery. A full battery will be indicated by a green LED. The supply voltage is 5V / 2A.

The massage icon pressed cyclically changes the intensity. Three massage strengths are available to suit the user's preferences.

Pressing the heating button activates the heating mode. There are 3 heating levels. The pressing cycles alternate in sequence on the display. The fourth press turns off the heating function.

Using the vibration button, as in the case of heating, we change its intensity cyclically. A fourth press turns off the vibration function.

When the heating or compression massage function is turned on, EMS mode can be activated. The EMS+ and EMS- buttons are used to control the intensity of electromassage. The vibration and EMS frequencies have been automatically adjusted by specialists for the best possible use.

The strength of the massage should be adjusted to your own preferences. The device must not be used by persons with injuries or who have undergone surgery in the last 3 months. Do not use the device in high temperatures and damp rooms. Absolutely keep away from water. The device cannot be used by people with hypertension, hypotension, during pregnancy, with a pacemaker or after severe strokes. People with skin problems should consult a doctor before use.

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Used electrical and electronic equipment and batteries and accumulators must not be mixed with other waste. This is indicated by the crossed-out bin symbol on the product packaging. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution. Proper handling of used electrical and electronic equipment and batteries and accumulators in the household contributes positively to the reuse and recovery, including recycling, of used equipment. Used electrical and electronic equipment and batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the used equipment is of the same type and has performed the same functions as the equipment sold. A distributor, when delivering household equipment to a purchaser, is obliged to take back used household equipment free of charge at the place of delivery of such equipment, provided that the used equipment is of the same type and has performed the same functions as the equipment delivered. A distributor running a retail outlet with a sales area within the meaning of art. 2 item 19 of the Act of 27 March 2003 on spatial planning and development of at least 400 m<sup>2</sup> devoted to the sale of household equipment, is obliged to accept, free of charge, in that outlet or in its immediate vicinity, used household equipment, none of the external dimensions of which exceeds 25 cm, without the need to purchase new household equipment. The product was introduced to the market after 13 August 2005.

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### **Short Declaration of Conformity**

The manufacturer declares, that the product Medivon Knee Pro complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility, as well as with Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU with regard to the list of restricted substances. The full text of the Declaration of Conformity is available at: [www.medivon.pl/deklaracje](http://www.medivon.pl/deklaracje)

**Importer:** Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska | Model: Medivon Knee Pro  
**Marka:** Medivon | **Made in China** | SN:01265904119281587 | [www.medivon.pl](http://www.medivon.pl) | mail: [serwis@medivon.pl](mailto:serwis@medivon.pl)

**Producent:** Xiamen Healthpal Electronic Co., Ltd.  
NO.148 BINGLANGXILI, SIMINGDISTRICT, XIAMEN,FUJIAN, CHINA  
E-mail: [vilin@xmhealthcare.com](mailto:vilin@xmhealthcare.com)

Detailed warranty conditions, how to submit complaints, what happens when you find a defect, what are exceptions to the warranty, when the warranty expires or is extended, and what are the rules regarding the processing of personal data, you can find on our website:  
<http://medivon.pl/warranty>