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User manual

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Medivon Gaia

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Everything for your health

The accupressure mat, also known as a spiked mat, is a modern and more delicate version of an acupuncture bed. This is currently one of the best-selling models on the European market.

It has helped thousands of people with ailments such as back pain, sciatica, muscle strain, insomnia, depression and fatigue. After the first use, you will find out about its beneficial properties without the use of pharmacological agents that destroy the body.

The Medivon mat has thousands of acupuncture points that release endorphins - a neuropeptide of happiness and well-being that also blocks pain. The result is a release from pain and tension and a feeling of deep relaxation.

It can be used when you lack energy, suffer from high levels of stress, back pain, insomnia, muscle tension, lumbago, sciatica, migraine, digestive problems, menstrual cramps, depression, etc.

Functionality and application

- Reduces pain in the back, neck, shoulders, hips and joints
- Improves blood circulation
- Increased energy levels
- Gives a deep sense of relaxation
- Reduces muscle tension
- Increases the level of endorphins (the feel-good hormone)
- Increases oxygenation of body cells
- It improves sleep and helps with insomnia
- Works well for improving digestion
- It gives a positive mental attitude

Q&A

Can I wash the mat in the washing machine if the mat is dirty?

Don't do this as the washing machine will destroy the acupressure spikes.

How to use it, I'm afraid of the pain because the mat is quite sharp?

You can put the fabric on the mat and then lay it on top of it.

How long will it take me to get a good result?

About 40 minutes, but time is different for different users.

Can I still lie on it if I feel pain at first?

If you have a problem with your skin, don't use the mat. If not, you can still lie on it, the pain starts for about 5 minutes, then slowly fades away and turns into blissful relaxation with a warm feeling.

After use, I have marks on my skin - is this normal?

The red marks are generated by acupuncture points, but don't worry, they will disappear automatically after a short time.

Can I use an acupressure mat if I have had heart disease?

You can not. Because the acupuncture point will needle the acupuncture points when you lie on the mat, it will then increase blood circulation and speed up the flow of blood.

How to clean the acupuncture mat when it is dirty?

Use a brush to remove dust, clean the mat with laundry detergent, then hang in the fresh air to dry. If you hang it in the sun, the fabric may fade.

Material:

Mat: 100% cotton

Filling: sponge

Acupressure points: polypropylene (ABS)

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The crossed-out wheeled bin symbol indicates that in the European Union after the end of the product's useful life it must be disposed of in a separate, dedicated collection point. Do not dispose of these products together with unsorted general household waste

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English

Hereby, Medivon Sp. z o. o. declares that the device equipment Medivon Gaia is in compliance with Directive. The full text of the EU declaration of conformity is available at the following internet address:

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TERMS OF WARRANTY

1. The product warranty is valid for a period of 24 months.
2. The warranty covers the period from the date of receipt of the device from Medivon Sp. z o. o. and covers defects resulting from causes inherent in the received device.
3. The device will be repaired as soon as possible, not exceeding 21 working days from the date of acceptance of the device for repair by an Authorized Service Center.
The customer has the right to apply for replacement of equipment free from defects, if the
4. Authorized Service finds that the removal of the defect is impossible or requires excessive costs.
5. If only part of the product is defective and can be detached from the part of the product working properly, the Customer's rights resulting from this warranty limit only to the defective part of the product.

In order to benefit from the warranty conditions, please send a complaint notification together with the warranty card to the e-mail address serwis@rgmedia.pl. The exact warranty conditions are under the address:

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I represent that I have read and accepted the terms and
conditions specified in the Warranty Certificate.

.....