medivon®

User manual

Hand massager

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User Manual

The massager is very easy to use - just one press of a button to start a relaxing hand massage, which greatly improves your blood circulation and relaxes the tense muscles of your hand after a hard day's work. The device is ready for immediate use after purchase. Make sure it's charged, insert your hand into the massager, then turn it on and choose the appropriate operating mode according to your preferences. When the massager is off, press and hold the POWER button for 2 seconds to turn on the device. To turn off, press the POWER button again.

Modes and intensity of the massager - the device has 6 modes to choose from, which can be changed with the (M) button, and 6 intensity levels adjusted by the (L) button.

Fan - to turn on the internal fan in the device, just press the button with the fan icon. To turn it off, press the button again.

Heating function - press the heating button to choose the level: low, high, or off.

Timer function - the device has built-in 3 timers: 10, 15, and 20 minutes, which can be changed by pressing the timer button. The device will automatically turn off when the timer expires.

Kneading function - press the button briefly to activate the kneading. The green light indicates kneading in the direction of the clock hands. The blue light indicates kneading opposite to the direction of the clock hands.

Safety air valve switch - if the massager's airbag cannot be normally deflated, switch the valve located at the bottom of the device to safely release the air in an emergency.

Caution: Adjust the massage strength to your preferences. The device should not be used by individuals with injuries or those who have undergone surgery in the last 3 months. Individuals with high blood pressure, low blood pressure, or after severe strokes should also refrain from using the device. People with skin problems should consult a doctor before use. This product is not a toy and is not suitable for children. Store the device at room temperature, in a dry and shaded place, out of the reach of children.

Maintenance and care: do not use the device in high temperatures or humid rooms. Always unplug the power adapter and let the device cool down after each use or before cleaning. To clean the device, wipe it with a soft, damp cloth. Keep strictly away from water.

Battery: only use the charger provided with the device to charge the battery. Before scrapping, remove the battery from the device. The massager should only be used with the charger provided with the device.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment. The device was placed on the market after 13.08.2005.

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The manufacturer hereby declares that the Medivon Revit device complies with Directive 2014/35/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits, Directive 2014/30/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65/EU of the European Parliament and of the Council of June 8, 2011 on restrictions on the use of certain hazardous substances in electrical and electronic equipment. The full content of the declaration of conformity is available at the following internet address: www.medivon.pl/deklaracje

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Importer: Armare S.A. | ul. Grzybowska 87, 00-844 Warszawa, Polska