medivon.

User manual
EYE MASSAGER

PL | EN | DE | FR | ES | PT | IT | CZ | NL
SV | SK | SL | EE | LT | LV | HU | RO | RU
NR | FI | EL | HR | DA | BG | TR | AR

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1. Construction of the device

1. Adjustment buckle
2. Elastic band
3. ON/OFF key/mode button
4. USB-C input
5. LED indicator

2. How to use

1. Before first use, please fully charge the device using the included USB cable.
2. Put the massager on your head and adjust the headband for comfortable use.
3. Press and hold the power button to turn on the product. Once activated, the button will act as a mode change function button. Each short press will change the operating mode:
 - Mode 1 Comprehensive Mode (default at power on)
 - Mode 2: Eye protection mode (vibration)
 - Mode 3: Eye protection mode (rhythmic vibration)
 - Mode 4: Sleep mode (gentle vibration)
4. Once the massage is complete, hold down the power button to turn off the device.

3. Charging and removing batteries

If you hear three consecutive beeps while using the device, it means that the built-in lithium battery has discharged. To charge it, you need to use a USB cable. While charging, the charging indicator flashes. It takes approximately 2-3 hours to fully charge. Battery life will be shortened if the battery is not used for a long time. It is recommended to charge the battery at least once a month.

The built-in lithium battery of this product contains substances that may pollute the environment. Please remove the built-in lithium battery before disposing of this product. Follow illustration 2 and the instructions:

1. Turn off or disconnect the power.
2. Remove the protective flap.
3. Remove the lithium battery.
4. The removed lithium battery should be safely disposed of at a location designated by a government or recycling facility.

4. Notes, warnings, contraindications

Please clean the product's eye shield before use. To clean, wipe the product with a slightly damp, soft cloth or towel. Do not use any corrosive substances.

People with sensitive skin or insensitive to heat should not use the product. Pregnant women, people suffering from heart disease, circulatory disorders, malignant tumors, cerebrovascular diseases and other acute diseases, as well as people under the care of a doctor, should not use the product. If you experience any discomfort or abnormal skin reaction while using the product, stop using it immediately and consult your doctor. People with eye problems should avoid using the product or use it under medical supervision. People who have undergone eye surgery and patients suffering from glaucoma, cataracts, retinal detachment and other eye diseases should not use the product. People wearing contact lenses should remove them before using the product. Store the product in a dry and airy place, away from children. Do not expose it to direct sunlight, high temperature or moisture. This is not a toy. Children and people unable to express their needs should not use the product. Do not use the product in places with high temperature, flammable, humid or places with electromagnetic radiation. Do not disassemble the product without permission; Please make sure the product is free from impact and water to avoid failure.

5. Product specification

Name: Eye massager

Model: Horizon Travel 2

Power supply: 5V 1000mA

Lithium battery: 3.7V 700mAh

Input power: 5W

Package size: 200x91x85mm

Net weight: about 190g

Product normal operating conditions:

Ambient temperature: +5°C ~ +40°C

Relative humidity: ≤ 80%

Atmospheric pressure: 86 kPa ~ 106 kPa

Product transport and storage conditions:

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Used electrical and electronic equipment and batteries and accumulators must not be mixed with other waste. This is indicated by the crossed-out bin symbol on the product packaging. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution. Proper handling of used electrical and electronic equipment and batteries and accumulators in the household contributes positively to the reuse and recovery, including recycling, of used equipment. Used electrical and electronic equipment and batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the used equipment is of the same type and has performed the same functions as the equipment sold. A distributor, when delivering household equipment to a purchaser, is obliged to take back used household equipment free of charge at the place of delivery of such equipment, provided that the used equipment is of the same type and has performed the same functions as the equipment delivered. A distributor running a retail outlet with a sales area within the meaning of art. 2 item 19 of the Act of 27 March 2003 on spatial planning and development of at least 400 m² devoted to the sale of household equipment, is obliged to accept, free of charge, in that outlet or in its immediate vicinity, used household equipment, none of the external dimensions of which exceeds 25 cm, without the need to purchase new household equipment. The product was introduced to the market after 13 August 2005

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The manufacturer hereby declares that the Medivon Hand device complies with Directive 2014/30/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, Directive 2011/65/EU of the European Parliament and of the Council of June 8, 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, and Commission Delegated Directive (EU) 2015/863 of March 31, 2015 amending Annex II to Directive 2011/65/EU. The full content of the declaration of conformity is available at the following address: www.medivon.pl/deklaracje

Model: Horizon Travel 2
SN: 06255904119283949
Marka: Medivon
Voltage: DC5V | Power 5W
Importer: Armare S.A.
ul. Grzybowska 87
00-844 Warszawa,
Poland
Made in China
e-mail: serwis@medivon.pl