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- www.beautifyly.eu
- www.youtube.com/beautifylyeu

Adsorption pressure:	60KP
Capacity:	30ml
Packed:	5V/1A
Power:	6W
Battery Capacity:	900mA
Landing	standard USB charging
Included:	device, USB cable, water tank, user manual, 6 tips in different sizes, non-woven fabric pads

Security:

Before using the product, read this user manual and follow the recommendations contained therein when using the device.

Device description [page 2-f1]:

1. Replaceable head
2. Indicator light
3. Heat control/heat control button
4. Power button
5. Charger input
6. Water reservoir

Hydrogen purification (hydrabrasion)

Hydrabrasion is a gentle treatment that cleanses, moisturizes and exfoliates the skin. It is definitely more gentle than microdermabrasion, so even people with sensitive, mature, vascular or acne-prone skin can use it. In order to achieve maximum cleansing effects, use special tips dedicated to specific parts of the face. Thanks to them, you can cleanse the skin even in hard-to-reach places and achieve the expected results much easier

Device operation:

1. Press and hold the power button for 3 seconds to turn your device on or off.
2. Briefly hold the power button to increase or decrease the intensity level.
3. The first intensity level turns on one green light, the second intensity turns on two green lights, and the third intensity turns on three green lights.
4. While charging, the charging indicator lights red. When the battery is full, the red light turns off and the green light turns on.
5. Charging time (full charge) of the device is approximately 3.5 hours.

Precautions:

1. Do not place the device in an environment with a strong magnetic field.
2. Be careful not to tilt the device more than 45°C when using, otherwise water will spray out of the exhaust port.
3. This device can only use clean water. Do not use other liquids that are not suitable for this device.
4. Do not use the device while it is charging.
5. The inside of the water tank can be disinfected and sterilized.
6. This device should not be soaked or immersed in water for cleaning. Please be careful when using the device and make sure that no water gets inside when using or cleaning it.
7. After each use of the device, empty the water tank of used water and then dry the device thoroughly.

How to use:

1. Remove makeup from your face using a makeup remover. We recommend that you also use a mask or pore-opening liquid before the treatment.
2. Fill 2/3 of the tank with clean water.
3. Attach the appropriate cleaning tip. Hold the device in a vertical position. Make smooth movements and do not stop the head in one place. After the procedure, remove the cylindrical filter from the hose with tweezers or a needle.
4. After cleansing, wash your face with tonic and apply a moisturizer or mask to your skin. Remember to protect yourself with SPF. Hydrabrasion treatment can be performed all year round. For best results, it should be performed regularly, 1-2 times a week, depending on the skin's needs. The duration of the treatment is approximately 5-10 minutes.

Function description:

1. Cleansing tip - removal of blackheads, dead skin, calluses, cleansing pores. Area of application: cheeks, forehead
2. Cleaning tip (increased suction power) - removing blackheads, cleaning pores. Area of use: T zone, nose part
3. Spiral cleaning tip - removal of dead skin, cleansing pores. Area of application: cheeks, forehead
4. Silicone cleaning tip - removing acne and calluses, cleansing pores. Area of application: cheeks, forehead. Tips made of medical silicone.

Device working angle: See page 2-f2.

More information about the company and products [page 2 - f3]:

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- www.beautify.eu; - www.youtube.com/beautifyeu

EN: Electrical and electronic devices, as well as batteries and accumulators, must not be mixed with other household waste. This is indicated by the crossed-out trash can symbol on the packaging. Used electrical and electronic devices, batteries, and accumulators must be taken to local selective waste collection points or to the retailer, and more information can be found in your municipality. Electrical and electronic equipment, batteries, and accumulators may contain hazardous substances, mixtures, and components harmful to

the environment and human health. Improper handling or damage to them can lead to health hazards or environmental pollution during further disposal or recycling. The device was introduced to the market after August 13, 2005.

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[EN] Summary declaration of conformity: The manufacturer hereby declares under his sole responsibility that the above product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of legislation. Member States in the field of electromagnetic compatibility and Directive 2011/65/EC of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment and Commission Delegated Directive (EU) 2015/863 of 2011. On 31 March 2015, amends Annex II to Directive 2011/65/EU of the European Parliament and of the Council on the list of restricted substances.

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Shortened warranty terms

This Guarantee is granted by a Polish company, Beautifly Sp. z o. o. based in Warsaw. Plac Bankowy 2, 00-095 Warsaw, Poland. You can find us in the official Polish register of enterprises under the KRS number: 0000887306, NIP: 5252854218 and REGON: 388382021 We guarantee that the equipment you have purchased from us (hereinafter referred to as the "Product") will work properly. If something goes wrong for a reason inherent in the Product, this warranty protects you. The warranty starts from the day you receive the Product and lasts:

24 months, if you are a consumer or a natural person running a business and the purchase of the Product is not of a professional nature for you. Otherwise, the warranty period is 12 months. If the law in your area requires a longer warranty, the warranty period will be extended to meet local law requirements.

In addition, if the Product does not comply with the contract, you have the right to legal protection from and at the expense of the seller. This warranty does not affect these rights - it does not limit, suspend or exclude them.

Detailed warranty conditions, how to submit complaints, what happens when you find a defect, what are exceptions to the warranty, when the warranty expires or is extended, and what are the rules regarding the processing of personal data, you can find on our website:

www.beautifly.eu/warranty-beautifly