medivon.

Neck massager

User manual

PL | EN | DE | FR | ES | PT | IT | CZ | NL
SV | SK | SL | EE | LT | LV | HU | RO | RU
NR | FI | EL | HR | DA | BG | TR | AR

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Device operation

The device is very easy to use - just hold the "-" button for a short moment to start a relaxing neck massage, which will greatly improve your blood circulation and relax tense muscles after a hard day's work.

Changing the massage intensity - to change the massage level, simply press the appropriate + or - button on the control panel. The intensity will vary from 1-15

Changing the massage mode - the massager has 4 defined automatic massage modes. By pressing the "M" button, we can select cyclically changing programs, the one that currently suits us best.

Warming up function - built-in heads act as a warming compress, this function is activated with heating button on the left side of the control panel. After activating the warm-up function, the temperature in the device will start to increase gradually. To increase the heating temperature or turn it off completely, press the heating button repeatedly.

To fully enjoy the complete neck massage, it is recommended to gently wipe the skin surface with a damp cloth before the massage. It is not recommended to massage the neck area for more than 30 minutes a day.

It is strictly prohibited to use the massager in the following conditions:

1. Embedded medical electronic devices such as pacemaker;
2. Medical electronic life support devices such as artificial heart and lungs;
3. Medical electronic devices such as electrocardiograph;
4. Flammable places such as a gas station, a workshop with high-density dust or a gas plant.

This will result in incorrect operation of medical electronic devices, which may pose a threat to the life and health of the user or people around him. If the massager is used with high-frequency surgical devices, the area around the electrode site may burn or damage the massager.

The performance of the massager may be unstable if the product is used near a shortwave or microwave therapy device (within 1 meter). During the massage, keep the electrode plate away from wounds and scars. When massaging, the

electrode plate should not be near the chest or heart, as this may cause cardiac fibrillation, modify without authorization. Please do not immerse it in water or throw it into fire, otherwise the product may explode.

Warnings and dangers

Before using the massager, please follow the following safety and use guidelines: Use only parts approved by the manufacturer. Servicing the product requires a permit and should be performed by professionals. If the power cord is damaged, it must be replaced by qualified personnel to prevent a risk. Avoid using or storing the product in humid or hot places. Once you have finished the massage, unplug the power cord and do not use the massager for more than 15 minutes at a time. If any adverse reactions occur, stop using the massager immediately. If the massager begins to emit an unpleasant or burning smell, unplug it and stop using it. Do not move the product by pulling the cord to avoid loosening. The product is not intended for people with physical, sensory or mental limitations, including children. People who are not sensitive to heat should use the heating function with caution. Consult your doctor before use if you are pregnant, have skin problems, use a pacemaker, suffer from heart disease, cancer, osteoporosis, or have blood pressure problems. Before cleaning the product, turn it off and unplug the power cord. Clean the surface with a dry, soft cloth and avoid corrosive detergents. Store in a cool, dry place, away from direct sunlight and moisture to prevent color fading or damage.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment.

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The manufacturer hereby declares that the Medivon Ultra EMS device complies with Directive 2014/35/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits, Directive 2014/30/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65/EU of the European Parliament and of the Council of June 8, 2011 on restrictions on the use of certain hazardous substances in electrical and electronic equipment. The full content of the declaration of conformity is available at the following internet address: www.medivon.pl/deklaracje

Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska
Model: Medivon Ultra EMS