medivon®

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Neck Electronic Muscle Stimulator

Pain Relief and Body Recovery

PL | EN

1. Product Introduction

The design of Neck Electronic Muscle Stimulator is mainly based on the principle of nourishing qi and blood according to traditional Chinese medicine. It helps to promote blood circulation to remove meridian obstruction, dissipate wound, relieve pain and other types of muscle discomfort. It adopts the summarized long-term clinical practice experience as well as combines with modern electronic science and technology, developing a new generation of high-tech products neck prevention. The streamlined ring design is accord with body neck physiological curvature. The device also adopts LF electrotherapy, to generate an efficient composite energy field and through human body's deep tissues and marrows. In the whole treatment process, its physiological positional always maintain the correct posture , providing a scientific way to the treatment of cervical vertebra disease. It will make the user enjoy multifunction system treatment at home with relaxation, efficiency and pleasure.

2. Product Intended use

1) Intended use: Promote neck blood circulation and relieve pain.

Contact with skins of neck. Use Time :15-30min.

2) Intended user: Adults aged between 18 and 60 need to relieve muscle fatigue and neck pain and promote blood circulation.

3) Benefit: Relieve pains; Promote blood circulation; The neck disability index decreased effectively; Enhance cervical range of motion (ROM).

4) Attention:

In case of any adverse events or serious complaints in the EU market, please contact EC Representative and the manufacturer as soon as possible, so that the manufacturer can confirm, handle and reply in time.

3. Contraindication

1) the device is contraindicated for use on patients with cardiac pacemakers fitted.


2) Treatment parts have the haemorrhage after acute trauma or fractures, disabling wound healing stage after the operation.

3) It is contraindicated for use by those who has a skin perceptual disturbance or is not sensitive about the heat.

4) It is prohibited to use when bathing, sweating and sleeping.

- 5) The patient with cerebral hemorrhage: It should be contraindicated in unsteady phase; the person who has sequela must be used under supervision of doctors.
- 6) It is contraindicated for use by those who has purulent inflammation, acute blood poisoning and continuous hyperpyrexia.
- 7) It is contraindicated for use by those who has acute cardiovascular and cerebrovascular diseases.
- 8) Have medical electronic equipment, such as cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device. Such use could cause electric shock, burns, electrical interference, or death.
- 9) Do not use on the open wounds, rashes, swollen, red, infected or inflamed areas or skin eruptions (e. g., phlebitis, thrombophlebitis, varicose veins); On top of, or in proximity to, cancerous lesions.
- 10) Metal allergies are contraindicated.


4. Primary Structure of Product


Remote control

- ① Electric pulse
- ② Heat moxibustion
- ③ Intensity adjustment
- ④ Switch /timing selector
- ⑤ Electrical pulse massage
- ⑥ Increase intensity
- ⑦ Decrease intensity
- ⑧ Heat moxibustion
- ⑨ LCD display

Host

- ① Power/charging switch
- ② Power connector
- ③ Conductive/heating plate
- ④ Red light irradiation device

USB cable

Instruction Manual
+with warranty card

5. Product Main Features

■ The product is designed on the basis of ergonomics which is in line with the neck curve and aesthetics and science.

■ It adopts the latest 3D intelligent mounting technology, and according to your neck curve adjust electrode easily, so that using is more comfortable;

■ There are two automatic combinational massage modes and four manual massage modes to be fit for wider range of people, and it can be selected according to different needs.

■ It contains a LED indicator light displayer for simple clarified operation, which is a humanized design.

6. Product Requirements and Major Parametric Description

6.1 Product Power Requirement

a) The mainframe is charged with DC 5V1000mA power supply
Lithium battery specification: 3.7V DC:1500mAh

Size: 8.6X33X43mm

Weight: about 25g

Power supply of remote control: CR2032 3V DC of Button battery b)

Device safety classes: class II type BF

6.2 Product Main Technical Parametric description

a) Impulse frequency: 1Hz~400Hz.

(Pulse modulation Modulation frequency 4KHZ-10KHZ)

b) Pulse width: 20μs~200μs.

c) Impulse waveform: square wave

d) Pulse duration: ≤3 seconds

e) The product output has no the DC component.

f) The individual pulsing electric quantity when in the maximum output amplitude: >7μC.

g) The maximum output energy of individual pulse: ≤300mJ.

h) The effective value of maximum output amplitude: ≤25V (50mA).

i) In the open measurement, the output voltage peak: ≤500V.

j) Influence of output end open circuit and short circuit: it is able to support the influence of output end open and short circuit, and its performance will not be allowed to weaken.

- k) Adjustment of output amplitude: continuous and balanced, the minimum output is not more than 2% of maximum output.
- l) Rated load impedance: 500Ω , $\pm 10\%$ franchise.
- m) The timer of therapeutic equipment: 15 minutes, $\pm 10\%$ franchise.
- n) Thermal therapy temperature: it should be in the range of $38^{\circ}\text{C} \sim 48^{\circ}\text{C}$. (under the $23^{\circ}\text{C} \pm 2^{\circ}\text{C}$ environment)
- o) Safety classification: internal electric source class and type BF equipment
- p) Boundary dimension:
Host (length 155mm x width 139mm x height 43mm) Controller (length 91mm x width 42mm x height 12 mm)
- q) Product software version No.: A/0
- r) Service life: ≥ 3 years.

6.3 Product Environmental Requirement

- a) Normal work environmental requirements: Environment temperature: $+5^{\circ}\text{C} \sim +40^{\circ}\text{C}$; Environment humidity: 15%-93%RH; Atmospheric environment conditions: 700hPa-1060hPa.
- b) Storage environment requirements: Environment temperature: $-25^{\circ}\text{C} \sim +55^{\circ}\text{C}$; Environment humidity: 0-93%RH; Atmospheric environment conditions: 700hPa-1060hPa.
- c) Transport environment requirements: Environment temperature: $-10^{\circ}\text{C} \sim 40^{\circ}\text{C}$; Environment humidity: 15%-93%RH; Atmospheric environment conditions: 700hPa-1060hPa.

7. Method of application

7.1 Charging of the built-in rechargeable battery.

- 1) The built-in rechargeable battery can serve as the power supply when the external power supply is not used.
- 2) Power supply with the built-in rechargeable battery (applicable to rechargeable ones).when it is not used as an external power supply ,it could be used for power supply.
- 3) The remote control uses the CR2032 Button battery, the positive pole upwards when is installed.

7.2 Functional Operation

1) The cervical part and conductive plate with a wet towel wipe wet, then therapeutic apparatus adheres closely to the neck, if there is a perspiration in cervical part, must be clean to use.

2) Press the "⏻" button by manual remote control, the host power indicator is lightened and working routine is started and the LED light of controller is lightened, that has the following setting will enter an automatic mode, that is:

A. MED-heat;

B. mode 1, the intensity 0; C. time for 15 minutes.

If you want to stop using, you can press "⏻" button to shut down.

When the remote control is operated, its infrared launch window must be corresponding with an infrared launch window of therapeutic apparatus, the apparatus will have a hint of "DI" while is receiving an infrared signal successfully.

3) Intensity adjustment

press the "+"、"-" button to switch the intensity of electric pulse or vibration, there are 16 gears intensity of electric pulse.

4) Heating mode switching

press the "🔥" button, switch the heating mode, there ishi-heat (3), med-heat (2), low-heat (1) and no heating (0).

5) Electric pulse mode switch

press the "⚡" button switch the electric pulse mode, press the "+", "-" button switch the electric pulse mode, and there are six modes, respectively as follows:

Electric pulse mode 1: as the auto-cycle mode, cyclical scart, massage, malaxation mode; Electric pulse mode 2: as the auto- cycle mode, cyclical massage, scart, thump mode;

Note: Electric pulse mode 1 and mode 2 are circulation the follow four modes waveform. Electric pulse mode 3: as the scart mode; Electric pulse mode 4: as the massage mode;

Electric pulse mode 5: as the malaxation mode; Electric pulse mode 6: as the thump mode;

5) The adjustment of timing

press the ⌚ button switch the timing, the adjustable range is 5 to 30 minutes, it can increase or decrease 5 minutes each pressing.

Precautions

It is suggested that using the therapeutic apparatus 2 times a day, vibration and electric pulse 10~15 minutes every time, heating 15~30 minutes every time.

7.3 Use of silicon electrode slices

1. Insert one end of the 2-core connector wire into the electrode port and connect the other end to the electrode slices; then, use the massager as specified above.

If they are used together with the conductive plates, there will be a little difference.

2. The electrode slices should be clean on their surfaces without any dust, oily or sticky substance; otherwise, the viscosity will be reduced. If their surfaces become dirty, please flush with clean water and then dry them in the air for repeated use.

3. The electrode slices could be used for about 80 times (20 minutes for each time).

If their viscosity still cannot be recovered after several times of cleaning, please consider to purchase new ones from your local dealer or the manufacturer.

4. Only the electrode slices provided for the massager, each with a size of 50mm wide, 62mm high and 2.5mm thick, could be used; and others are not allowed.

5. Please dispose the waste electrode slices according to your local environmental.

7.4 The use of rechargeable batteries

1) In the course of using the therapeutic apparatus, it shows the electric quantity of internal lithium battery is not enough when you heard "DI" buzzing 3 times in row, you can charge the battery using adapter and turn the toggle switch to "charge" position, then the green indicator of the power supply is flashing; after the battery is charging fully, the green indicator is normal on. The charging time is about 3~4 hours.

2) It can be used about one week (30 minutes a day) after fully charged.

3) Under a normal using condition, the rechargeable batteries life $\geq$ 500 times charge and discharge, can use about 3 years (30 minutes a day).

4) When the rechargeable batteries are not used for a long time, it will reduce the service life of rechargeable batteries, it is suggested that at least have a full charging every month.

5) The discarding method of the post-batteries should be deal with according to the urban environmental protection.

7.5 The use of button battery to remote control

1) The remote control adopts the“CR2032”button battery, can use about 60days (according to use 30 minutes a day).

2) When pressing the button, the brightness of LCE backlight of remote control is distinct weakening or can only remote close to the host infrared launch window, showing the electric quantity of battery is not enough and please replace the new battery.

3) If the product is not used for a long time, please take out the battery, otherwise it may cause the battery leakage and cause the machine error.

4) Do not use the outside of the specified battery.

5) Do not misplace the positive and negative pole of battery, the positive pole up, and the negative pole down.

6) Do not use this battery under the environment which more than 45°C , otherwise it will affect the performance and life of the battery.

7) The discarding method of the post-batteries should be deal with according to the urban environmental protection.

8. Safety Precautions

The manual said that the purpose of warning sign and legend is safe and proper using the product by yours, and prevents the harm to you and others.

The warning sign and legend as well as their meaning is as follows:

Warning sign	Meaning
 Contraindication	It shows that it will appear the dangerous of casualties or serious injury in the error use.
 Warning	It shows that it will appear the possibility of casualties or serious injury in the error use.
 Precautions	It shows that it will appear the possibility of personal injury or damaged goods in the error use.

Warning

- 1) The therapeutic apparatus cannot be used with the HF apparatus to avoid burns or damaged apparatus
- 2) If the patient himself uses the therapeutic apparatus and HF apparatus at the same time, the part of massage plate may be cause burns on apparatus, it may also damage the apparatus; If use the apparatus near (1 meter) the short wave or microwave therapeutic apparatus, the output of apparatus may be instability.
- 3) Do not modify this equipment without authorization of the manufacturer.
- 4) When need replacing the lithium batteries, please contact aftersales staff designated and authorized by the manufacturer directly, its accessories shall not be replaceable.
- 5) Product battery replacement should be replaced by professional maintenance staff, otherwise it will produce risk.
- 6) Shall not make personnel can no longer take care of themselves, infant, or not sensitive person use.
- 7) Simultaneous connection of a patient to a high frequency surgical ME equipment may result in burns at the site of the stimulator electrodes and possible damage to the stimulator.
- 8) Operation in close proximity (e.g. 1 m) to a shortwave or microwave therapy ME equipment may produce instability in the stimulator output.
- 9) Application of electrodes near the thorax may increase the risk of cardiac fibrillation.
- 10) Advice that stimulation should not be applied across or through the head, such as directly on the scalp, eyes, covering the mouth, on the front of the neck, (especially the carotidsinus).
- 11) It should be shut down then using again when the apparatus is moved or changed the therapeutic parts in the using process, otherwise there will be a strong stimulation.
- 12) It is not permitted to give children or people who were no ability to express consciousness to use.
- 13) If you feel unwell because of using the product, please stop using it immediately to the doctor consult.
- 14) Please pull up the power plugs when you are done or do not use the product.

15) Do not used in conjunction with other medical electronic devices, such as cardiac pacemaker, artificial heart and lung and other medical electronic devices with fuels life, electrocardiograph and other medical electronic devices, otherwise it will lead to danger.

17) Do not disassemble, repair and transform the therapeutic apparatus, otherwise cause failure or get an electric shock accident.

18) The therapeutic apparatus should be placed in the position which is easy to move the plugs to use, it is easy to move the plugs in emergency circumstances.

19) The long-term effects of chronic electrical stimulation are unknown, please do not use longer than 30 minutes every time.

20) Stimulation should not be applied over the carotid sinus nerves, particularly in patients with a known sensitivity to the carotidsinus reflex.

21) Stimulation should not be applied transthoracically in that the introduction of electrical current into the heart may cause cardiac arrhythmias.

22) Stimulation should not be applied transcerebrally.

23) Stimulation should not be applied over swollen, infected, or inflamed areas or skin eruptions, e.g., phlebitis, thrombophlebitis, varicose veins, etc.

Stimulation should not be applied over, or in proximity to, cancerous lesions.

Contraindication

1) Patients wired with cardiac pacemaker are prohibited to use!

2) Patients with haemorrhage after acute trauma or fractures or disabling wound after operation and those at healing stage after the operation should not use.

3) It is contraindicated for those who has a skin perceptual disturbance or cannot feel the heat.

4) Neither the patients who has the symptom of conscio-
usness disturbance onor children should not use it.

5) It is prohibited to use when bathing, sweating and sleeping.

6) The patient with cerebral hemorrhage which is in an unsteady phase should not use. The person who has sequela must be used under the supervision of doctors.

- 7) It is forbidden for those who has purulent inflammation, acute blood poisoning and continuous hyperpyrexia.
- 8) It is forbidden for those who has acute cardiovascular and cerebrovascular diseases to use.
- 9) Those who has medical electronic equipment implanted in their bodies, such as cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic devices Such could cause electric shock, burns, electrical interference, or even death.
- 10) Do not use on the open wounds, rashes, swollen, red, infected or inflamed areas or skin eruptions (e. g., phlebitis, thrombophlebitis, varicose veins); On top of, or in proximity to, cancerous lesions.
- 11) People with metal allergies are forbidden to use. Metal allergies are disabled.

Precautions

- 1) Safety of stimulators for use during pregnancy has not been established.
- 2) Caution should be used for patients with suspected or diagnosed heart problems.
- 3) Caution should be used for patients with suspected or diagnosed epilepsy.
- 4) Caution should be used in the presence of the following:
- 5) When there is a tendency to hemorrhage following acute trauma or fracture;
- 6) Following recent surgical procedures when muscle contraction may disrupt the healing process;
- 7) Over the menstruating or pregnant uterus; and
- 8) Over areas of the skin which lack normal sensation.
- 9) Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium. The irritation can usually be reduced by using an alternate conductive medium, or alternate electrode placement.
- 10) Electrode placement and stimulation settings should be based on the guidance of the prescribing practitioner.
- 11) The device should be kept out of the reach of children.

12) The device should be used only with the leads and remote control recommended for use by the manufacturer. 13) Don't use while driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user at undue risk of injury.	
Adverse	
Skin irritation and burns beneath the electrodes have been reported with the use of stimulators.	

9. Product maintenance

- 1) Before opening the power switch, please check the power supply and power plugs whether is in the normal working condition, each connecting parts whether is connected in good;
- 2) If this apparatus there is no evidence in the charging status, you should check whether there is a power supply, the power plugs are plugged in good, host toggle switch is in the "charging" position;
- 3) If the apparatus runs normally, show disorderly and button invalidly after charging. To turn it off after a few seconds and then opening, if still not normal, check the voltage whether is more than or less than the specified value;
- 4) If the remote control displays normally, but the host cannot run, please check the infrared launch window to remote control whether is corresponding with infrared launch window to therapeutic apparatus when operating the remote control, the apparatus will have a "DI" hint while the apparatus receives infrared signal successfully;
- 5) If all operation of the apparatus are normal but without output (the user without feelings) after charging. Please check the apparatus' electrode whether direct contacted with cervical skin, the contacting parts whether has hair, clothes and so on;
- 6) While cleaning, shutdown first, let the power line stay away from the electric supply socket. You can wipe by using soft cloth or towel with a little water and after dry.
- 7) After each use, please clean conducting plate, you can wipe by using soft cloth or towel with a little water and after dry, when conducting plate is particularly dirt, you can wipe by using soft cloth with a little medicinal alcohol (alcohol concentration as 75%).

- 8) The apparatus should be placed in the place which is dry, ventilated and good insulation use.
- 9) When the apparatus is moved, must be handled with care, avoiding the shake.
- 10) The apparatus should be checked the power plugs, to check whether the output is good, if has abnormal conditions, you should pay attention to it and service it in time.
- 11) The apparatus' shell board should be pay attention to protection, avoiding the abrasion.
- 12) The non-professional cannot disassemble the host to avoid getting an electric shock or damaging the apparatus, an accident can happen because of it, all at your peril.
- 13) The product has a free warranty for one year and provided by the dealer, and there is a compensable maintenance for the goods which expired for the warranty period.

14) Common faulted elimination:

Fault phenomenon	The items need to be affirmed	Handling method
The apparatus cannot boot	1. the host toggle switch whether is in"ON" position	Confirm whether the host toggle switch whether is in the"ON" position
	2. the electric quantity of inner battery whether is not enough	If the electric quantity is not enough please use the external power source
	3. the power line whether has been connected	Check the power line whether connected to the host;
	4. the power supply whether is normal	Check the power plugs supply whether is normal
	5. when operated the remote control, the infrared launch window to remote control whether is corresponding with the infrared launch window to therapeutic apparatus	The infrared launch window to remote control is corresponding with the infrared launch window to therapeutic apparatus

Cannot charging or without charging display	1. whether the adapter has connected good	Check the power line whether connected to the host
	2. the host toggle switch whether is in "charge" position	Confirm the host toggle switch whether is in "Charge" position
Display normally, but without response or feeling	1. the apparatus' electrode whether direct contacted to cervical skin	Check the apparatus' electrode whether direct contacted to cervical skin
	2. the parts contacted with electrode whether has hair, clothes and so on	Check the parts contacted with electrode whether has hair, clothes and so on

10. EMC Declaration

1. This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.

2. Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!

3. Caution: this machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

The use of ACCESSORIES, transducers and cables other than those specified, with the exception of transducers and cables sold by the MANUFACTURER of the AMBULATORY Neck Electronic Muscle Stimulator as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the Neck Electronic Muscle Stimulator.

Table 1

Guidance and manufacturer's declaration – electromagnetic emissions		
The MODEL P7 is compliance for each EMISSIONS test specified by the standard, e.g. EMISSIONS class and group.		
Emissions	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The P7 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Table 2

Guidance and manufacturer's declaration – electromagnetic immunity		
The P7 is compliance for each IMMUNITY test specified by the standard, e.g. IMMUNITY test level.		
Immunity test	IEC 60601-1-2 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Radiated RF EM fields IEC 61000-4-3	10V/m 80MHz-2.7GHz 80% AM at 1kHz	10V/m 80MHz-2.7GHz 80% AM at 1kHz
Electrical fast transient/burst IEC 61000-4-4	±2 kV 100 kHz repetition frequency	±2 kV 100 kHz repetition frequency
Surge IEC 61000-4-5	±0.5 kV, ±1 kV line-to-line; ±0.5 kV, ±1 kV and ±2 kV line-to-ground;	±0.5 kV, ±1 kV line-to-line; ±0.5 kV, ±1 kV, and ±2 kV line-to-ground;

Conducted disturbances induced by RF fields IEC 61000-4-6	3V 0.15MHz-80MHz 6V in ISM and amateur radio bands between 0.15MHz and 80MHz 80% AM at 1kHz	3V 0.15MHz-80MHz 6V in ISM and amateur radio bandsbetween 0.15MHz and 80MHz 80% AM at 1kHz
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT : 0.5 cycle a) At 0o, 45o, 90o, 135o, 180o, 225o, 270o, and 315o.	0% UT : 0.5 cycle a) At 0o, 45o, 90o, 135o, 180o, 225o, 270o, and 315o .
	0% UT : 1 cycle 70% UT : 25/30 cyclesb) Single phase: at 0o	0% UT : 1 cycle 70% UT : 25/30 cyclesb) Single phase: at 0o
	0% UT : 250/300 cyclesb)	0% UT : 250/300 cyclesb)
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz
NOTE a) UT is the a.c. mains voltage prior to application of the test level; b) E.g. 25/30 means 25 periods at 50 Hz or 30 periods at 60 Hz.		

Table 3 – Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment

Test Frequency (MHz)	Band (MHz)	Service a)	Modulation b)	Maximum Power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)	Compliance level
385	380-390	TETRA 400	Pulse Modulation 18 Hz	1.8	0.3	27	27
450	430-470	GMRS 460 FRS 460	FMc) ± 5 kHz deviation 1 kHz sine	2	0.3	28	28
710	704-787	LTE-Band 13, 17	Pulse Modulation 217 Hz	0.2	0.3	9	9
745							
780							
810	800-960	GSM 800/900, TETRA 800, Iden 820, CDMA 850, LTE Band 5	Pulse Modulation 18 Hz	2	0.3	28	28
870							
930							
1720	1700 -1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse Modulation 217 Hz	2	0.3	28	28
1845							
1970							
2450	2400 -2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse Modulation 217 Hz	2	0.3	28	28
5240	5100 -5800	WLAN 802.11 a/n	Pulse Modulation 217 Hz	0.2	0.3	9	9
5500							
5785							

NOTE:

- a) For some services, only the uplink frequencies are included.
- b) The carrier shall be modulated using a 50 % duty cycle square wave signal.
- c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

11. Product Scrapping Disposal

The therapeutic equipment belongs to medical device, if the complete machine is already aging and is not able to use, the process mode Haccording to the local legislation carry through scrap processing.

12. Accessories and Parts

Manual 1pcs

Warranty Card 1pc

Controller 1pc

Charging cable 1pc

Electrode Pads 2pc

Electrode Pads connection cable 1pc

13. The Paraphrase of Graphic Symbol

Graphic symbol	Paraphrase
	Product batch number tag
	Product serial number tag
	Manufacturer
	EU representative
	Production date
	Precautions
	Warning against danger
	Applied part of type BF
	Class II equipment
	This product only for indoor use
	It means prohibition (the thing was not permitted to do)
	It means enforcement (the thing must be observed)

	Keep out of the sun
	Keep away from rain
	Contact its local authorities to determine the proper method of disposal of potentially bio hazardous parts and accessories.
	CE label and notified body code
IP22	The dustproofing and waterproofing grade, that is, to prevent the solid object is greater than 12mm to invade, the inclination of 15 degrees can prevent dripping invasion still not cause harmful effects.
	This symbol is used to direct the user to refer to documentation for additional information regarding the system use or description.
	Please refer to the instructions for use

14. Executive Standards

The product conforms to the following standards and laws:

1. IEC 60601-1:2005+A1:2012 Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.
2. EN60601-1-11:2015 Medical electrical equipment General requirements for basic safety and essential performance- Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment.
3. EN60601-2-10 (2nd edition) Medical electrical equipment - Part 2-10: Particular requirements for the safety of nerve and musclestimulators.
4. IEC60601-1-2:2014 Medical electrical equipment — Part 1-2: General requirements for safety — Collateral standard: Electromagnetic compatibility — Requirements and tests.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out basket on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may lead to health damage or environmental pollution during further disposal or recycling of the equipment.

Shenzhen Pango Medical
Electronics Co., Ltd.
Main Site: No.25 1st Industry
Zone, Fenghuang Road, Xikeng
Village, Henggang Town,
Longgang District, Shenzhen,
Guangdong China

Lotus NL B.V.
Address: Koningin
Julianaplein 10, 1e
Verd, 2595AA, The
Hague, Netherlands.
Tel: +31644168999

We, the manufacturer, herewith declare under our sole responsibility that the above-mentioned products meet the provisions of the Regulation (EU) 2017/745 on Medical Devices (MDR). All supporting documentations are retained under the premises of the manufacturer. The full text of the declaration of conformity is available at the internet address: www.medivon.pl/deklaracje

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Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska