PL EN DE FR ES PT IT CZ SK SL EE LT LV
NB FI HU RO AR RU EL HR SV DA NL TR BG

BEAUTIFLY

SereneEye

PRO

series

Glasses with
red light

User manual

/beautifuly.global

/beautifuly_global

/@beautifuly.global

/@beautifuly.global

f1

f2

facebook.com/beautifly.global

youtube.com/@beautifly.global

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Material:	ABS + silicone
Charging method:	USB-C charging
Charging time:	approx. 24 hours
Battery capacity:	350 mAh
Voltage:	3.7 V
Output power:	1.2 W

Safety:

Read this user manual before using this product and follow the instructions contained therein when using the device.

Device description [page 2-f1]:

On/Off switch (1) Elastic band clasp (2) USB port (3).

Light sensitivity test:

Before using the device for the first time, check if you are allergic to light: hold it against your forearm for 3 minutes. If your skin becomes red or itchy, stop using it and consult a dermatologist. Always clean your face before use to avoid a reaction to cosmetics. If you experience a burning sensation, discontinue use immediately.

Usage:

1. Cleanse your face and pat it dry.
2. Fasten and adjust the elastic headband.
3. Press the power button to turn on the device.
4. Press the mode button to select the desired mode (light therapy, microvibration, light therapy + microvibration).
6. After three minutes, the device will automatically turn off. Remove the mask and clean it.
7. Apply serum or moisturizer immediately after treatment.

Charging Method:

The product should not be used while charging. (Note: the device will not function and the buttons will not respond while charging.) Before first use, it is recommended to fully charge the device. Connect the USB-C cable to the device, then connect the USB plug to the adapter. Once the power is on, charging will begin. The device is fully charged when the indicator light changes from red to green. A fully charged device will last for approximately 30 minutes of use and 10 treatments.

Device Care and Maintenance:

Store this product in a dry place at room temperature and avoid exposure to dust, direct sunlight, or humid conditions. Clean the mask before and after use.

Precautions:

Keep this product away from small children and animals to prevent it from being treated as a toy. Handle the product with care during use. This device is designed as a household electronic product and is not suitable for diagnosis or treatment. It should not be used by pregnant or breastfeeding women, individuals who are photosensitive or taking photosensitizing medications, or those with genetic skin conditions or medical conditions such as lupus, photosensitive eczema, or albinism. Do not use it while charging. Do not use the device if it is damaged. Do not repair, disassemble, or modify the device yourself. Any abnormalities should be assessed by a service center. Do not drop or hit the device to avoid malfunctions. Wipe the product with a soft, damp cloth after each use. Do not wash the device under water or immerse it. Do not use caustic agents to clean the device. Store it in a cool, dry place.

More information about the company and products [page 2 - f2]:

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[EN] Rules for handling electrical and electronic equipment: Waste electrical and electronic equipment as well as batteries and accumulators cannot be mixed with other waste. This is indicated by the crossed out bin symbol on the product packaging. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage may result in damage to health or environmental pollution. Proper household handling of waste electrical and electronic equipment and batteries and accumulators positively contributes to the reuse and recovery, including recycling, of waste equipment. Waste electrical and electronic equipment as well as batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the used equipment is of the same type and performed the same functions as the equipment sold. When delivering equipment intended for households to the buyer, the distributor is obliged to collect used equipment from households free of charge at the place of delivery of the equipment, provided that the used equipment is of the same type and performed the same functions as the delivered equipment. A distributor running a retail unit with a sales area within the meaning of Art. 2 point 19 of the Act of March 27, 2003 on spatial planning and development of at least 400 m² devoted to the sale of equipment intended for households, is obliged to accept waste equipment from households free of charge in this unit or in its immediate vicinity, none of the external dimensions exceeds 25 cm, without the need to purchase new household equipment. The product was introduced to the market on August 13, 2005

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[EN] Summary declaration of conformity: The manufacturer declares, under its sole responsibility, that the SereneEye Pro product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Directive Commission delegated act (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of substances subject to restriction and Directive (EU) 2017/2102 of the European Parliament and of the Council of November 15, 2017 amending Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment. The full content of the declaration of conformity is available at www.beautifly.eu/deklaracje

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