BERDSEN

BD-101

INFRARED
**FOREHEAD
THERMOMETER**

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EN

Instructions manual

Read and follow the instructions before installation and use. Keep the instructions for future reference.

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


www.berdsen.pl

CE 0123

Please read these instructions carefully and thoroughly before using this product.

The BD-101 infrared forehead thermometer has been specially designed for safe use on the forehead. It is a device capable of measuring body temperature by detecting the intensity of infrared light emitted from the forehead area. It works by converting the recorded heat into a temperature reading, which is then displayed on an LCD screen. With proper operation, the thermometer will quickly and accurately assess the user's body temperature. It is a practical device that enables accurate temperature monitoring safely and conveniently. The device is dedicated to the non-invasive measurement of human body temperature, carried out through the skin of the forehead. It is ideal for both home and clinical use, enabling measurements for people of all ages. Thanks to its safe and convenient application on the forehead area, this thermometer provides an effective tool for monitoring body temperature in a variety of environments.

Indications for use: the infrared forehead thermometer is used to measure body temperature by measuring the temperature on the forehead.

This device is compatible with the following standards:

- ISO 80601-2-56 Medical electrical equipment - Part 2-56: Detailed requirements for basic safety and essential performance of clinical thermometers for measuring body temperature,
- IEC 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical

electrical systems used in the home healthcare environment and complies with the requirements of IEC 60601-1-2(EMC), IEC 60601-1(Safety). The manufacturer is ISO 13485 certified. The product complies with the requirements of MDD Directive 93/42/EEC, "0123" is the identification number of the notifying body.

WARNING

- Use of the infrared thermometer is not restricted by gender or age.
- Do not touch the temperature probe with your hands.
- Use of this forehead thermometer is not a substitute for consultation with a doctor.
- Do not allow children to take the temperature unsupervised, as some parts are small enough to be swallowed.
- Never immerse the device in water or other liquids.
- Do not modify the device without the manufacturer's approval.
- Do not expose the thermometer to extreme temperatures (below -25°C/-13°F or above 55°C/131°F) or excessive humidity (>95% relative humidity).
- Keep the battery away from children.
- Remove the battery if the device is not used for an extended period of time.
- Do not place the thermometer in direct sunlight or with cotton wool, otherwise the accuracy of the measurement will be affected.
- Portable and mobile RF (Radio Frequency) communications may affect the performance of this device. It is therefore necessary to take special precautions and follow the recommendations regarding electromagnetic compatibility (EMC). Detailed information on

the application of EMC precautions can be found at the end of the manual included in the device documentation.

- ME equipment should not be cleaned or wiped down during use.
- Any serious incident related to the equipment should be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is located.
- The ME equipment probe should not be serviced or maintained during use with the patient.

DEVICE

Design of the device

1. Probe
2. Test button
3. Sound switch
4. Settings button
5. Battery cover
6. Display

Icons on the display

- a. Sound sign
- b. Battery sign
- c. Unit sign: °C or °F
- d. AM and PM sign in time mode
- e. Time and date
- f. Temperature result
- g. Memory mode
- h. Front temperature mode
- i. Object mode

BASIC FUNCTIONS

CLOCK

The real-time clock will function in memory and help support the measurement results.

Refer to the Clock Settings section to learn how to configure the time when using it for the first time.

FOREHEAD MODE

The thermometer has been designed for practical use. It is not a substitute for visiting your doctor. You should also remember to compare the measurement result with your regular body temperature.

To learn how to measure your body temperature, refer to the Features section.

OBJECT MODE

The object mode shows the actual uncorrected surface temperature, which differs from the body temperature. This can help monitor whether the object temperature is appropriate for a baby or patient, such as baby milk.

Measurement range of object mode: 0°C~100°C (32°F~212°F).

Laboratory accuracy of object mode: $\pm 4\%$ or $\pm 2^\circ\text{C}$ (4°F), whichever is greater.

To learn how to measure the object temperature, refer to the Features section.

MEMORY MODE

There are 30 memory sets available for forehead and object measurements. Each memory also stores the measurement date/time/mode icon.

BEEP

One short beep will sound when testing is complete.

°C/ °F SWITCH

To learn how to switch between degrees Celsius and Fahrenheit, refer to the Clock Settings section.

SWITCHING THE SOUND ON/OFF


Use the button to select whether the sound should be switched on or off.

ADVANTAGES OF THE THERMOMETER

The infrared forehead thermometer accurately measures body temperature, focusing on the area of key organs (see Figure 1). Its use involves measuring the temperature of the skin surface over the temporal artery, which is the main artery of the head.

The temporal artery, which is a direct extension of the aorta - the main trunk of the arterial system - connects to the heart via the carotid artery. It provides a constant flow of blood, making it a particularly representative site for overall body temperature. Therefore, changes in body temperature are reflected more rapidly in the forehead than in other areas such as the mouth, rectum or armpit.

Vital organs


SETTING THE CLOCK

When using the thermometer for the first time, set the device parameters. With the thermometer off, press and hold the **settings button** to enter the settings mode.

1. Setting the unit

Press the **test button** to select the desired unit.


Once the unit is set, press the **settings button**, the time format will be displayed.

2. Setting the time format

The unit can display the time in 12-hour (AM/PM) or 24-hour format.

Press and release the **test button** to select the format. Once the preferred time format has been selected, press the **settings button**, the hour digit will start flashing automatically.

3. Setting the time

Press and release the **test button** to advance one hour until the correct time is displayed.

Once the hour is set, press the **settings button** and the minute value will start flashing automatically.

4. Setting the minute

Press and release the **test button** to advance one minute until the correct minute appears.

Once the minute is set, press the **settings button** and the year digit will start flashing automatically.

5. Setting the year

Press and release the **test button** to advance one year until the correct year appears.

After setting the year, press the **setting button**, the month digit will appear.

6. Setting the month

Press and release the **test button** to advance one month until the correct month appears.

After setting the month, press the **settings button** to display the Date digit.

7. Setting the date

Press and release the **test button** to advance one day until the correct month appears.

Once the day has been set, press the **setting button** to exit the setting mode.

TEMPERATURE MEASUREMENT TIPS

To ensure an accurate reading of body temperature with an infrared forehead thermometer, it is useful to consider the following factors that can affect the accuracy of the measurement:

1. Knowing each person's normal body temperature when they feel well is crucial for an accurate diagnosis of fever. To determine this baseline temperature, it is advisable to take several measurements when the patient is well. To confirm the readings, it is useful to take an additional measurement with a standard digital thermometer.
2. Users must be in the room for 30 minutes before taking a measurement.
Note: Users and the thermometer should be at the same ambient temperature for at least 10 minutes before taking a reading.
3. Users should not drink, eat or be physically active, such as bathing, showering, shampooing and drying their hair before/during the measurement. Remove cap and wait 10 minutes before taking a reading.
4. Oils or cosmetics on the forehead may give a lower temperature reading than the actual reading. Remove dirt from the forehead before taking the measurement. Wait at least 10 minutes after washing the forehead before taking a measurement.
5. Keeping your hand on your forehead for a long time will affect the temperature reading.
6. Do not take the temperature on scars, open wounds or abrasions.
7. Do not use the thermometer on a sweaty or perspiring forehead as this may affect the reading.

8. Do not measure during or immediately after feeding your baby.
9. Do not measure the temperature with this thermometer near areas that are very hot, such as fireplaces and cookers.
10. The probe window of the thermometer is the most delicate part of the device. Do not touch the probe window. Damage or dirt to the probe window can affect the accuracy of the reading.
11. If the thermometer is stored in an environment significantly different from the test site, place it in the test site for approximately 30 minutes before use.
12. It is not intended for use in oxygen-rich environments or in the presence of flammable anaesthetic mixtures with air, oxygen or nitrous oxide.

FUNCTIONS

FOREHEAD TEMPERATURE MEASUREMENT

1. Press the test button, the display will be activated to show all segments. After the self-test, the display will show the icons as in figure 2 along with beeps to start a new measurement.

2. Point the thermometer at the centre of the forehead at a distance of less than 5 cm (see figure 3) and then press the test button.
Note: Do not remove the thermometer from the forehead before the beep sounds.
3. Read the temperature on the display.
4. The unit will automatically switch off if left idle for 30 seconds.

TURN THE SOUND ON/OFF


You can flick the sound switch to turn the sound on or off.

SWITCHING FRONT AND OBJECT MODE

Press and hold the settings button to switch between forehead mode and object mode.

MEASURING OBJECT TEMPERATURE

1. Press the test button to switch the thermometer on, and once in measurement mode you can measure the temperature of an object (see figure 4).
2. Point the thermometer at the centre of the object whose temperature you want to measure, from a distance of less than 5 cm.
3. Press the test button and then read the temperature on the display.
4. The device will automatically switch off if left idle for 30 seconds.

AFTER MEASUREMENT

1. Switch off the power. After 30 seconds of inactivity, the device will automatically enter power saving mode and switch off. This is a feature designed to extend battery life.
2. Clean the probe after each use to ensure accurate readings and avoid cross contamination.
(See Care and cleaning section for details).

BACKLIGHTING

In forehead mode:


1. The display will illuminate GREEN for 3 seconds when the unit is ready to measure, and the display will illuminate GREEN for 3 seconds when the measurement is completed with a reading below 37.3 °C (99.1 °F).
2. The display will light up YELLOW for 3 seconds when the measurement is completed with a reading below 37.8 °C (100.0 °F).
3. The display will light RED for 3 seconds after the measurement is completed with a reading equal to or higher than 37.8 °C (100.0 °F).

In object mode:

The display will illuminate GREEN for 3 seconds when the instrument is ready for measurement and the measurement is complete.

MEMORY MODE

1. The memory mode can be accessed in the forehead mode or in the object mode:
After switching on the thermometer and following the steps shown in Figures 2-4 or completing the test, press the settings button. The letter M will appear in the lower left corner of the display (see figure 5).
2. The thermometer will automatically store the last 30 temperature readings. A measurement mode icon is also stored in each memory. Each time the settings button is pressed, the screen displays the previous readings corresponding to numbers 1 to 30. The number 1 reflects the most recent reading, while the number 30 indicates the oldest reading stored in memory (see Figure 6).
3. In memory mode, the "▲" or "▼" marker will not change. The user can press the test button to take new measurements.



CARE AND CLEANING

1. To ensure accurate readings, the probe window must always be kept clean, dry and free of damage. The accuracy of temperature measurements can be affected by the presence of dirt, fingerprints, dust and other debris on the surface of the probe window. Damage to this area can significantly affect the performance of the instrument, leading to erroneous readings or other problems.
2. The probe window is the most delicate part of the thermometer. Use a soft cloth lightly moistened with a 75% isopropyl alcohol solution to disinfect the probe window and thermometer. Do not use abrasive cleaners. After thorough cleaning, it is recommended to wait at least 10 minutes for it to dry before taking a temperature measurement.
Note: Do not use chemicals other than isopropyl alcohol to clean the probe window.
3. Use a soft, dry cloth to clean the thermometer display and its exterior.
4. Do not put the thermometer directly into water.
5. Store the thermometer in a dry place, free from dust and dirt and away from direct sunlight.
6. Place the thermometer back in its original packaging after use.

BATTERY REPLACEMENT

1. Replace the battery when the "  " appears in the upper right corner of the LCD display.
2. Slide the battery cover down as shown in Figure 7.
3. Remove the battery and install 2 new AAA alkaline batteries as shown in Figure 8.
4. Slide the battery cover back.

SPECIFICATIONS

Measuring range	Front mode: 34.0°C~43.0°C (93.2°F~109.4°F)
Place of measurement	Forehead (Forehead mode)
Manufacturer model	DET-3017 (Ref. BD-100)
Reference site	Oral (This thermometer converts the reference forehead temperature to display its "oral equivalent").
Operating mode	Forehead mode (adjustment mode)
Laboratory accuracy	Front mode: $\pm 0.2^{\circ}\text{C}$ (0.4°F) during 35.5°C ~ 42.0°C (95.9°F ~ 107.6°F) over the operating temperature range 15°C ~ 35°C (59.0°F ~ 95.0°F) $\pm 0.3^{\circ}\text{C}$ (0.5°F) for other measurements and operating temperature range
Display resolution	0.1°C or 0.1°F
Measurement time	Approximately one second
Operating temperature range	5°C ~ 40°C (41°F ~ 104°F), 15%~85% relative humidity, non-condensing. Atmospheric pressure: 70 kPa ~ 106 kPa
Storage and transport temperature range	-25°C ~ 55°C (-13°F ~ 131°F), 15% ~ 95% relative humidity, non-condensing Atmospheric pressure: 70 kPa ~ 106 kPa
Clinical accuracy	0-1 year: Clinical deviation: -0.19°C (-0.34°F) ; Clinical repeatability: 0.06°C (0.11°F) ; Compliance limits: 0.68°C (1.22°F); 1- 5 years: Clinical bias: -0.25°C (-0.45°F) ; Clinical repeatability: 0.07°C (0.13°F) ; Compliance limits: 0.73°C (1.31°F) Over 5 years: Clinical deviation: -0.25°C (-0.45°F) ; Clinical repeatability: 0.17°C (0.31°F) ; Compliance limits: 0.60°C (1.08°F)
Shock resistance	to a fall from a height of 3 feet
Dimension	155*39.2*49.4mm
Weight	approx. 97 grams (with batteries)
Battery	DC3V (2×AAA battery)
Battery life	Approximately 6000 readings
Part used	Probe
Expected service life	Three years
Degree of protection	IP22
Accessories	Battery
Contraindications	Non

TROUBLESHOOTING

Komunikat	Problem	Solution
	The thermometer is not working properly.	Discharge the battery, wait 1 minute and power it up again. If the message reappears, contact your dealer for assistance.
	The ambient temperature is not between 5°C and 40°C (41°F–104°F).	Place the thermometer indoors for at least 30 minutes at a room temperature between 5°C and 40°C (41°F–104°F).
	In forehead mode: The measured temperature is higher than 43.0 °C (109.4 °F). In object mode: the measured temperature is higher than 100 °C (212 °F).	Carefully read the temperature measurement instructions, then take a new temperature measurement.
	In front mode: The measured temperature is lower than 34.0 °C (93.2 °F). In object mode: the measured temperature is lower than 0 °C (32 °F).	Carefully read the temperature measurement instructions, then make sure the lens filter is clean and take a new temperature measurement.
	The thermometer is functioning properly.	The thermometer can be used normally
	When the battery indicator flashes, it means that the battery level is low, but you can continue measuring.	The thermometer will take a correct measurement, but the batteries should be replaced soon.
	The thermometer cannot operate due to low battery level.	Replace with two new AAA size alkaline batteries.

SYMBOL EXPLANATION

	Manufacturer
	General symbol for recovery/ recyclable
	Authorized representative in the European Community
	Date Of manufacture
	Batch Code
	Unique device identifier
	Storage and Transportation Temperature Limit:
	TYPE BF APPLIED PART
	Refer to instruction manual/ booklet
	Atmospheric pressure limitation
	Storage and Transportation Humidity limitation:
	Disposal of this product and used batteries should be carried out in accordance with the national regulations for the disposal of electronic products.
IP22	The first num.2: Protected against solid foreign objects of 12,5mm Ø and greater. The second num.2: Protection against vertically falling water drops when ENCLOSURE tilted upto 150
 0123	The product conforms to the requirements of the EC Directive MDR(EU) 2017/745 on medical devices.
	Medical Device
	Caution

INFORMATION FOR USERS ON HOW TO DEAL WITH WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT.

If the device, packaging, instruction manual and the like are marked with the symbol of a crossed-out wheeled waste container, it means that the device is subject to separate waste collection in accordance with Directive 2012/19/EU of the European Parliament and of the Council. The used device should not be disposed of with other household waste and should not be treated as municipal waste. By disposing of electro-waste in the container you pose a threat to the environment. Take the used device to a separate collection point for electrical and electronic equipment organized by the public administration. By segregating and submitting used electrical and electronic equipment for processing, recovery, recycling and disposal, you protect the environment from pollution and contamination, contribute to reducing the use of natural resources and lowering the cost of producing new equipment. Proper disposal and scrapping helps eliminate the adverse effects of scrapped products on the environment and human health. For details on the recyclability of this product, please contact your local municipality, city cleaning services or the store where you purchased the product.

INFORMATION ON ELECTROMAGNETIC COMPATIBILITY

The device complies with the requirements of the international standard IEC 60601-1-2 for electromagnetic compatibility. The requirements are met under the conditions described in the table below. The device is an electrical medical product and is subject to special precautions with regard to EMC, which must be published in the operating instructions. Portable and mobile HF communication equipment can affect the device. Using the device in combination with unapproved accessories may adversely affect the device and alter electromagnetic compatibility. The device should not be used directly near or between other electrical equipment.

Table 1

Guidelines and manufacturer's declaration - electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or user of the device should ensure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - guidance
Conducted emission CISPR 11	Not applicable	The device uses RF energy exclusively for its internal functions. Therefore, RF emissions are very low and should not cause any interference to nearby electronic equipment.
Radiated emission CISPR 11	Group 1 Class B	The device uses RF energy exclusively for its internal functions. Therefore, the RF emissions are very low and should not cause any interference to nearby electronic equipment.
Acoustic emissions IEC 61000-3-2	N/A	
Voltage fluctuation / Flicker emissions IEC 61000-3-3	N/A	

Table 2

Guidance and manufacturer's declaration - electromagnetic immunity		
The device is intended for use in the electromagnetic environment specified below. The customer or user of the equipment should ensure that it is used in such an environment.		
Immunity test	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4kV, ±8 kV air	±8 kV contact ±2 kV, ±4kV, ±8 kV, ±15 kV air
Electrostatic transients/ impacts IEC 61000-4-4	±2 kV for supply lines Repetition rate 100 kHz ±1kV for input/output lines	N/A
Overvoltages IEC 61000-4-5	±0,5kV, ±1kV in line-to-line differential mode	N/A
Voltage drops, short interrup- tions and voltage fluctuations on power input lines IEC 61000-4-11	0% UT (100% UT drop) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT (100% UT decrease) for 1 cycle at 0° 70% UT (30% drop in UT) for 25/30 cycles at 0° 0% UT (100% UT drop) for 250/300 cycles at 0°	N/A
Power frequency magnetic field (50/60 Hz) IEC 61000-4-8	30 A/m, 50/60Hz	30 A/m, 50/60Hz
Proximity magnetic fields	IEC 61000-4-39	See Table 3
NOTE: UT is the mains voltage before the test level is applied.		

Table 3

ENCLOSURE PORT IMMUNITY test specifications for proximity magnetic fields		
Test frequency	Modulation	IMMUNITY TEST LEVEL (A/m)
30 kHz ^{a)}	CW	8
134,2 kHz	Pulse modulation ^{b)} 2,1 kHz	65 ^{c)}
13,56 MHz	Pulse modulation ^{b)} 50 kHz	7,5 ^{c)}
a) This test is only applicable to MEDICAL DEVICES AND SYSTEMS intended for use in HOME CARE ENVIRONMENTS. (b) The carrier shall be modulated using a 50% duty cycle square wave signal. (c) r.m.s. before modulation is applied.		

Table 4

Guidance and manufacturer's declaration - electromagnetic immunity		
The device is intended for use in the electromagnetic environment specified below. The customer or user of the equipment should ensure that it is used in such an environment.		
Immunity test	IEC 60601 test level	Compliance level
RF conducted IEC 61000-4-6	3 Vrms 150 kHz up to 80 MHz 6 Vrms 150 kHz to 80 MHz outside the ISM bands	N/A
RF radiated IEC 61000-4-3	10 V/m 80 MHz to 2,7 GHz	10 V/m

NOTE 1 At 80 MHz and 800 MHz the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. The propagation of electromagnetic waves is affected by absorption and reflection from structures, objects and people.

a The ISM (industrial, scientific and medical) bands between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. Amateur radio bands from 0.15 MHz to 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 1.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.

b The compliance levels in the ISM frequency bands 150 kHz to 80 MHz and in the frequency range 80 MHz to 2.7 GHz are intended to reduce the likelihood that mobile/portable communication equipment could cause interference if accidentally brought into patient areas. For this reason, an additional factor of 10/3 has been included in the formulae used to calculate the recommended distance for transmitters in these frequency bands.

c Field strengths from fixed transmitters, such as mobile/wireless phone base stations and fixed mobile radios, amateur radio stations, AM and FM radio stations and TV stations, cannot be accurately predicted theoretically. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic survey of the site should be considered. If the measured field strength at the site where the device is used exceeds the relevant RF compliance level above, the device should be observed to check normal operation. If abnormal operation is observed, additional measures such as changing the orientation or location of the device may be necessary.

d In the frequency range 150 kHz to 80 MHz, the field strength should be less than 3 V/m.

Table 5

Recommended distances between RF wireless communication devices				
The device is intended for use in an electromagnetic environment where RF interference is controlled. The customer or user of the device can help prevent electromagnetic interference by maintaining a minimum distance between the RF wireless communication equipment and the device, according to the following recommendations, depending on the maximum output power of the communication equipment.				
Frequency MHz	Maximum power W	Distance	Test level IEC 60601	Compliance level
385	1.8	0.3	27	27
450	2	0.3	28	28
710	0.2	0.3	9	9
745				
780				
810	2	0.3	28	28
870				
930				
1720	2	0.3	28	28
1845				
1970				
2450	2	0.3	28	28
5240	0.2	0.3	9	9
5500				
5785				
Note 1: These guidelines may not apply in all situations. The propagation of electromagnetic waves is affected by absorption and reflection from structures, objects and people.				

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JOY TECH Healthcare Co., Ltd.
**EC REP Shanghai International Holding Corp. GmbH (Europe)