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GB Instruction Manual
Personal Scale PS-72E

- 1 Battery Compartment (on underside)
- 2 Unit key (on underside) for kilograms (kg), pounds (lb) stones (st)
- 3 Display
- 4 Standing area

Explanation of symbols

IMPORTANT
Non-observance of these instructions can result in serious injury or damage to the device.

WARNING
These warning notes must be observed to prevent any injury to the user.

CAUTION
These notes must be observed to prevent any damage to the device.

NOTE
These notes give you useful additional information on the installation and operation.

LOT LOT number

Manufacturer

Recycling symbols/codes:
These are used to provide information about the material and its proper use and recycling.

GB Safety instructions

Read the instruction manual carefully before using this device, especially the safety instructions and keep the instruction manual for future use. Should you give this device to another person, it is vital that you also pass on these instructions for use.

- The device is intended for indoor private use only and should not be employed for commercial purposes.
- Only use the device for the according the instructions for use. The guarantee is invalidated if the device is misused.
- Never place the scales on an uneven surface. Never stand on just one corner of the scales. Do not stand on one side or on the edge of the scales when weighing.
- Do not stand on the scales with wet feet. Do not stand on the scales wearing socks.
- The scale is designed for domestic use. It is not intended for commercial use in hospitals or other medical facilities..
- The scales have a measuring range up to 150 kg (330 lb / 24 st). Do not overload the scales.
- Handle the scales with care. The scales must not be subjected to impacts or vibration. Do not drop it. Keep the scales away from water.
- Do not attempt to repair the device yourself in the event of a malfunction since this will invalidate the warranty. Repairs should only be carried out by authorised service centres.
- If the device is not going to be used for a long period, please remove the battery.
- Children must not play with the appliance.

BATTERY SAFETY INSTRUCTIONS

- Do not disassemble the battery!
- Replace the battery when the battery symbol "Lo" appears on the display.
- Remove weak batteries from the battery compartment immediately, as they can leak and damage the appliance!
- Increased risk of leakage, avoid contact with skin, eyes and mucous membranes! In the event of contact with battery acid, rinse the affected areas immediately with plenty of clean water and consult a doctor immediately!
- If a battery has been swallowed, consult a doctor immediately!
- Insert the battery correctly, observe the polarity!
- Keep batteries away from children!
- Do not recharge, short-circuit or dispose of batteries in fire! There is a risk of explosion!
- Do not dispose of used batteries and rechargeable batteries in household waste, but in hazardous waste or at a specialised battery collection point!

Scope of delivery

Please check first of all that the device is complete and is not damaged in any way. In case of doubt, do not use the appliance and contact your dealer or your service centre.

The following parts are included:

- 1 **ecomед** personal scale **PS-72E**
- 1 Battery Lithium button cell 3V type CR2032
- 1 Instruction manual

Packaging is recyclable or can be returned to the raw material cycle. Please dispose of packaging material that is no longer required properly. If you notice any transport damage when unpacking, please contact your dealer immediately.

WARNING
Make sure that the packaging films do not get into the hands of children! There is a danger of suffocation!

Insert / replace battery

Before you can operate the scales, you must remove the insulating strip from the battery compartment **1** on the underside of the appliance. Replace the battery when the battery replacement symbol "Lo" appears in the display **3** or if nothing is shown in the display after you have entered the appliance. Pay attention to the polarity (positive pole upwards; pointing towards the cover).

Change the unit of weight

You can display your weight in kilograms (kg), stones (st) or pounds (lb). Press unit key **2** until the desired unit is displayed.

Weighing ("step-on" function)

1. Make sure the scales are on a firm, level surface. If you only want to weigh yourself but have moved the scales just before, you must first initialise the scales. To do this, press your foot on the middle of the surface of the scales. "0.0 kg" appears in the display **3**. **When the scale has switched off, it is ready for the "Step-on" function.** This procedure is not necessary if you have not moved the scales.
2. Step onto the scales and stand still. The device switches on automatically and the display **3** briefly shows "0.0 kg".
3. Your weight is measured, flashes twice and is then finally displayed.
4. Step off the scales. The display changes to "0.0 kg". The scales switch off automatically after approx. 10 seconds.

Error messages

Err = The scales are overloaded.

Lo = The battery is weak and needs to be replaced.

If your personal scales do not function to your satisfaction, check the following before contacting the service centre:

- Check that you have inserted the battery correctly.
- Check that you have selected the correct weight unit.
- Check that the scales are standing freely on a firm and level surface. It must not touch a wall or other object.
- Repeat the weighing process.

Cleaning and care

Never use aggressive cleaning agents or strong brushes. Only clean the bathroom scales with a soft, slightly damp cloth. Never use harsh cleaning agents or alcohol. Do not allow water to penetrate the appliance. Only use the appliance again when it is completely dry.

Technical specifications

Name and model:	ecomед personal scale PS-72E
Power supply:	3V==, 1 x 3V Lithium button cell type CR2032
Measuring range:	up to 150kg, 330lb, 24st
Graduation:	0,1kg, 0,2lb, 1lb
Autom. switch-off:	after approx. 10 seconds
Dimensions (LxWxH):	approx. 30 x 30 x 2 cm
Weight:	approx. 1.5 kg
Art. number:	23511
EAN number:	40 15588 23511 5

Disposal

This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Consult your municipal authority or your dealer for information about disposal.

Warranty and repair terms

Your statutory warranty rights are not restricted by our guarantee below. Please contact your dealer or the service centre in case of a claim under the warranty. If you have to return the unit, please enclose a copy of your receipt and state what the defect is.

The following warranty terms apply:

1. The warranty period for **ecomед** products is 2 years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Defects in material or workmanship will be removed free of charge within the warranty period.
3. Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
4. The following is excluded under the warranty:
 - a. All damage which has arisen due to improper treatment, e.g. non-observance of the user instructions.
 - b. All damage which is due to repairs or tampering by the customer or unauthorised third parties.
 - c. Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
 - d. Accessories which are subject to normal wear and tear.
5. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

medisana GmbH, Carl-Schurz-Str. 2, 41460 NEUSS, GERMANY

In accordance with our policy of continual product improvement, we reserve the right to make technical and visual changes without notice.

The current version of this instruction manual can be found under <https://docs.medisana.com/23511>.

Information about service can be found here:
www.medisana.com/servicepartners

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