**medivon.**

## **User manual**

PL | EN | DE | FR | ES | PT | IT | CZ | SK  
SL | EE | LT | LV | HU | RO | RU | NB | FI  
EL | HR | SV | DA | NL | TR | BG | AR

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The Medivon Seraph massager is a device designed to improve microcirculation, reduce cellulite, promote muscle relaxation, and care for the skin. Massage should be performed according to the instructions, avoiding prolonged pressure on one spot. People with heart disease, thrombosis, cancer, skin infections, osteoporosis, fever, as well as pregnant or menstruating women, should not use the device without consulting a doctor.

The device is equipped with four interchangeable attachments: a flat one for smoothing the skin, a ring one for wave massage, a wave one for breaking down fat, and an anti-cellulite head. Before beginning, select the appropriate attachment and attach it to the bottom of the device.

To turn on the massager, press and hold the power button for approximately 2 seconds. Short presses allow you to change the intensity modes from lowest to highest (4 levels in total). Massage should be performed by moving the device over the skin, without applying pressure – the pressure applied depends on the area being massaged and should not cause discomfort.

The Medivon Seraph is charged via a USB Type-C cable. When the device's LED flashes red, it means the battery is below 10% and requires charging. Do not use the device while charging. A full charge takes approximately 3-4 hours. To protect the battery, do not leave the device discharged for extended periods.

After use, wipe the device with a dry cloth – do not use water or detergents. Do not immerse it in water or use it in humid environments. After use, turn off the device and store it in a dry place.

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Used electrical and electronic equipment and batteries and accumulators must not be mixed with other waste. This is indicated by the crossed-out bin symbol on the product packaging. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution. Proper handling of used electrical and electronic equipment and batteries and accumulators in the household contributes positively to the reuse and recovery, including recycling, of used equipment. Used electrical and electronic equipment and batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the used equipment is of the same type and has performed the same functions as the equipment sold. A distributor, when delivering household equipment to a purchaser, is obliged to take back used household equipment free of charge at the place of delivery of such equipment, provided that the used equipment is of the same type and has performed the same functions as the equipment delivered. A distributor running a retail outlet with a sales area within the meaning of art. 2 item 19 of the Act of 27 March 2003 on spatial planning and development of at least 400 m<sup>2</sup> devoted to the sale of household equipment, is obliged to accept, free of charge, in that outlet or in its immediate vicinity, used household equipment, none of the external dimensions of which exceeds 25 cm, without the need to purchase new household equipment. The product was introduced to the market after 13 August 2005

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The manufacturer declares under its sole responsibility that the Medivon Seraph product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances. The full text of the declaration of conformity is available on the website [www.medivon.pl/deklaracje](http://www.medivon.pl/deklaracje)

Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska  
Model: Medivon Seraph | SN:07255904119286599 | Marka: Medivon | email: [serwis@medivon.pl](mailto:serwis@medivon.pl)

Producent: **Xiamen Healthpal Electronic Co., Ltd.**  
NO.148 BINGLANGXILI, SIMINGDISTRICT. XIAMEN,FUJIAN, CHINA  
email: [vilin@xmhealthcare.com](mailto:vilin@xmhealthcare.com)

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#### Shortened warranty terms

This Guarantee is provided by Armare S.A., a company governed by Polish law, with its registered office in Warsaw, ul. Grzybowska 87, 00-844, Poland. You can find us in the official Polish register of companies under the KRS number: 0000730590, NIP: 5213833614 and REGON: 380721860. We guarantee that the equipment you have purchased from us (hereinafter referred to as the "Product") will work properly. If something goes wrong for a reason inherent in the Product, this warranty protects you. The warranty starts from the day you receive the Product and lasts: 24 months, if you are a consumer or a natural person running a business and the purchase of the Product is not of a professional nature for you. Otherwise, the warranty period is 12 months. If the law in your area requires a longer warranty, the warranty period will be extended to meet local law requirements.

In addition, if the Product does not comply with the contract, you have the right to legal protection from and at the expense of the seller. This warranty does not affect these rights - it does not limit, suspend or exclude them.

Detailed warranty conditions, how to submit complaints, what happens when you find a defect, what are exceptions to the warranty, when the warranty expires or is extended, and what are the rules regarding the processing of personal data, you can find on our website: [www.medivon.pl/warranty](http://www.medivon.pl/warranty)