medivon®

User manual

PL | EN | DE | FR | ES | PT | IT | CZ | SK
SL | EE | LT | LV | HU | RO | RU | NB | FI
EL | HR | SV | DA | NL | TR | BG | AR

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Device manual

The device is ready to use right out of the box. Make sure the massager is fully charged – a full charge takes approximately 3 hours. Then, apply the desired amount of oil to your skin, select one of the 6 interchangeable attachments, and turn the device on. To turn on the massager, hold the power button for a few seconds – you will hear a beep, and the display will show "1" or "2," indicating the selected mode.

The HEAT button (flame icon) activates the heating function. Each press increases the temperature, up to a maximum of 55°C. The "0" setting means the heating is off. The Suction button adjusts the vacuum massage intensity – you can choose from 9 suction levels and an intermittent mode, which provides temporary relief and comfort during the treatment.

You can perform the treatment on different body parts simultaneously – thanks to the 6 interchangeable attachments, you can attach multiple cups at the same time, similar to traditional cupping therapy. Each attachment features a convenient release button that allows for quick release of pressure from the cup after the massage, making removal quick, painless, and convenient.

After completing the massage, hold the POWER button again to turn off the device. The device will automatically turn off after 15 minutes of inactivity.

IMPORTANT:

- Adjust the massage intensity and temperature to your individual needs.
- Do not use the device on damaged skin, post-operative skin (up to 3 months after surgery), or in the event of inflammation.
- Avoid use in areas of high humidity, in the bathtub, shower, or near water.
- It is not recommended for use by individuals with varicose veins, hypertension, hypotension, stroke, or skin conditions – in such cases, consult a doctor before use.

The device is ideal for individuals experiencing muscle tension, fatigue, circulation problems, or cellulite. Regular use supports regeneration, improves skin condition, and improves well-being.

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Used electrical and electronic equipment and batteries and accumulators must not be mixed with other waste. This is indicated by the crossed-out bin symbol on the product packaging. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution. Proper handling of used electrical and electronic equipment and batteries and accumulators in the household contributes positively to the reuse and recovery, including recycling, of used equipment. Used electrical and electronic equipment and batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the used equipment is of the same type and has performed the same functions as the equipment sold. A distributor, when delivering household equipment to a purchaser, is obliged to take back used household equipment free of charge at the place of delivery of such equipment, provided that the used equipment is of the same type and has performed the same functions as the equipment delivered. A distributor running a retail outlet with a sales area within the meaning of art. 2 item 19 of the Act of 27 March 2003 on spatial planning and development of at least 400 m² devoted to the sale of household equipment, is obliged to accept, free of charge, in that outlet or in its immediate vicinity, used household equipment, none of the external dimensions of which exceeds 25 cm, without the need to purchase new household equipment. The product was introduced to the market after 13 August 2005

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The manufacturer declares, under its sole responsibility, that the Medivon Eclipse device complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances. The full content of the declaration of conformity is available at www.medivon.pl/deklaracje

Importer: Armare S.A. | ul Grzybowska 87, 00-844 Warszawa, Poland.
Model: Eclipse | **SN:** 07255904119286582 | **Marka:** Medivon
e-mail: serwis@medivon.pl

Producent: Zhejiang Xilian Electronic Technology Co., Ltd.
No.2, Yihuan Road, Xinxin Village, Kunyang Town, Pingyang County,
Wenzhou, 325499 Zhejiang, China.
e-mail: export@clavcare.com

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<http://medivon.pl/warranty>

Shortened warranty terms

This Guarantee is provided by Armare S.A., a company governed by Polish law, with its registered office in Warsaw, ul. Grzybowska 87, 00-844, Poland. You can find us in the official Polish register of companies under the KRS number: 0000730590, NIP: 5213833614 and REGON: 380721860.

We guarantee that the equipment you have purchased from us (hereinafter referred to as the "Product") will work properly. If something goes wrong for a reason inherent in the Product, this warranty protects you. The warranty starts from the day you receive the Product and lasts:

24 months, if you are a consumer or a natural person running a business and the purchase of the Product is not of a professional nature for you. Otherwise, the warranty period is 12 months. If the law in your area requires a longer warranty, the warranty period will be extended to meet local law requirements.

In addition, if the Product does not comply with the contract, you have the right to legal protection from and at the expense of the seller. This warranty does not affect these rights - it does not limit, suspend or exclude them.

Detailed warranty conditions, how to submit complaints, what happens when you find a defect, what are exceptions to the warranty, when the warranty expires or is extended, and what are the rules regarding the processing of personal data, you can find on our website: <http://medivon.pl/warranty>