RY200B/C

PL / BEZPRZEWODOWY IRYGATOR DO ZĘBÓW EN / CORDLESS ORAL IRRIGATOR
DE / KABELLOSE MUNDDUSCHE RU / БЕСПРОВОДНОЙ ИРРИГАТОР
ES / IRRIGADOR DENTAL INALÁMBRICO FR / HYDROPULSEUR SANS FIL
IT / IDROPULSORE SENZA FILI

Eldom Sp. z o.o. • ul. Pawła Chromika 5a • 40-238 Katowice, POLAND
tel: +48 32 2553340 • fax: +48 32 2530412 • www.eldom.eu

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/EN/ The disposal of used electrical and electronic equipment (applies to European Union countries and other European countries with separate waste-collection systems).

This symbol on the product or its packaging indicates that it should not be classified as household waste. It should be handed over to an appropriate company dealing with the collection and recycling of electrical and electronic equipment. The correct disposal of the product will prevent potential negative consequences for the environment and human health resulting from hazardous substances present in the product. Electrical devices must be handed over to restrict their re-use and further treatment. If the device contains batteries, remove them, and hand them over to a storage point separately. DO NOT THROW EQUIPMENT INTO THE MUNICIPAL WASTE BIN. Material recycling helps to preserve natural resources. For detailed information on how to recycle this product, please contact your local authority, the recycling company, or the shop where you bought it.

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SAFETY INSTRUCTIONS

- **Before first use, read this instruction manual carefully.**
- **This appliance may be used by children aged 8 years and above and by persons with reduced physical or mental capabilities, or lack of experience and knowledge, if they are supervised or instructed on how to use the appliance safely and understand the hazards involved. Children must not play with the appliance. Cleaning and maintenance must not be carried out by children without supervision.**
- **The appliance is intended for cleaning the oral cavity. Do not use it for any other purpose.**
- Do not place the appliance near heat sources.
- Always switch off the appliance before cleaning, maintenance, in case of malfunction, or after use.
- Do not use with damaged attachments, as this may cause injury.
- Do not use outdoors.
- Use only with the accessories supplied.
- Repairs may only be carried out by an authorized service center. Any modifications or the use of non-original spare parts or components are prohibited and may compromise safety.
- During charging, do not place the appliance where it could fall into a bathtub, washbasin, or shower.
- Do not touch the appliance if it has fallen into liquid during charging or while plugged in. In such a case, unplug it immediately.
- Eldom Sp. z o.o. is not liable for any damage resulting from improper use of the appliance;

SPECIAL INSTRUCTIONS

- **Persons with a pacemaker should consult a physician before using the appliance.**
- **Use the appliance only as described in this manual or as recommended by your dentist.**
- If you have undergone dental procedures within the last 2 months, consult your dentist before using the irrigator.
- Fill the appliance only with water at a temperature not exceeding 40°C.
- Do not direct the water jet under the tongue, into the ear, or towards other sensitive areas. The pressure generated by the appliance may cause injury to such areas.

- **Do not use mouth rinses containing iodine, bleach, or tea tree oil.**
- Remove oral jewelry before using the irrigator.
- Do not use the appliance for more than 5 minutes at a time or more frequently than every 2 hours.
- If bleeding occurs after using the irrigator, consult your dentist. Do not use the appliance if there are wounds in the oral cavity.
- Charge only with the charger specified in this manual.
- Do not operate the appliance without the nozzle (1) in place.
- Do not press the nozzle (1) too firmly against the gums, teeth, or tongue.

This appliance contains non-replaceable batteries.

- The battery must be removed from the appliance before disposal.
- The battery must be disposed of in a safe manner.
- Do not disassemble the product except for the purpose of battery disposal.

▪ **CLEANING AND MAINTENANCE**

- Clean the appliance after each use.
- Before cleaning, make sure the appliance is switched off.
- Do not immerse the cable (5) in water.
- Do not use strong or abrasive cleaning agents or products.
- Wipe the cable (5) and the housing (4) with a damp cloth.
- Rinse the nozzle (1) thoroughly. Do not squeeze it to avoid damage.
- To clean the inside of the tank, remove the bottom lid (11) by pressing the button (10).
- Remember to descale the appliance at least once every 3 months. For this purpose, use only descaling agents intended for dental irrigators.
- Do not clean the appliance by immersing it in water or placing it under running water, as this may cause irreversible damage to the appliance.

USER MANUAL

CORDLESS ORAL IRRIGATOR RY200B/C

EN

▪ GENERAL DESCRIPTION

1. Replaceable nozzle
2. Power switch
/ Mode selector
3. Mode indicator
4. Housing
5. USB cable
6. Socket
7. Nozzle lock
8. Charging port
9. Water inlet
10. Tank lid lock
11. Tank lid

▪ TECHNICAL DATA

Pressure: 100–160 PSI (7–11 Bar)
Pulsation: >1500 strokes/min
Capacity: 230 ml
Battery: DC 3.7V, 1500mAh, Li-ION
Waterproof rating: IPX7
Charging time: approx. 6 hours
Operating time: up to 15 days (2 times per day)
Charging method: wired - connect the cable (5) to a power adapter with the following parameters: **DC 5V, 0.4A (not included)**

▪ INTENDED USE

The irrigator is designed for daily oral hygiene. With different nozzles, it can be used to clean interdental spaces, hard-to-reach areas, gaps in orthodontic appliances, and the surface of the tongue. The appliance is also intended for massaging the gums. Please note that the irrigator is not a substitute for brushing teeth with a toothbrush and toothpaste but serves as a supplement. Regular use of the appliance may help reduce gum disease and significantly improve oral hygiene.

▪ OPERATING THE APPLIANCE

Before first use, read the entire content of this manual carefully. After unpacking the appliance, make sure it has not been damaged during transport. In case of any doubt, do not use the appliance until you have contacted the service center. Before first use, rinse the appliance with warm water.

Charging the appliance

- the appliance can only be used after charging,
- low battery is indicated by the control light (3) – after switching off the device, the indicator light will turn red,
- charge the appliance at room temperature,
- to extend battery life, the first 3 charging cycles should be completed without interruption, after the battery has been fully discharged,
- switch off the appliance before charging,
- connect one end of the cable (5) to the socket (8), and the other end to a power adapter (not included) with the parameters specified above, then connect it to the mains,
- during charging the indicator light (3) will be red, and it will turn green once fully charged,
- the appliance cannot be used while charging.

NOTE. In the case of a completely discharged battery, the indicator light (3) may turn on only after approx. 5 minutes. This is a normal phenomenon.

Types of nozzles and their applications

Nozzles (1) should be replaced every 3 months for hygienic reasons and immediately if damaged.

A – standard nozzle for irrigation.

B – tongue cleaning nozzle. Place the nozzle in the middle of the tongue, press gently, and then switch on the irrigator. It is recommended to start cleaning at the lowest level and gradually increase it, adjusting to your individual needs.

C – orthodontic nozzle. Designed for cleaning hard-to-reach areas of orthodontic appliances.

D – periodontal nozzle. Designed for cleaning periodontal pockets. It is also used for delivering medicinal liquids between the gum line and the periodontal pocket. This nozzle should only be used in gentle mode. To use it properly, place it on the tooth at an angle of about 45 degrees, then move it toward the base of the tooth (the point where the gum meets the tooth).

Nozzles (1) are consumable parts and are not covered by the warranty.

Tips (1) can only be replaced with original tips, available also in the online store www.eldom.eu.

▪ **USAGE**

- open the valve (9) or the lid (11) by pressing the lock (10),
- fill the tank with water or a specialized liquid – dental liquid should be used in a 1:1 ratio, when using other liquids follow the manufacturer's instructions, close the valve (9) or the lid (11).

NOTE. Fill the appliance only with water (liquids) at a temperature not exceeding 40°C.

- insert the appropriate nozzle (1) into the socket (6),
- to replace the nozzle, press and hold the button (7), then gently pull the nozzle (1) upwards.

NOTE. The nozzle can only be replaced when the appliance is switched off.

- use button (2) to select the operating mode:
 - **DIY mode** – after pressing button (2) the indicator will flash every 0.5 seconds for the first 3 seconds. By pressing the button, the intensity can be adjusted from minimum to maximum, and after releasing the button the setting will be saved automatically.
 - **strong** – for standard use,
 - **gentle** – for sensitive teeth,
 - **pulse** – for gum massage,
- place the nozzle (1) in your mouth, then press button (2) to switch on the appliance and select the operating mode,
- to switch off the appliance press button (2) again,
- after irrigation fill the tank with warm water and rinse the appliance – direct the nozzle (1) towards the washbasin/bathtub/shower.

▪ **TIPS**

- during irrigation, hold the appliance so that the nozzle (1) is directed upwards,
- keep your mouth open so that the water can freely flow out into the washbasin,
- place the nozzle (1) against the teeth at a 90-degree angle and move it clockwise,
- clean sensitive areas between the teeth carefully,
- move the water jet along the gum line – this will effectively rinse away food residues from both the surface and between the teeth, while at the same time providing a gentle gum massage.

PROBLEM	POSSIBLE CAUSE	SOLUTION
Appliance does not work	Battery discharged	Charge the battery
Incorrect pressure	Nozzle (1) clogged with scale	Replace the nozzle (1)
Leak from the appliance	Valve (9) or lid (11) not closed properly	Close the valve (9) or lid (11) correctly
Leak around socket (2)	Nozzle (1) not mounted correctly	Remove and insert the nozzle (1) properly into socket (2)

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