Panasonic

Operating Instructions
(Household use) Sonic Vibration Toothbrush
Model No. EW-DCA3

Series 300

English	3	Norsk	190
Deutsch	23	Suomi	210
Français	45	Polski	230
Italiano	66	Česky	251
Español	87	Slovensky	272
Nederlands	108	Magyar	293
Português	129	Română	315
Svenska	150	Türkçe	336
Dansk	170		

Operating Instructions

(Household use) Sonic Vibration Toothbrush

Model No. EW-DCA3

Contents

Safety precautions.....	4
Important information	8
Parts identification	9
Before cleaning your teeth	10
Features.....	10
Charging.....	10
Basic usage.....	12
Tips on usage	14
After use	14
Changing the brush head (Sold separately).....	16
FAQ	16
Troubleshooting	17
Battery life	20
Removing the built-in rechargeable battery.....	20
Specifications	21

Thank you for purchasing this Panasonic product.

Before operating this unit, please read these instructions completely and save them for future use.

WARNING

- Do not immerse the USB power adaptor and USB cable in water or clean them with water.
Doing so may cause electric shock or fire due to a short circuit.
- Always unplug the adaptor from a household outlet when cleaning the USB power adaptor and USB cable.
Failure to do so may cause electric shock or injury.
- Toothbrushes can be used by children and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Cleaning and user maintenance shall not be made by children without supervision.
- Children shall not play with the appliance.
- Do not use anything other than the supplied USB power adaptor or USB cable for any purpose. Also, do not use any other product with the USB power adaptor or USB cable. (See page 9.)
Doing so may cause burns or fire due to a short circuit.
- The following symbol indicates that a specific detachable power supply unit is required for connecting the electrical appliance to the supply mains.
The type reference of power supply unit is marked near the symbol.

Safety precautions

Make sure to follow these instructions.

To prevent accidents, injuries or damage to property, please follow the instructions below.

- The following chart indicates the degree of damage caused by improper operation.

 DANGER	Indicates a potential hazard that will result in serious injury or death.
 WARNING	Indicates a potential hazard that could result in serious injury or death.
 CAUTION	Indicates a potential hazard that could result in minor injury or property damage.

WARNING

► Power supply

- Do not use in a way that would exceed the rating of the household outlet or the wiring.
 - Exceeding the rating by connecting too many plugs to one household outlet may cause fire due to overheating.
- Fully insert the power plug of the USB power adaptor.
 - Failure to do so may cause fire or electric shock.

► About the USB power adaptor and USB cable

- Never use them if they are damaged or if the power plug fits loosely in a household outlet. Do not damage them.
- Do not place them over or near water filled sink or bathtub.
- Do not deface, modify, forcefully bend or pull them.
 - Also, do not place anything heavy on, twist or pinch them.
 - Doing so may cause electric shock or fire due to a short circuit.
- Do not connect or disconnect them to a household outlet with a wet hand.
 - Doing so may cause electric shock or injury.
- Clean regularly to prevent dust from accumulating.
 - Failure to do so may cause fire due to insulation failure caused by humidity.
 - Disconnect the USB power adaptor and wipe it with a dry cloth.
- Do not charge this product using any USB cable other than the supplied USB cable. Also, do not use any other product with the supplied USB cable.
 - Doing so may cause fire, electric shock, or injury.
 - The USB cable provides the overheat protection function.
 - If charging is abnormal or the Type-C plug heats up, disconnect the power supply, check the appliance socket, and remove any water or foreign objects from the appliance socket.

► In case of an abnormality or malfunction

- Immediately stop using this product and disconnect USB power adaptor if there is an abnormality or malfunction.
 - Failure to do so may cause fire, electric shock, or injury.
- <Abnormality or malfunction cases>
 - The main unit, USB power adaptor or USB cable is deformed or abnormally hot.
 - The main unit, USB power adaptor or USB cable smells of burning.
 - There is abnormal sound during use or charging of the main unit.
- Immediately request inspection or repair at an authorized service centre.

WARNING

► This product

- This product has a built-in rechargeable battery. Do not throw into flame or heat. Do not charge, use, or leave in high temperature environments.
 - Doing so may cause overheating, ignition, or explosion.
- Do not modify or repair.
 - Doing so may cause fire, electric shock, or injury.
 - Contact an authorized service centre for repair (battery change etc.).
- Never disassemble except when disposing of the product.
 - Doing so may cause fire, electric shock, or injury.
- Unplug the USB cable from the main unit when cleaning.
 - Failure to do so may cause electric shock or injury.

► Preventing accidents

- Do not place it within the reach of children or infants. Do not let them use it.
 - Doing so may cause an accident or injury due to accidental ingestion of removable parts.

CAUTION

► Protecting your teeth and gums

- Do not use the brush with toothpaste that contains a high abrasive content, such as toothpaste for removing tar or for whitening. Please use general toothpaste.
- Do not apply the brush head strongly against your teeth or gums, or use it continuously against a same spot.
- Do not apply the plastic part of the main unit or brush head to your teeth.
- People who are unable to operate should not use this appliance.
- People with weak oral senses should not use this appliance.
 - Doing so may cause injury or damage to teeth and gums.
- People who may have gum disease, have teeth being treated, or are concerned about symptoms within their mouth should consult a dentist before use.
 - Failure to do so may cause injury or damage to teeth and gums.

► This product

- Never use for a purpose other than cleaning your teeth.
 - Doing so may cause injury.
- Do not share the brush head with your family or other people.
 - Doing so may result in infection or inflammation.

CAUTION

► Note the following precautions

- Do not drop or hit against other objects.
 - Doing so may cause injury.
- Do not allow needles or any other foreign objects to adhere to the power plug of the USB power adaptor or the plugs of the USB cable.
 - Doing so may cause electric shock or fire due to a short circuit.
- Do not charge the appliance while the appliance socket is wet.
 - Doing so may cause fire or electric shock.
- Disconnect the USB power adaptor from the household outlet when not charging.
 - Failure to do so may cause electric shock or fire due to electrical leakage resulting from insulation deterioration.

Handling of the removed battery when disposing

DANGER

- The rechargeable battery is exclusively for use with this appliance. Do not use the battery with other products.
- Do not charge the battery after it has been removed from the product.
- Do not throw into fire or apply heat.
- Do not hit, disassemble, modify, or puncture with a nail.
- Do not let the positive and negative terminals of the battery get in contact with each other through metallic objects.
- Do not carry or store the battery together with metallic jewellery such as necklaces and hairpins.
- Do not use or leave the battery where it will be exposed to high temperatures, such as under direct sunlight or near other sources of heat.
- Never peel off the tube.
 - Doing so may cause overheating, ignition, or explosion.

WARNING

- After removing the rechargeable battery, do not keep it within reach of children and infants.
 - The battery will harm the body if accidentally swallowed.
If this happens, consult a physician immediately.
- If the battery fluid leaks out and comes in contact with the eyes, do not rub the eyes, and thoroughly rinse with clean water, such as tap water.
 - Failure to do so may cause eye injury.
Consult a physician immediately.


CAUTION

- If the battery fluid leaks out and comes in contact with the skin or clothes, rinse off with clean water, such as tap water.
 - Failure to do so may cause inflammation.

Important information

- Do not clean this product by using paint thinner, benzene, alcohol, or other chemicals. Otherwise, it may lead to product fracture or color fading. Soft cloth slightly dipped in tap water or suds can be used to clean this product's shell.
- After use, this product shall be kept away from the pool, shower or other areas exposed to water or humidity.
- This product shall be kept away from area of high temperature or direct sunlight.
- Charge the appliance correctly according to these Operating Instructions. (See pages 10-12 "Charging".)
- This appliance contains a battery that is only replaceable by skilled persons. Please contact an authorized service centre for the details of repair.

A Main unit

- ① Toothbrush shaft
- ② Power switch
(power off/on/mode selection)
- ③ Mode indicator (White)
- ④ Charge status indicator () (Red)
- ⑤ Appliance socket (Type-C)
- ⑥ Socket cover

B Extra-fine brush

- ⑦ Identification ring

C USB cable

- For charge only
- ⑧ Cord
- ⑨ USB plug (Type-A)
- ⑩ USB plug (Type-C)
- With overheat protection function
- ⑪ Overheat protection function mark

D USB power adaptor (ES-RU1-02)

- ⑫ Power plug
- ⑬ USB port (Type-A)

Before cleaning your teeth

► About this product

The Sonic Vibration Toothbrush is an electric toothbrush with sonic vibration*¹ of approximately 31 000 brush strokes/minutes*².

*¹ Sonic vibration in acoustic wave field.

*² Only using SOFT mode.

Brush strokes/minutes may vary depending on the mode selected.

► Bleeding of the gums

When you first begin using the toothbrush, you may experience slight gum bleeding even if your gums are healthy. This is because your gums are receiving stimulation by the toothbrush for the first time. The bleeding should stop after one or two weeks of use. Bleeding that persists for more than two weeks may be indication of an abnormal gum condition, such as gum disease. In this case, stop using the toothbrush and consult your dentist.

Features

► About soft start function

To get accustomed to the high speed vibration, the toothbrush will gradually reach the selected speed and vibration approximately 3 seconds after it starts operating.


Charging

1 Open the socket cover and insert the USB plug firmly into the main unit and USB power adaptor.

- Wipe any water or toothpaste residue off the appliance socket area.



2 Insert the USB power adaptor firmly into a household outlet.

- Charging is completed after approx. 2 hours.

3 Disconnect the USB power adaptor from the household outlet when the charge status indicator () lights off.

- Close the socket cover.

► When to charge

- When the appliance is switched on: the charge status indicator () will blink (once per second).
- When the appliance is switched off: the indicator lights off after 5 seconds.

► Even when the appliance battery has run out of charge

- Charging for 2 minutes will supply power for approx. 2-minute use. (Quick charge.)
(The available use time may vary depending on the ambient temperature and the manner in which the appliance is used.)

► Battery runtime on a full charge

Mode	Battery runtime
CLEAN SOFT	Approx. 31 days*
SENSITIVE	Approx. 90 days*

(The available use time may vary depending on the ambient temperature and the manner in which the appliance is used.)

* Using for 2 min. 2x per day.

► Notes


- The charging time may vary when using for the first time, or when it has not been in use for more than 6 months. Also, the charge status indicator () may not glow for 5 to 10 minutes after charging starts, but it will glow if charging continues.
- The required charging time may vary depending on the starting charge level.

- Recommended ambient temperature for charging is 10 °C – 35 °C. When the ambient temperature is out of the recommended range, the battery may not charge properly or may take longer to be fully charged.
- If the charge status indicator () blinks quickly (twice per second), see page 17.

Basic usage


1 Attach the brush head to the main unit.

There is a small gap between the brush head and the main unit (a few mm), but this is not a malfunction.

2 Apply toothpaste.

Put a small amount of toothpaste on the brush head.

3 Put the brush head in your mouth and press the power switch to turn on the main unit.

► Select the desired mode

To switch the mode, press the power switch within 5 seconds.


- The last mode used for more than 5 seconds will be stored.

4 Lightly press and move the toothbrush against your teeth.


► Example brushing sequence

- ① front of the teeth
- ② biting plane of the teeth
- ③ back of the teeth

► About 30-second brush timer

- The toothbrush will pause once at 30-second intervals to indicate the brushing time. It pauses twice at 2-minute intervals then starts again.

5 Turn off the main unit before removing it from your mouth.

Notes

- The appliance may not operate in an ambient temperature approximately 5 °C or lower.
- Use the appliance in the appropriate mode.

Mode	Strength	Usage or desired care
CLEAN	Strong  Weak	For everyday cleaning at standard intensity.
SOFT		For first-time users or for gentle brushing.
SENSITIVE		For extra gentle brushing to prevent gum irritation

Tips on usage


Extra-fine brush

When you want to brush mainly the gum pockets or between the teeth and gums.

How to apply/move the toothbrush


Front surface and the biting surface of teeth

- Apply the brush head at a 90° angle.

Gum pocket and the boundary between the teeth and gums

- Place the brush head at a 45° angle.
- Apply the bristles lightly into the gum pocket and move the toothbrush slowly.

For the back of the front teeth or misaligned teeth

- Place the brush vertically.
- Move toothbrush to match the angle of the teeth.

After use

- Do not use thinner, benzene, alcohol, etc. to clean the appliance. Doing so may cause failure, cracking or discoloration of the parts.

Remove the brush head.

Brush head

Wash the brush head under running water.

- Do not wash with water hotter than 80 °C.
Doing so may damage the brush head.

Main unit


Close the cover and wash with water

Before washing the toothbrush, push the cover down to see it is tightly closed.

Wash the main unit under running water.

- Do not damage or rub the waterproof rubber seal (A) with hard materials.
- Do not leave the main unit under water.
- Do not wash with water hotter than 40 °C.
Doing so may damage the main unit.

Wipe off water from the brush head and main unit for storage. (This is to keep the toothbrush clean and to reduce malfunction risks.)

- Do not charge the toothbrush when it is wet.


Replacing the identification rings

► Removing

Slide up with your thumb.


► Attaching

Insert by covering the indentation (a) in the brush head with the inside of the identification ring (b).

Changing the brush head (Sold separately)

- Brushes are consumables.
- The brush should be replaced with a new one every 3 months.

Part number	Replacement Brush for EW-DCA3	
WEW0801	Extra-fine brush	

FAQ

	Question	Answer
Write using	Why does my toothbrush temporarily stop or reduce the brush action about every 30-second?	► The toothbrush will notify you of the elapsed brushing time every 30-second. (See page 13.)
Battery	Will the battery deteriorate if not used for an extended period of time?	► The battery will deteriorate if not used for over 6 months, so make sure to perform a full charge at least once every 6 months.
	Can I charge the appliance after every use?	► The battery used is a lithium-ion, so charging after every use will have no effect on battery life.

	Problem	Possible cause	Action
While charging	The charge status indicator () does not glow when charging starts.	▶ You have just purchased the appliance or have not used it for a while. ▶	The charge status indicator () may not glow for 5 to 10 minutes. It will eventually glow if kept connected. If the charge status indicator () does not glow after a while, consult with an authorized service centre for repair.
		▶ The USB cable is not correctly connected with the main unit and/or the USB power adaptor. ▶	Make all connections correctly and check that the charge status indicator () lights on. (See page 11.)
	The charge status indicator () is blinking rapidly. (Twice every second)	▶ The USB power adaptor is not correctly connected to the household outlet. ▶	
		▶ The temperature is too low or high while charging. ▶	Charge the appliance within the recommended charging temperature of 10 °C to 35 °C. If the charge status indicator () blinks even when charging within the recommended temperature range, please contact an authorized service centre for repair.

	Problem	Possible cause	Action
While charging	▶ Toothbrush does not charge.	▶ The USB cable is not correctly connected with the main unit and/or the USB power adaptor. ▶ The USB power adaptor is not correctly connected to the household outlet.	▶ Make all connections correctly and check that the charge status indicator () lights on. (See page 11.)
	▶ The toothbrush can only be used for a few minutes, even if it is fully charged.	▶ The wet appliance socket activated the overheat protection feature of the USB cable to suspend the charging action.	▶ Allow the appliance socket to dry before recharging.
		▶ The rechargeable battery has reached its operating life. (Approx. 3 years)	▶ Replace the battery. Contact the authorized service centre where purchased.

	Problem	Possible cause	Action
While using	The vibration increases gradually when pressing the power switch.	▶ The “Soft start function” is in operation.	▶ To avoid exposing the user to rapid mechanical vibration all at once, the appliance takes 3 seconds after power-on to increase its vibration action to the brushing level. Afterward, the appliance will be ready for use.
	The charge status indicator () lights off without blinking.	▶ The rechargeable battery has reached its operating life. (Approx. 3 years)	▶ Replace the battery. Contact the authorized service centre where purchased.
	The appliance does not operate even after charging.	▶ The appliance is used with approximately 5 °C or lower.	▶ Do not use the appliance with approximately 5 °C or lower.

If the problems still cannot be solved, contact the store where you purchased the unit or a service centre authorized by Panasonic for repair.

Battery life

Approximately 3 years. The battery may have reached its operating life if the usage frequency decreases drastically, even after a full charge. (The operating life of the battery may differ depending on usage or storage conditions.)


Removing the built-in rechargeable battery

This product has a built-in rechargeable battery. Be sure to remove the battery before disposing of this product.

- The battery is to be disposed safely.
- This figure must only be used when disposing of this product and must not be used to repair it.
- If you dismantle this product yourself, it will no longer be waterproof, which may cause it to malfunction.
- We recommend seeking the advice of a dealer or Panasonic's authorized service centre before deciding to dispose of this product, so as to extend the service life of this product as much as possible.
- If you decide to dispose of this product but do not want to disassemble this product by yourself to remove the battery, we recommend taking this product to Panasonic's authorized service centre which will remove and recycle the battery.

► How to remove the built-in rechargeable battery

- Power off this product before removing the battery. Turn on the power and keep the power on until the battery is completely discharged.
- Perform steps ① to ⑤ to disassemble the main unit and remove the battery with a screwdriver.
 - ③ While pushing at the toothbrush shaft (See page 9.), push the two battery holder tabs (a).
 - ④ Insert the tip of the driver inward (b) and disconnect the battery terminal.
 - ⑤ Remove the battery by rotating it.
- Insulate and protect the ⊕ and ⊖ terminals of the battery by applying tape to them. (To prevent short circuit of terminals and to avoid surface cuts on your fingers caused by the cut sections of the terminals).
- If you cannot remove the battery, please contact an authorized service center.

For environmental protection and recycling of materials

This product contains a Li-ion battery.

Please make sure that the battery is disposed of at an officially assigned location, if there is one in your country.

Specifications

Power source	See the name plate on the product (Automatic voltage conversion)
Charging time	Approx. 2 hours (Room temperature: 10 °C – 35 °C)
Airborne Acoustical Noise	54 (dB (A) re 1 pW)

This product is intended for household use only.

For more information about the energy efficiency of the product, please visit our website, www.panasonic.com, and enter the model number in the search box.

Authorized Representative in EU:

Panasonic Marketing Europe GmbH
Panasonic Testing Centre
Winsbergring 15, 22525 Hamburg, Germany

Importer:

Panasonic Marketing Europe GmbH
Hagenauer Strasse 43, 65203 Wiesbaden, Germany

Manufactured by:

Panasonic Corporation
1006, Oaza Kadoma, Kadoma City, Osaka, Japan

Disposal of Old Equipment and Batteries

Only for the European Union and countries with recycling systems

These symbols on the products, packaging, and/or accompanying documents mean that used electrical and electronic products and batteries must not be mixed with general household waste. For proper treatment, recovery and recycling of old products and used batteries, please take them to applicable collection points in accordance with your national legislation.

By disposing of them correctly, you will help to save valuable resources and prevent any potential negative effects on human health and the environment.

For more information about collection and recycling, please contact your local authority. Penalties may be applicable for incorrect disposal of this waste, in accordance with national legislation.

Note for the battery symbol (bottom symbol):

This symbol might be used in combination with a chemical symbol. In this case, it complies with the requirement set by the Directive for the chemical involved.

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© Panasonic Corporation 2026
AB Yetkili Temsilcisi: